



[Manuals.plus](#) /

> [Oral-B](#) /

> Oral-B iO Series 10 Rechargeable Electric Toothbrush User Manual

Oral-B iO M10.4B4.2A

Oral-B iO Series 10 Electric Toothbrush User Manual

Model: iO M10.4B4.2A | Brand: Oral-B

INTRODUCTION

The Oral-B iO Series 10 Electric Toothbrush represents the pinnacle of oral care technology, designed to provide a professional clean feeling and enhanced gum health. This manual provides comprehensive instructions for setup, operation, maintenance, and troubleshooting to ensure optimal performance and longevity of your device.

IMPORTANT SAFETY INFORMATION

Please read all instructions carefully before using this product. Failure to follow these instructions could result in injury or damage to the product.

- **Not suitable for children under 3 years old.**
- Replace brush head every 3 months or sooner if brush head becomes worn.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

WHAT'S INCLUDED

- Oral-B iO Series 10 Electric Toothbrush Handle
- iO Sense Smart Charger with Real-time Feedback
- 4 Replacement Toothbrush Heads
- Charging Travel Case



Figure 1: Complete Oral-B iO Series 10 package contents, showing the Cosmic Black toothbrush handle, iO Sense charger, four brush heads, and charging travel case.

SETUP

Initial Charging

Before first use, fully charge your Oral-B iO Series 10 toothbrush. Place the toothbrush handle onto the iO Sense Smart Charger. The display on the handle and the charger will indicate the charging progress. A full charge typically takes approximately 3 hours.



Figure 2: The Oral-B iO Series 10 toothbrush and iO Sense charger as packaged, ready for initial setup.

Attaching the Brush Head

Align the brush head with the magnetic shaft on the toothbrush handle and gently push it into place until it clicks securely. The iO Series 10 uses a unique magnetic attachment for easy and secure brush head changes.

iOTM
SERIES 10



Oral-B's iconic round brush head plus gentle micro-vibrations

La brosette arrondie emblématique de Oral-B combinée aux microvibrations douces



SMART PRESSURE CONTROL

• Light ring changes color to help you optimize brushing pressure:

- Perfect
- Too Much
- Not Enough

• Customize your brush's light color in the Oral-B App

CONTRÔLE DE PRESSION INTELLIGENT

• L'anneau de lumière change de couleur pour vous aider à optimiser la pression de votre brossage :

- Pression parfaite
- Trop de pression
- Pas assez de pression

• Personnalisez la lumière de votre brosse à dents dans l'application Oral-B



INTERACTIVE COLOR DISPLAY

• 7 brushing modes
• Real-time personalized coaching

ÉCRAN INTERACTIF EN COULEURS

• 7 modes de brossage
• Conseils personnalisés en temps réel



ARTIFICIAL INTELLIGENCE

• 3D teeth tracking & position detection
• Bluetooth connected app

INTELLIGENCE ARTIFICIELLE

• Reconnaissance des dents en 3D et détection de la position
• Application à connexion Bluetooth



iO senseTM

Personalized live-coaching base

Base avec conseils personnalisés en temps réel



INTUITIVE GUIDING LIGHTS

• Real-time feedback on pressure & brushing coverage

VOYANTS LUMINEUX INTUITIFS

• Rétroaction sur la pression et le brossage en temps réel



WIFI CLOCK AND TIMER

HORLOGE ET MINUTEUR WIFI



MAGNETIC CHARGER

• Just ~3 hours for full charge

CHARGEUR MAGNÉTIQUE

• Seulement 3 heures environ pour une recharge complète



Trust the brush Dentists use the most worldwide

Faites confiance à la brosse la plus utilisée par les dentistes dans le monde entier

Figure 3: The toothbrush handle with different brush heads and the iO Sense charger, illustrating the magnetic attachment system.

Powering On/Off and Mode Selection

Press the power button to turn the toothbrush on. Use the mode selection button to cycle through the 7 available brushing modes: Daily Clean, Sensitive, Gum Care, Intense Clean, Super Sensitive, Tongue Clean, and Whitening. The interactive color display on the handle will show the selected mode.



Figure 4: The interactive display on the toothbrush handle, showing the 2-minute timer and mode selection options.

Pressure Sensor Guidance

The visible pressure sensor protects your gums by alerting you if you brush too hard or too softly. The light ring around the brush head will illuminate: green for just right pressure, red for too much pressure, and white for too light pressure. Adjust your brushing force accordingly.



Figure 5: Visual representation of the pressure sensor feedback, indicating optimal and excessive pressure.

iO Sense Smart Charger and App Integration

The iO Sense charger provides real-time feedback on brushing coverage, time, and pressure directly via the charger's display, even without opening the Oral-B app. For a more detailed experience, sync your toothbrush with the free Oral-B app via Bluetooth. The app offers live coaching, tracks your brushing style across 6 zones, and helps you set long-term oral health goals.

DENTIST-INSPIRED ROUND BRUSH HEAD



Figure 6: The iO Sense charger and Oral-B app working together to provide comprehensive brushing data.

[Official Product Video: Oral-B iO Series 10 Electric Toothbrush Overview](#)

Your browser does not support the video tag.

Video 1: This official video from Oral-B provides a comprehensive overview of the iO Series 10 Electric Toothbrush, highlighting its key features such as the iO Sense smart charger, interactive display, and micro-vibrating bristles for a professional clean.

MAINTENANCE

Cleaning the Toothbrush

After each use, rinse the brush head thoroughly under running water with the handle switched on. Remove the brush head and clean both the handle and the brush head separately. Wipe the handle with a damp cloth. Do not immerse the charger in water.

Replacing Brush Heads

For optimal cleaning results, replace your Oral-B iO brush head every 3 months, or sooner if the bristles become splayed or worn. The toothbrush display will also indicate when it's time for a replacement.



Figure 7: A pack of Oral-B iO Series 10 replacement brush heads, essential for maintaining optimal oral hygiene.

Charging and Travel

The iO Series 10 features a 3-hour fast charger. The charging travel case allows for convenient charging on the go, ensuring your toothbrush is always ready for use. The display on the handle always tracks your battery level.

Your browser does not support the video tag.

Video 2: This short preview video demonstrates the charging process of the Oral-B iO Series 10 toothbrush on its iO Sense charger, highlighting the quick charge capability.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Toothbrush not charging.	Improper placement on charger, power outlet issue.	Ensure toothbrush is correctly seated on the iO Sense charger. Check power connection to the outlet.
App not connecting to toothbrush.	Bluetooth is off, app not updated, interference.	Ensure Bluetooth is enabled on your device. Update the Oral-B app to the latest version. Restart both the toothbrush and your device.
Pressure sensor always red/white.	Brushing too hard/soft consistently.	Adjust brushing pressure until the light ring turns green. Refer to operating instructions for proper technique.

SPECIFICATIONS

- **Product Dimensions:** 7.01 x 3.94 x 9.96 inches; 1.72 Pounds
- **Item Model Number:** iO M10.4B4.2A
- **Batteries:** 1 Lithium Ion battery required (included)
- **Manufacturer:** Procter & Gamble
- **Age Range (Description):** Adult
- **Special Feature:** Pressure Sensor
- **Power Source:** Battery Powered
- **Item Firmness Description:** Soft

WARRANTY AND SUPPORT

Oral-B offers a 30-day money-back guarantee for carefree switching to their electric toothbrushes. For detailed warranty information, product support, or to register your product, please visit the official Oral-B website or contact their customer service. Contact details can typically be found on the product packaging or the Oral-B brand store.

Oral-B Store: [Visit the Oral-B Store on Amazon](#)