

BNF KD-7910

BNF® Kd-7910 Blue Light Acne Pen User Manual

Model: KD-7910 | Brand: BNF

1. INTRODUCTION AND PRODUCT OVERVIEW

The BNF® Kd-7910 Blue Light Acne Pen is a specialized device designed for skin treatment, focusing on acne, scars, and wrinkles. It utilizes a specific 415nm blue light in the visible spectrum, combined with 42°C heat and BIO micro electric current, to address various skin concerns.

Key functions include:

- Effectively reducing the grease of pimples and accelerating their natural resolution.
- Tightening loose skin and shrinking enlarged pores.
- Promoting lymph detoxification.
- Diminishing inflammation and reducing pain.

This device is applicable for oily skin, severe skin allergy, and red spots. It can also assist with spider veins, Rosacea, post-operative bruising, and age-related Purpura.



Figure 1: Overview of the BNF Kd-7910 Blue Light Acne Pen.

2. PRODUCT COMPONENTS

Familiarize yourself with the different parts of your device:



Figure 2: Labeled diagram of the device components.

Probe: The treatment head that emits blue light and heat, and applies micro electric current.

Lid: Protective cover for the probe.

LED Indicator: Illuminates to show device status (e.g., on/off, active treatment).

Power Switch: Button to turn the device on or off.

Conductive Strip: Area on the device body that needs to be in contact with your skin for the micro electric current to function.

Battery Case: Compartment for the battery.

3. SETUP

3.1 Battery Installation

1. Unscrew the bottom cap of the device to access the battery compartment.

2. Insert one (1) AA battery (not included) into the battery case, ensuring the correct polarity (+/-) as indicated inside the compartment.
3. Screw the bottom cap back on securely.

3.2 Initial Check

After installing the battery, press the Power Switch. The LED Indicator should light up, and the probe should begin to warm slightly, indicating the device is ready for use.

4. OPERATING INSTRUCTIONS

4.1 Preparation

- Ensure your skin is clean and dry before use.
- Remove the protective lid from the probe.

4.2 Usage Steps

1. Press the Power Switch to turn on the device. The LED indicator will illuminate, and the probe will start to heat up.
2. Place the probe directly onto the affected area of your skin (e.g., acne, scar, wrinkle). Ensure your fingers are in contact with the conductive strip on the device body for the micro electric current to activate.
3. Hold the device gently on the area for approximately 2-5 minutes per session. You may feel a slight warmth and tingling sensation, which is normal.
4. Move the device to the next target area and repeat the process.
5. After completing your treatment, press the Power Switch again to turn off the device. The LED indicator will turn off.

4.3 Recommended Usage

For optimal results, use the device once a day. Consistent use over several weeks can lead to noticeable improvements in the visibility of acne, scars, and thread veins.



Figure 3: The device emitting blue light during operation.

EFFECTIVE KILL BACTERIA IN PORES

NO Acne



before



after

NO Red Spot



before



after



Figure 4: Illustrative results of skin treatment (Note: Individual results may vary).

5. MAINTENANCE

5.1 Cleaning

- After each use, ensure the device is turned off.
- Wipe the probe and the device body with a soft, dry cloth. For thorough cleaning, a cloth lightly dampened with rubbing alcohol can be used, ensuring no liquid enters the device.
- Do not rinse the device under water or immerse it in any liquid.

5.2 Storage

- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- Replace the protective lid on the probe when not in use.
- If the device will not be used for an extended period, remove the battery to prevent leakage.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is dead or incorrectly inserted.	Replace with a new AA battery, ensuring correct polarity.
Blue light is not visible or weak.	Low battery or device malfunction.	Replace battery. If problem persists, contact customer support.
No warmth or tingling sensation.	Fingers not in contact with conductive strip, or device malfunction.	Ensure firm contact with the conductive strip. Check battery. If problem persists, contact customer support.
Skin irritation or discomfort.	Overuse, sensitive skin, or allergic reaction.	Discontinue use immediately. Consult a healthcare professional if irritation persists.

7. SPECIFICATIONS

Feature	Detail
Model	KD-7910
Item Model Number	BNF62030462
Voltage	DC 1.5V
Current	0.36A
Power	0.9W
Battery	1 x AA Battery (not included)
Material	Plastic
Device Dimensions (Length x Diameter)	13.5 cm x 2 cm
Product Dimensions (LxWxH)	8 x 5 x 3 cm
Country of Origin	China



Figure 5: Device dimensions (135mm length, 20mm diameter).

8. SAFETY INFORMATION

- This device is for external use only. Do not use internally.
- Keep out of reach of children.
- Do not use on broken, irritated, or infected skin.
- Avoid direct contact with eyes.
- Consult a healthcare professional before use if you are pregnant, have a pacemaker, or have any underlying medical conditions, especially skin-related conditions.
- Discontinue use if you experience any adverse reactions such as excessive redness, swelling, or pain.
- Do not attempt to disassemble or repair the device yourself.

9. WARRANTY AND SUPPORT

Specific warranty information for this product is not detailed in the provided documentation. For any product inquiries, technical support, or warranty claims, please contact the manufacturer/importer directly:

Manufacturer/Importer: NALINISHUBHRA ENTERPRISE LLP

Customer Care: +91-9873384373

Email: nalinishubhra@gmail.com

© 2023 BNF. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How can I request a refund for the BNF KD-7910 Blue Light Acne Pen if it is not working for me?

Manuals+ is a documentation repository and does not sell products or process refunds directly. To request a refund for your BNF KD-7910 Blue Light Acne Pen, you must contact the retailer where you purchased...

[Read full answer](#)