

[Manuals.plus](#) /

› [NCVI](#) /

› NCVI Hands Free Breast Pump User Manual

NCVI 8791

NCVI Hands Free Breast Pump User Manual

Model: 8791 | Brand: NCVI

PRODUCT OVERVIEW

The NCVI Wearable Hands-Free Breast Pump is designed for convenience and discretion, allowing you to pump anywhere, anytime. This electric breast pump features multiple modes and intensity levels to optimize milk production, mimicking a baby's natural suckling. Its lightweight and slim design ensures comfortable wear inside your bra, and its low noise operation allows for discreet use without disturbing your baby or others. The pump is USB rechargeable for on-the-go use and is easy to clean with detachable accessories.



Image: Two NCVI Hands Free Breast Pumps, showcasing their compact and wearable design.

WHAT'S INCLUDED

Your NCVI Hands Free Breast Pump package includes the following components:

- 2 x NCVI Hands Free Breast Pumps (main units)
- 1 x USB Charging Cable
- 2 x 28mm Flanges (pre-installed)
- 2 x 21mm Flange Inserts
- 2 x 24mm Flange Inserts
- Instruction Manual

SETUP INSTRUCTIONS

Initial Preparation

1. **Wash Hands:** Before handling any breast pump components, thoroughly wash your hands with soap and water.
 2. **Clean Components:** Although the pump comes fully assembled, it is recommended to disassemble and clean all parts that come into contact with breast milk before first use. Refer to the "Maintenance and Cleaning" section for detailed instructions.
 3. **Charge Battery:** Ensure the breast pump is fully charged before initial use. Connect the provided USB cable to the charging port located on the side of the pump and plug it into a suitable USB power source.
 4. **Select Flange Size:** The pump comes with 28mm flanges pre-installed, along with 21mm and 24mm inserts. Choose the flange insert size that best fits your nipple for optimal comfort and milk expression. To change the insert, gently pull out the existing one and press the desired size into place within the flange.
-

OPERATING THE BREAST PUMP

Powering On/Off

- To turn on the pump, press and hold the power button (usually indicated by a circular icon) on the top of the unit until the indicator light illuminates.
- To turn off the pump, press and hold the power button again until the indicator light turns off.

Mode Selection

The NCVI breast pump offers multiple modes to suit your needs. Press the 'M' button to cycle through the available modes:

- **Massage Mode:** Designed to stimulate milk let-down with gentle, rapid suction patterns.
- **Normal Mode:** Provides a steady, consistent suction for efficient milk expression.
- **Pacifier Option (Simulation Mode):** Mimics a baby's natural suckling rhythm and pattern.
- **Combination Mode:** Alternates between light and stronger suction patterns.

Intensity Adjustment

Each mode has 9 adjustable intensity levels. Use the '+' and '-' buttons to increase or decrease the suction strength to your comfort level.

Milk Collection

The transparent milk collection cup has measurement markings on the side. When the milk reaches the maximum fill line, the pump should be paused, and the milk transferred to a storage container to prevent overflow.

CARE AND MAINTENANCE

Proper cleaning and maintenance are crucial for the longevity and hygiene of your breast pump. All parts that come into contact with breast milk are BPA-free and easy to clean.

Disassembly for Cleaning

1. Separate the pump motor from the milk collection cup by pressing the release button on the front of the pump.
2. Twist off the lid of the milk collection cup.
3. Carefully pull out the diaphragm and flange insert from the main flange unit.
4. Disassemble any other removable parts as indicated in the detailed manual.

Cleaning Guidelines

- Wash all disassembled parts (excluding the pump motor and USB cable) with warm, soapy water immediately after each use.
- Rinse thoroughly with clean water.
- Air dry on a clean drying rack or use a clean cloth. Ensure all parts are completely dry before reassembly.
- The pump motor and USB cable should not be submerged in water. Wipe them clean with a damp cloth if necessary.

Reassembly

Once all parts are dry, reassemble the pump by reversing the disassembly steps. Ensure all components are securely fitted to maintain proper suction.

TROUBLESHOOTING COMMON ISSUES

If you encounter any issues with your NCVI Hands Free Breast Pump, please refer to the following general advice. For more detailed troubleshooting or persistent problems, consult the full instruction manual or contact customer support.

- **No Power:** Ensure the pump is fully charged. Check the USB cable and power source connection.
 - **Weak Suction:** Verify all pump components are correctly assembled and securely sealed. Check that the flange insert size is appropriate for your nipple. Ensure the diaphragm is properly seated and not damaged.
 - **Leaking:** Confirm all parts are correctly assembled and tightened. Ensure the milk collection cup is not overfilled.
 - **Unusual Noise:** Disassemble and reassemble the pump to ensure all parts are correctly aligned. If the noise persists, contact customer support.
-

PRODUCT SPECIFICATIONS

Product Dimensions: 8.9 x 6.5 x 3.4 inches; 1.87 Pounds

Item Model Number: 8791

Batteries: 2 Lithium Ion batteries required (included)

Manufacturer: NCVI

Operation Mode: Electric

Special Features: Double Pump, Ergonomic, Hands Free Usage, Portable, Wearable

Material: Polypropylene

Power Source: Battery Powered

IMPORTANT SAFETY INFORMATION

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

Always consult with a healthcare professional before using any breast pump, especially if you have pre-existing medical conditions or concerns. Ensure proper hygiene practices are followed to prevent contamination.

WARRANTY AND CUSTOMER SUPPORT

For warranty information, product support, or any inquiries not covered in this manual, please contact NCVI customer service through the retailer's platform or the official NCVI website. Keep your proof of purchase for warranty claims.

© 2024 NCVI. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question