

[Manuals.plus](#) /

› [JFMkaer](#) /

› JFMkaer Compression Sock Aid Device Instruction Manual

JFMkaer B0CCYPY28J

JFMkaer Compression Sock Aid Device Instruction Manual

Model: B0CCYPY28J

1. INTRODUCTION

This manual provides detailed instructions for the safe and effective use of your JFMkaer Compression Sock Aid Device. Please read this manual thoroughly before initial use and retain it for future reference. This device is designed to assist individuals in putting on compression stockings without excessive bending or strain.

The JFMkaer Compression Sock Aid Device is suitable for calves with a circumference of **16 inches or less**. It features a patented hard plastic sock puller and ergonomic handles for ease of use.



Image 1.1: The JFMkaer Compression Sock Aid Device shown with a pair of donning gloves.

2. COMPONENTS AND FEATURES

- **Main Frame:** The core structure that holds the sock.
- **Adjustable Handles:** Extendable handles for comfortable use while sitting or standing.
- **Sock Puller:** Hard plastic component designed to hold the sock in place.
- **Donning Gloves:** Included to assist with handling compression socks and reducing friction.

EASY TO USE



Image 2.1: Detailed view of the device highlighting its foldable handle, high-quality hard plastic, and ergonomic tilt angle design. The device features foldable handles for convenient storage and a robust hard plastic construction for durability and effective sock application.

3. SETUP AND SIZING

3.1 Unboxing and Initial Inspection

Upon receiving your JFMkaer Compression Sock Aid Device, carefully remove all components from the packaging. Inspect the device for any signs of damage. Ensure all parts, including the main frame and donning gloves, are present.



TIPS: Includes Instruction Manual & Gloves
Please Check And Read Before Use

Image 3.1: The product packaging, showing the device and included gloves. An instruction manual is also typically included.

3.2 Sizing Compatibility

This device is designed for individuals with a calf circumference of **16 inches (approximately 40.6 cm) or less**. To determine if this device is suitable for you, measure the widest part of your calf. Using the device with a calf circumference exceeding this limit may result in difficulty or damage to the device or stockings.

SIZING INSTRUCTIONS



1. Please read the attached instructions before use



2. For Calf Circ. Under 16 inches

Image 3.2: Visual guide for measuring calf circumference to ensure compatibility with the device.

4. OPERATING INSTRUCTIONS

Follow these steps to effectively use the JFMkaer Compression Sock Aid Device:

1. **Prepare the Device:** Place the sock aid device on a stable surface, such as the floor or a low table. Ensure the handles are extended to a comfortable height if you are standing, or folded if you prefer to use it while seated.
2. **Load the Sock (Step 1):** Gather your compression sock and, using the provided donning gloves, stretch the sock over the hard plastic sock puller. Pull the sock down until the heel pocket is aligned with the curved part of the device, and the toe of the sock is at the bottom. The sock should be inverted over the device.
3. **Position Your Foot (Step 2):** While seated or standing, carefully insert your toes into the opening of the sock. Ensure your foot is centered.
4. **Pull the Sock On (Step 3):** Using the handles, gently pull the device upwards, guiding your foot into the sock. As you pull, the sock will slide smoothly onto your foot and up your leg. Continue pulling until the sock is fully on your leg and the device can be easily removed.

5. **Adjust the Sock:** Once the sock is on, use your hands (with gloves) to smooth out any wrinkles or bunched areas, ensuring a comfortable and proper fit.

Easy to use:

Flexible - Even if you have weak hand strength



step 1



step 2



step 3

NOTE : Follow the video a few times and you will be able to use the JFMkaer sock puller.

Image 4.1: Visual demonstration of the three main steps for putting on a compression sock using the device.

This device is suitable for various types of hosiery, including compression stockings, sports stockings, and silk stockings, minimizing the need to bend.

No Bending - Suitable For Compression Stockings, Sports Stockings And Silk Stockings



Image 4.2: The device positioned next to a bed, illustrating its use for various stockings without bending.

5. MAINTENANCE AND CARE

5.1 Cleaning

To clean the device, wipe it with a damp cloth and mild soap. Avoid using abrasive cleaners or harsh chemicals, as these may damage the plastic. Ensure the device is completely dry before storing.

5.2 Storage

When not in use, fold the handles inward for compact storage. Store the device in a cool, dry place away from direct sunlight and extreme temperatures to prolong its lifespan.

6. TROUBLESHOOTING

- **Difficulty Loading the Sock:** Ensure the sock is fully inverted over the device and the heel pocket is correctly aligned. Use the donning gloves to get a better grip on the sock material.
- **Sock Not Sliding Up Easily:** Verify that your calf circumference is within the recommended 16-inch limit. If the sock is too tight for the device, it may not slide smoothly. Ensure your foot is pointed and relaxed when inserting it.
- **Device Feels Unstable:** Ensure the device is placed on a flat, stable surface. If using while standing, maintain balance and use the handles for support, not as a primary balance aid.
- **Cracking Sounds or Damage:** If you hear cracking sounds or notice any damage, discontinue use immediately. This may indicate stress on the plastic components.

7. SPECIFICATIONS

- **Model:** B0CCYPY28J
- **Brand:** JFMkaer
- **Product Dimensions:** 4.49 x 9.02 x 12.99 inches (11.4 x 22.9 x 33 cm)
- **Weight:** 2.12 Pounds (approximately 0.96 kg)
- **Material:** High-quality hard plastic
- **Calf Circumference Compatibility:** Up to 16 inches (40.6 cm)

8. WARRANTY AND SUPPORT

8.1 Warranty Information

Specific warranty details for the JFMkaer Compression Sock Aid Device are not provided in this manual. Please refer to the product packaging or contact the retailer/manufacturer for warranty information.

8.2 Customer Support

For further assistance, questions, or concerns regarding your JFMkaer Compression Sock Aid Device, please contact the seller or manufacturer directly. You may also find additional resources by scanning the QR code provided on the product packaging or visiting the manufacturer's website.

Additional resources may be available at: <https://qr17.cn/EYU12V>

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for

conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question