

[Manuals.plus](#) /

› [Lenmen](#) /

› Lenmen Graphene-Infused Heating Pad User Manual

## Lenmen HP006

# Lenmen Graphene-Infused Heating Pad User Manual

Model: HP006

### INTRODUCTION

The Lenmen Graphene-Infused Heating Pad is designed to provide targeted heat therapy for deep pain relief and muscle relaxation, particularly for the neck and shoulders. Its weighted design ensures optimal contact, while advanced graphene technology promotes deeper heat penetration. This manual provides comprehensive instructions for safe and effective use, setup, operation, and maintenance of your heating pad.



Figure 1: Lenmen Graphene-Infused Heating Pad with controller and gift box.

## KEY FEATURES

- **Extended Heating Coil:** Features a 354" heating coil, 150% longer than ordinary heating pads, for faster and wider heat distribution across the entire pad within a minute.
- **Graphene-Infused Technology:** Utilizes AdvInfuse technology to incorporate graphene material, allowing infrared rays to penetrate subcutaneous tissues for deeper muscle relaxation and more effective soothing.
- **Ergonomic Weighted Design:** Designed with vertical and horizontal wiring to conform to body contours. Filled with keep-warm ceramic beads at the edges, ensuring the pad drapes naturally and stays in place, providing constant warmth.
- **Customizable Heat & Timer Options:** Offers 6 temperature settings (108°F ~ 144°F) and 4 auto-off timer options (0.5h ~ 2h) for personalized comfort.
- **Soft & Machine-Washable Fabric:** Made from micro plush flannel fabric for a luxurious and soft feel. The pad is machine-washable after removing the controller for easy cleaning.

- **Versatile Application:** Suitable for use on the neck, shoulders, upper back, trapezius, waist, abdomen, legs, and arms.

## UPGRADED HEATING WIRE

**OURS**



**354"** ✓

Heating Faster

Denser Heating Wiring

Heating Covers All Pad Area



**Safer Ptc/Ntc Heating Wire**

**OTHERS**



**236"** ✗

Heating Slowly

Sparse Heating Wiring

Heating Just Part Of The Pad



**Single-Layer Heating Wire**

Figure 2: Comparison of Lenmen's upgraded heating wire for faster and denser heat distribution.

# GRAPHENE INFUSED

Relax deeper muscles, Improve blood circulation, Boost metabolism



Figure 3: Graphene infusion allows heat to penetrate deeper into muscle tissues for enhanced relief.

# MORE THOUGHTFUL DESIGN

FURTHER COMFORTABLE EXPERIENCE




Figure 4: The thoughtful design with ceramic beads ensures optimal comfort and fit.

# MORE TEMPERATURE & TIMER CHOICES

## CUSTOMIZED WARMTH!



Figure 5: The LED controller allows easy selection of 6 heat levels and 4 auto-off timer options.

### SAFETY INFORMATION

To ensure safe operation of your Lenmen Graphene-Infused Heating Pad, please adhere to the following safety guidelines:

- Do not use while sleeping.
- Keep away from water and damp areas.
- Avoid covering the heating pad with heavy blankets or other materials that may cause overheating.

- Do not use on sensitive skin areas, open wounds, or areas with poor circulation.
- Consult a healthcare provider before use if you are pregnant or have a medical condition.
- Unplug the heating pad from the power outlet when not in use.
- Ensure the voltage matches your power supply (120 volts).

## SETUP

Before first use, unpack the heating pad and ensure all components (heating pad, controller, power cord) are present and undamaged. Lay the heating pad flat on a clean, dry surface.

1. Connect the controller to the heating pad. Ensure the connection is secure.
2. Plug the power cord into a standard 120V electrical outlet.
3. Position the heating pad on the desired body area (neck, shoulders, back, etc.). The weighted edges will help it conform and stay in place.

# ENVELOPED IN WARMTH

## EVERY COLD DAY



Figure 6: The heating pad can be used on various body parts for targeted warmth.

### OPERATING INSTRUCTIONS

To use your Lenmen Graphene-Infused Heating Pad:

1. Press the power button on the LED controller to turn on the heating pad.
2. Select your desired heat setting using the temperature control button. There are 6 temperature levels ranging from 108°F to 144°F.
3. Choose an auto-off timer option (0.5h, 1h, 1.5h, or 2h) using the timer button. The pad will automatically shut off after the selected duration.

4. Ensure the heating pad is positioned to allow even distribution of heat to the desired area.
5. Avoid prolonged use or falling asleep with the heating pad on, even with the auto-off timer.

Your browser does not support the video tag.

Video 1: Demonstration of the Lenmen Heating Pad's features and operation, including heat and timer settings.

## MAINTENANCE

Proper maintenance ensures the longevity and hygiene of your heating pad:

- Always disconnect the controller from the heating pad before cleaning.
- The micro plush flannel fabric is machine-washable. Wash in cold water on a gentle cycle.
- Do not bleach, iron, or dry clean the heating pad.
- Tumble dry on low heat or air dry. Ensure the pad is completely dry before reconnecting the controller and storing or using.
- Store the heating pad in a cool, dry place when not in use.

## TROUBLESHOOTING

If you encounter issues with your heating pad, refer to the following common troubleshooting steps:

Problem	Possible Cause	Solution
Heating pad does not turn on.	Not properly plugged in; power outlet issue; controller malfunction.	Ensure power cord is securely plugged into both the pad and the outlet. Test the outlet with another device. Contact customer support if the issue persists.
Uneven heating or no heat in certain areas.	Internal wiring damage; improper placement.	Adjust the pad's position to ensure full contact. If the issue persists, discontinue use and contact customer support.
Pad shuts off unexpectedly.	Auto-off timer activated; overheating protection; loose connection.	Check the timer setting. Ensure the pad is not covered by heavy materials. Re-secure all connections.

## SPECIFICATIONS

Feature	Detail
Brand	Lenmen

Feature	Detail
Model Number	HP006
Color	Blue
Material	Cotton (Micro Plush Flannel)
Package Dimensions	8.62 x 7.95 x 4.17 inches
Item Weight	1 Kilograms (approx. 2.2 Pounds)
Voltage	120 volts
Heat Settings	6 (108°F ~ 144°F)
Auto-Off Timer	4 options (0.5h, 1h, 1.5h, 2h)

## WARRANTY AND SUPPORT

Lenmen products are manufactured with quality and customer satisfaction in mind. For any warranty claims, technical support, or inquiries regarding your Graphene-Infused Heating Pad, please contact Lenmen customer service directly through their official channels. Refer to the product packaging or the Lenmen brand store on Amazon for specific contact information.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is intended for personal comfort and temporary relief of muscle pain and stiffness.