

INIA T033KW

INIA HAUTE IPL Laser Hair Removal Device Instruction Manual

Model: T033KW | Brand: INIA

INTRODUCTION

This manual provides comprehensive instructions for the safe and effective use of your INIA HAUTE IPL Laser Hair Removal Device. Please read all instructions carefully before operation to ensure optimal results and user safety. The device is designed for at-home hair reduction, offering advanced features like Double-Pulse technology and Sapphire Ice-Cooling for a comfortable and efficient experience.

IMPORTANT SAFETY INFORMATION

- For external use only. Do not use on broken or irritated skin.
- If any irritation or adverse reactions occur, discontinue use and consult a healthcare professional.
- Keep out of reach of children.
- Not suitable for tattooed skin, wounded areas, and pregnant individuals.
- Always wear protective goggles during use to shield your eyes from intense light.
- Ensure the device is used on clean, shaved skin, free from deodorants, lotions, or other products.
- Perform a patch test on a small skin area before full treatment to check for any adverse reactions.

SKIN TONE & HAIR COLOR COMPATIBILITY

The INIA HAUTE IPL device is most effective on lighter skin tones with darker hair. It is not suitable for all skin and hair types. Please refer to the charts below to determine your eligibility for safe and effective IPL use.



Image: Skin Tone and Hair Color Chart. This chart illustrates suitable skin tones (Fair, Medium, Olive/Light, Brown) and hair colors (Black, Dark Brown, Brown, Light Brown) for effective IPL treatment. It also indicates unsafe skin tones (Dark Brown, Black) and ineffective hair colors (Red, Light Blonde, White/Gray).

- **Suitable Skin Tones:** Fair, Medium, Olive/Light, Brown.
- **Unsuitable Skin Tones:** Dark Brown, Black.
- **Suitable Hair Colors:** Black, Dark Brown, Brown, Light Brown.
- **Ineffective Hair Colors:** Red, Light Blonde, White/Gray.

KEY PRODUCT FEATURES

- **Double-Pulse Technology:** Emits two pulses within milliseconds for ultra-fast hair removal. This maximizes energy absorption by hair follicles, offering superior effectiveness compared to single-pulse mode while minimizing skin harm.
- **Rapid 5 Flashes per Second:** Features an Auto-Glide mode for quick and effortless full-body treatment, reducing missed spots.
- **Ultra Sapphire Ice-Cooling:** Maintains a cool 50°F (10°C) for up to 20 minutes, providing a refreshing icy sensation for ultimate comfort and protecting the skin from heat injury during treatment.
- **5 Intensity Levels:** Allows for customization to target various hair types and skin sensitivities.
- **Safety Certifications:** Rigorously tested and holds prestigious safety certifications for peace of mind.

SETUP INSTRUCTIONS

1. **Prepare Skin:** Clean and shave the area to be treated. Ensure the skin is dry and free from any deodorants, lotions, or other products.
2. **Connect Power:** Plug the device into a power outlet using the provided adapter. The device is corded electric and must be plugged in to operate.
3. **Wear Goggles:** Put on the protective IPL goggles included in the package to shield your eyes from the intense light flashes.
4. **Power On:** Long press the power button for 2 seconds to turn on the device. The fan inside will activate, and you may hear its sound. The ice cooling symbol will light up by default.

Powerful 21J Output for Quicker,
More Efficient Hair Removal
with Fewer Sessions



Image: Visual guide for device setup, showing steps to prepare skin, wear goggles, plug in the device, and begin hair removal.

OPERATING INSTRUCTIONS

1. **Select Intensity Level:** Short press the level button to adjust the treatment level. Level-1 is the lowest intensity, and Level-5 is the highest. Choose a suitable level appropriate for your skin and hair type. Start with Level-1 and gradually increase as comfortable.
2. **Position Device:** Place the light-emitting window of the device firmly and vertically against the skin area you wish to treat. Ensure full contact with the skin; the indicator lights will illuminate when properly positioned.
3. **Initiate Flash:**
 - **Manual Flash Mode:** Press the flash button to emit a single pulse. This mode is ideal for smaller, precise areas like the upper lip or underarms.
 - **Auto-Glide Mode:** Long press the flash button for 2 seconds to activate Auto-Glide mode. The device will continuously emit flashes as you glide it over the skin, ideal for larger areas like legs or back.
4. **Treatment Frequency:** For the first few weeks, use the device approximately 3 times a week. After initial hair reduction, reduce frequency to once every 2-4 weeks for maintenance.
5. **After Treatment:** After finishing the treatment, keep pressing the power button for 2 seconds to turn off the device. Allow the device to cool to room temperature.

Your browser does not support the video tag.

Video: Demonstration of the Blex Laser Hair Removal device in Black, showing how to shave, plug in, turn on, wear goggles, and apply the device to the skin for hair removal. This video illustrates the Auto-Drift and Dual-Pulse modes, and emphasizes holding the device vertically against the skin.

MAINTENANCE

- After each use, ensure the device has cooled down.
- Use a dry and clean cloth to gently wipe the device, especially the light-emitting window, to remove any residue.
- Store the device and its accessories in the original gift box or a safe, dry place at room temperature.

TROUBLESHOOTING

Error reminders may be encountered, indicated by the level indicator light when the device is running. Refer to the following:

- **ERROR 1:** Indicates a reminder of input voltage abnormality. When the voltage is lower than the value, or exceeds the highest value, the first light (far away from light-emitting window) of the level indicator light will light up. Customers need to contact the after-sales department.
- **ERROR 2:** Indicates a reminder of fan abnormality. When the device is running, the second light of the level indicator light will light up. Customers need to contact the after-sales department.
- **Device Not Flashing:** If the white light isn't on when positioned, adjust the device's position to ensure full contact with the skin.

PRODUCT SPECIFICATIONS

- **Model Name:** T033KW
- **Power Source:** Corded Electric
- **Special Feature:** Sapphire Ice-Cooling
- **Recommended Uses:** Hair Removal
- **Product Dimensions:** 10.43 x 6.88 x 4.44 inches
- **Weight:** 2.56 Pounds
- **Manufacturer:** INIA

WARRANTY AND SUPPORT

Your INIA HAUTE device comes with a **1-year warranty**, extendable to **2 years upon registration**. Free replacement handsets are provided if you run out of flashes within the warranty period. For any inquiries, concerns, or to contact the after-sales department, please refer to the contact information provided in your product packaging or visit the official INIA website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Email

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Details

Submit question