

[Manuals.plus](#) /

› [LifePro](#) /

› LifePro Portable Red Light Therapy Torch User Manual

LifePro LumiCure Pro

LifePro LumiCure Pro Light Therapy Torch

User Manual

INTRODUCTION

The LifePro LumiCure Pro Light Therapy Torch is a portable and lightweight device designed for convenient light therapy sessions. This compact device utilizes specific wavelengths of light to support various applications, offering flexibility for use at home or while traveling. Regular short sessions are recommended for optimal results.



Image: The LumiCure Pro Light Therapy Torch alongside its zippered carrying case.

WHAT'S IN THE BOX

Upon opening your LumiCure Pro package, please ensure all components listed below are present:

- LifePro LumiCure Pro Light Therapy Torch
- Zippered Carrying Case
- Torch Pouch
- Wrist Strap
- Cleaning Cloth
- Battery Charger
- Rechargeable Batteries (2)
- User Manual (this document)

WHAT'S IN THE BOX



Image: An overview of all items included in the LumiCure Pro product box.

SETUP

1. **Charge Batteries:** Before first use, fully charge both rechargeable batteries using the provided battery charger. Insert the batteries into the charger according to the polarity indicators. The charger's indicator light will show when charging is complete.
2. **Insert Batteries:** Unscrew the cap at the base of the LumiCure Pro torch. Insert the charged batteries into the battery compartment, ensuring correct polarity (+/-). Screw the cap back on securely.
3. **Attach Wrist Strap (Optional):** If desired, attach the wrist strap to the designated loop on the torch for added security during use.

OPERATING INSTRUCTIONS

The LumiCure Pro offers multiple light modes for various applications. Familiarize yourself with the device's operation for effective use.

1. **Power On/Off:** Press and hold the power button located on the side of the torch for approximately 2 seconds to turn the device on or off.
2. **Change Light Mode:** Once powered on, short press the power button to cycle through the available light modes. The LumiCure Pro features:
 - **Red Light Therapy (630nm & 660nm):** Effective for various body applications.
 - **Near-Infrared Light (850nm & 900nm):** Provides natural support for different applications.
 - **Combined Red and Near-Infrared Light:** Utilizes both wavelengths simultaneously.
 - **10Hz Pulse Mode:** This mode introduces short pauses between light pulses, which may encourage cellular regeneration and recovery. Activate this mode by cycling through the light options.
3. **Session Duration:** For general use, a 2- to 3-minute session per area is recommended. Adjust duration based on individual needs and professional advice.
4. **Application:** Hold the torch directly over the desired area of the body. The device is designed for direct skin contact or close proximity.

IMPROVE YOUR SESSIONS WITH 10HZ PULSE MODE

Research Shows short pauses between each light pulse encourage cell regeneration & recovery!



Image: A user demonstrating the application of the LumiCure Pro Light Therapy Torch on their arm.

MAINTENANCE

- **Cleaning:** Use the provided cleaning cloth or a soft, dry cloth to wipe the surface of the torch after each use. Do not use abrasive cleaners or immerse the device in water.
- **Battery Care:** Recharge batteries when the device's performance diminishes. Store the device and batteries in a cool, dry place when not in use for extended periods. Remove batteries if storing for several months.
- **Storage:** Store the LumiCure Pro in its zippered carrying case to protect it from dust and damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Batteries are not charged or incorrectly inserted.	Ensure batteries are fully charged and inserted with correct polarity. Press and hold the power button for 2 seconds.
Light output is dim.	Low battery charge.	Recharge the batteries.
Device turns off unexpectedly.	Low battery charge or loose battery connection.	Recharge batteries. Check that the battery cap is securely tightened.

SPECIFICATIONS

- **Brand:** LifePro
- **Model:** LumiCure Pro
- **Light Source Type:** Infrared (IR)
- **Wavelengths:** 630nm, 660nm (Red Light); 850nm, 900nm (Near-Infrared Light)
- **Power Source:** Battery Powered (Rechargeable)
- **Color:** Black Lumicure Pro
- **Special Feature:** Portable, Lightweight
- **Package Dimensions:** 9.2 x 6.75 x 2 inches
- **Weight:** 1.37 Pounds
- **Manufacturer:** LifePro

WARRANTY AND SUPPORT

The LifePro LumiCure Pro Light Therapy Torch comes with a **Lifetime Warranty**, reflecting our commitment to quality and customer satisfaction. For any questions, technical support, or warranty claims, please contact LifePro customer service or visit the official LifePro store.

LifePro Official Store: [Visit the LifePro Store on Amazon](#)

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How should I use the LifePro LumiCure Pro Light Therapy Torch to help improve my sleep?

While the LifePro LumiCure Pro manual does not provide a specific protocol exclusively for sleep, red light therapy is often used to support the body's natural rhythms. To use the device effectively, follow these...

[Read full answer](#)