

[Manuals.plus](#) /

› [Oklar](#) /

› Oklar Blood Pressure Monitor User Manual

Oklar W1101L

Learn how to use your Oklar W1101L Blood Pressure Monitor with this quick start guide. Includes steps for charging, wearing the cuff, taking measurements, and understanding readings, plus dual-user mode info.

Oklar Blood Pressure Monitor User Manual

Brand: Oklar

Model: W1101L

INTRODUCTION

Thank you for choosing the Oklar Blood Pressure Monitor. This device is designed for convenient and accurate blood pressure and heart rate measurement at home. Its user-friendly features, including a large LED display and voice broadcast, make it accessible for all users. Please read this manual thoroughly before use to ensure proper operation and accurate readings.



Image: Oklar Wrist Blood Pressure Monitor with LED display showing readings.

SETUP

Unboxing and Initial Setup

Your Oklar Blood Pressure Monitor package includes:

- Blood Pressure Monitor
- Quick Guide
- Carrying Box
- Power Cord (Type-C cable)
- Instruction Manual

LED DISPLAY

Clearly visible at the night time

VOICE BOARDCAST

Read your result correct



Image: Oklar Blood Pressure Monitor and its accessories including carrying case and cable.

Charging the Device

The monitor is equipped with a built-in rechargeable battery. Use the provided Type-C cable to charge the device. A charging adapter is not included. Ensure the device is fully charged before first use for optimal performance.



Image: Oklar Blood Pressure Monitor connected to a Type-C charging cable.

OPERATING INSTRUCTIONS

Wearing the Cuff

To ensure accurate measurements, proper cuff placement is crucial:

1. Put the wristband on your left wrist.
2. Ensure that the distance from the wristband to the bottom of the palm is approximately 1 cm.
3. The adjustable cuff fits wrists ranging from 4.92 to 8.46 inches, providing a comfortable and secure fit.



DUAL USERS SUPPORTED



120 READING
RECORDS



120 READING
RECORDS

switch the user data by the **MEM** button



Image: Person wearing the Oklar Blood Pressure Monitor cuff on their wrist.

Taking a Measurement

Follow these steps for a fast and precise reading:

1. Sit up straight and relax.
2. Keep your heart and the machine level.
3. Press the "START/STOP" button. The measurement process takes approximately 35 seconds.
4. The device will automatically turn off after 60 seconds of inactivity.

EASY TO USE

- Wear the cuff on the left wrist
- Keep the monitor at the heart level
- Press the **START** button



Image: Hand pressing the START button on the Oklar Blood Pressure Monitor.

Understanding Readings

The large LED backlit display ensures clear readability. The device also features a voice broadcast function, ideal for individuals with poor eyesight. You can adjust the volume or turn off the sound using the SET and MEM buttons.

- **SYS (Systolic Pressure):** The top number, indicating the pressure in your arteries when your heart beats.
- **DIA (Diastolic Pressure):** The bottom number, indicating the pressure in your arteries between beats.
- **PULSE:** Your heart rate in beats per minute.



Image: Close-up of the Oklar Blood Pressure Monitor's LED display showing blood pressure and pulse readings.

Dual User Mode and Memory

The monitor supports a convenient dual user mode, allowing two users to store their readings separately. It can store up to 240 readings (120 for each user), enabling easy tracking and comparison of measurements over time.



DUAL USERS SUPPORTED



120 READING
RECORDS



120 READING
RECORDS

switch the user data by the **MEM** button



Image: Oklar Blood Pressure Monitor showing dual user icons on its display.

MAINTENANCE

Cleaning and Storage

To maintain your device, wipe it with a soft, dry cloth. Avoid using abrasive cleaners or submerging the device in water. Store the monitor in its included carrying case when not in use to protect it from dust and damage.



Image: Oklar Blood Pressure Monitor stored neatly in its compact carrying case.

TROUBLESHOOTING

Common Issues and Solutions

If you encounter any issues with your device, please refer to the following general guidelines:

- **Inaccurate Readings:** Ensure the cuff is correctly placed, you are relaxed, and the device is at heart level. Avoid talking or moving during measurement.
- **Device Not Turning On:** Check if the battery is charged. Connect the Type-C cable to recharge.
- **Display Issues:** If the display is unclear or shows error codes, try restarting the device. If the problem persists, contact customer support.

SPECIFICATIONS

Product Specifications

Feature	Detail
---------	--------

Product Dimensions	3.03 x 2.55 x 2.63 inches
Item Model Number	W1101L
Batteries	1 Lithium Ion battery (included)
Item Weight	3.59 Ounces
Display Type	LED
Use For	Wrist
Cuff Range	4.92 to 8.46 inches

Important Safety Information

Safety Guidelines

Keep out of reach of children. Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

Warranty and Support

Customer Support

For any questions, technical support, or warranty inquiries, please visit the official Oklar Store on Amazon or contact their customer service directly. Refer to the product packaging or the Quick Guide for specific contact details.

Visit the [Oklar Store](#) for more information.

Ask a Question About This Manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How do I read the charging indicator on the Oklar W1101L Blood Pressure Monitor?

To charge the Oklar W1101L Blood Pressure Monitor, connect the provided Type-C cable to the device and a power source. While the manual does not explicitly detail specific light patterns for the charging indicator,...

[Read full answer](#)

What does the "LO" reading on the screen of the Oklar W1101L Blood Pressure Monitor mean if the battery is already fully charged?

On the Oklar W1101L Blood Pressure Monitor, a "LO" reading typically indicates one of two things: Low Battery: Although you mentioned the battery is fully charged, this is the most common reason for this indicator....

[Read full answer](#)

How do I reset the Oklar W1101L Blood Pressure Monitor and fix an Error 2 (ER2) code?

While the Oklar W1101L manual does not provide a specific button combination for a factory reset, you can effectively reset the device by allowing it to power down or by ensuring it is fully...

[Read full answer](#)

What does the ER2 error code mean on the Oklar W1101L Blood Pressure Monitor, how do I fix it, and what is the warranty policy?

While the manual for the Oklar W1101L does not explicitly list a code-by-code breakdown, an "ER2" (Error 2) on most digital blood pressure monitors typically indicates a signal error or measurement failure caused by...

[Read full answer](#)

How do I access the stored memory on the Oklar W1101L Blood Pressure Monitor? The MEM button only shows the last reading.

To access the full memory on your Oklar W1101L, follow these steps: Enter Memory Mode: While the device is off, press the MEM button. The display will typically show the average of the last three...

[Read full answer](#)