

VogueNow DermGlo

VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit

USER MANUAL - MODEL: DERMGLO

Introduction

Thank you for choosing the VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit. This kit is designed to provide a comprehensive skincare experience, combining advanced LED light therapy with complementary tools for cleansing, exfoliation, and nourishment. Please read this manual thoroughly before use to ensure proper operation and to maximize the benefits of your new skincare system.

Your browser does not support the video tag.

Video Description: This video provides an overview of the VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit, showcasing the LED mask and its three complementary products: the ultrasonic facial brush, ultrasonic scrubber, and green tea serum. It highlights the comprehensive spa experience offered by the kit.

Safety Information

- Read all instructions before using the device.
- Do not use if you are pregnant, have epilepsy, or are taking photosensitizing medications. Consult a physician if you have any medical conditions.
- Avoid direct eye exposure to the LED lights. Protective eyewear is recommended.
- Keep out of reach of children.
- Do not use on broken, irritated, or infected skin.
- Ensure the device is fully charged before use and disconnect from power when not in use.
- For external use only.

Package Contents

Your VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit includes the following items:

- VogueNow DermGlo LED Light Therapy Mask
- Ultrasonic Facial Cleansing Brush
- Ultrasonic Facial Skin Scrubber
- Glean Glow Rejuvenating Green Tea Serum
- Power Adapter and USB Cable for LED Mask
- User Guide

Visible Results

in 6 – 12 Weeks

- **Smoother**
Wrinkles
- **Tighter**
Skin
- **Clearer**
Skin
- **Enhanced**
Glow



Wrinkles



Redness



Dull Skin

Image Description: A flat lay image displaying all components of the VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit. The pink LED mask is prominently featured alongside the pink ultrasonic facial scrubber, the pink ultrasonic facial cleansing brush, and a small bottle of Glean Glow Rejuvenating Green Tea Serum. The items are arranged neatly, highlighting the complete set.

Product Overview

The DermGlo kit is designed for a comprehensive at-home spa experience, targeting various skin concerns with its multi-functional devices.

LED Light Therapy Mask

The core of the kit, this mask features 150 LEDs and 7 distinct light wavelengths to address different skin needs. It is designed with an eye mold for comfort during use.



Ultimate 4-Step Spa-Like Regime for 10x Better Results

Image Description: An infographic highlighting the performance features of the DermGlo LED Face Mask. Icons and text indicate that the mask is CE-cleared for safety, utilizes non-invasive technology, offers eye-mold comfort, is derm approved, and provides 7 effective LED colors.

Ultrasonic Facial Cleansing Brush

A soft silicone brush that uses ultrasonic vibrations for deep cleansing, removing makeup and impurities from the skin's surface.

7 Wavelengths

Targeted Results



Red Light

Reduces Wrinkles



Green Light

Brightens Skin



Blue Light

Soothes Skin



Purple Light

Calms Skin



Cyan Light

Reduces Skin Tension



Yellow Light

Detoxifies Skin



White Light

Rejuvenates Skin



Image Description: A close-up shot of the pink ultrasonic facial cleansing brush, showing its soft silicone bristles and ergonomic design. The VogueNow logo is visible on the handle.

Ultrasonic Facial Skin Scrubber

This device provides deeper exfoliation, helping to clear pores, blackheads, and whiteheads, which enhances the penetration of LED light and serums.



Dr. Merrick D. Elias
FAAD - Dermatologist

Dr. Elena Dinkollari
Dermatologist

Dr. Anna Chacon
FAAD - Dermatologist

Loved by Skincare Professionals

Image Description: A close-up shot of the pink ultrasonic facial skin scrubber, highlighting its metal spatula tip and sleek design. The VogueNow logo is visible on the body of the device.

Gleam Glow Rejuvenating Green Tea Serum

A specially formulated serum designed to be used with the LED mask to boost skin-rejuvenating results, providing hydration and brightening effects.

Engineered for Performance

Designed for Daily Use



CE-cleared for
Safety



Non-invasive
Tech



Eye-mold
Comfort



Derm
Approved



7 Effective LED
Colors



Image Description: A close-up image of the Gleam Glow Rejuvenating Green Tea Serum bottle, showing its pink packaging and dropper applicator. The label indicates its rejuvenating properties.

Setup

1. **Charge Devices:** Before first use, fully charge the LED Face Mask, Ultrasonic Facial Cleansing Brush, and Ultrasonic Facial Skin Scrubber using their respective charging cables. The LED mask connects to a control unit, which then connects to a power source.
2. **Prepare Skin:** Ensure your face is clean and free of makeup before beginning any treatment.

Operating Instructions: 4-Step Spa-Grade Routine

Follow this recommended 4-step routine for optimal results:



Image Description: A collage of four images demonstrating the 4-step spa-grade light therapy routine. Each image shows a woman performing one step: cleansing with the brush, exfoliating with the scrubber, applying serum, and using the LED mask for rejuvenation.

1. Step 1: Cleanse with Ultrasonic Facial Cleansing Brush

Apply your favorite cleanser to your face. Wet the silicone brush and turn it on. Gently move the brush in circular motions across your face for 1-2 minutes to deeply cleanse and remove surface impurities. Rinse your face and the brush thoroughly.

2. Step 2: Exfoliate with Ultrasonic Facial Skin Scrubber

Ensure your skin is damp. Turn on the ultrasonic scrubber and gently glide the metal spatula across your skin, particularly on areas with blackheads or clogged pores. The vibrations will help to dislodge impurities and exfoliate dead skin cells. Rinse your face and the scrubber.

3. Step 3: Nourish with Gleam Glow Rejuvenating Green Tea Serum

Apply a few drops of the Gleam Glow Green Tea Serum to your clean, dry face. Gently massage it into your skin until absorbed. This serum is formulated to enhance the effects of the LED light therapy.

4. Step 4: Rejuvenate with LED Light Therapy Mask

Connect the LED mask to its control unit and power source. Place the mask comfortably on your face, ensuring the eye mold aligns. Select your desired LED light color and intensity using the control unit. The recommended session duration is 15-30 minutes. It is advisable to use the mask while lying down for best comfort due to its weight. After the session, turn off the mask and gently remove it.

LED Light Color Functions

The DermGlo LED Face Mask offers 7 distinct light colors, each targeting specific skin concerns:



7 Wavelengths for Targeted Results

	Red Light Reduces Wrinkles		Cyan Light Reduces Skin Tension
	Green Light Brightens Skin		Yellow Light Detoxifies Skin
	Blue Light Soothes Skin		White Light Rejuvenates Skin
	Purple Light Calms Skin		

Image Description: An infographic detailing the 7 wavelengths of LED light therapy provided by the mask. Each color (Red, Green, Blue, Purple, Cyan, Yellow, White) is listed with its corresponding skin benefit, such as reducing wrinkles, brightening skin, soothing skin, calming skin, reducing skin tension, detoxifying skin, and rejuvenating skin.

- **Red Light:** Helps in making the skin firm and reduces the appearance of wrinkles.

- **Green Light:** Brightens skin and helps to reduce hyperpigmentation.
- **Blue Light:** Soothes irritated skin and targets acne-causing bacteria.
- **Purple Light:** Calms skin and promotes cell regeneration.
- **Cyan Light:** Reduces skin tension and inflammation.
- **Yellow Light:** Detoxifies skin and improves lymphatic drainage.
- **White Light:** Offers a complete 360° skin rejuvenation, penetrating deeply to tighten and improve skin vitality.

Maintenance

- **Cleaning:** After each use, wipe the LED mask, ultrasonic brush, and scrubber with a soft, damp cloth. Do not submerge the electronic devices in water.
- **Storage:** Store all components in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Serum Storage:** Keep the Green Tea Serum in a cool, dark place.

Troubleshooting

Problem	Possible Cause	Solution
Device not turning on	Low battery or not connected to power.	Ensure devices are fully charged or properly connected to a power source.
LED lights not working	Loose connection or incorrect setting.	Check all cable connections. Ensure the correct light mode and intensity are selected on the control unit.
Skin irritation	Overuse or sensitivity.	Reduce frequency or duration of use. Discontinue use if irritation persists and consult a dermatologist.
Scrubber/Brush not vibrating	Low battery or device malfunction.	Recharge the device. If the issue persists, contact customer support.

Specifications

- **Model Number:** DermGlo (VN-DermGlo)
- **UPC:** 197644922376
- **Package Dimensions:** 10.47 x 8.62 x 5.2 inches
- **Weight:** 2.69 Pounds
- **LED Mask Features:** 150 LEDs, 7 Wavelengths (Red, Green, Blue, Purple, Cyan, Yellow, White)
- **Included Components:** LED Face Mask, Ultrasonic Facial Cleansing Brush, Ultrasonic Facial Skin Scrubber, Green Tea Serum
- **Skin Type Compatibility:** All
- **Product Benefits:** Skin Brightening, Skin Rejuvenation, Wrinkle Reduction

Warranty and Support

Your VogueNow DermGlo LED Face Mask Light Therapy Premium Spa Kit comes with a **1-year warranty** from the date of purchase, in addition to Amazon's 30-day guarantee. This warranty covers manufacturing defects and malfunctions under normal use.

If you encounter any issues with your product or require assistance, please contact VogueNow customer support. Our team is dedicated to ensuring your satisfaction and will assist you promptly.

Contact Information: Please refer to the contact details provided on the product packaging or the official VogueNow website for the most up-to-date support information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question