

Phonak Power Smokey Dome Small

Phonak Power Smokey Dome Small (10mm) Hearing Aid Domes Instruction Manual

INTRODUCTION

This instruction manual provides essential information for the proper use, installation, and maintenance of your Phonak Power Smokey Dome Small (10mm) hearing aid domes. These genuine OEM replacement domes are manufactured by Sonova in Switzerland, ensuring quality and compatibility with specified Phonak hearing aid models.

PRODUCT OVERVIEW

Phonak Power Smokey Domes are designed to provide a secure and comfortable fit for your hearing aid receiver, optimizing sound delivery. This package contains 10 domes, each measuring 10mm (0.39 inches) in diameter. They are made from a durable, smokey-colored silicone material.



Image: Phonak Power Smokey Dome Small (10mm) domes shown in their transparent packaging, indicating a pack of 10 units.

Per pack Included 10 Domes Total

PHONAK



**GENUINE OEM
PHONAK DOMES**

Made in Switzerland

Image: A close-up of the Phonak Power Smokey Dome packaging, highlighting the 'Genuine OEM Phonak Domes' and 'Made in Switzerland' labels.

UNDERSTANDING DOME TYPES AND SIZES

Phonak offers various dome types and sizes to suit individual hearing needs and ear canal anatomies. The Power Smokey Domes are designed for specific amplification requirements.

Dome Sizes

Power Smokey Domes are available in Small (S), Medium (M), and Large (L) sizes. This product specifically provides the Small (10mm) size.



Power Smokey dome S

054-1993



Power Smokey dome M

054-1994



Power Smokey dome L

054-1995

PHONAK

Image: Visual comparison of Power Smokey Dome sizes Small (S), Medium (M), and Large (L), each with its corresponding reference number.

Detailed Size Chart for Various Dome Types

Refer to the chart below for a comprehensive overview of different Phonak dome types and their respective dimensions. This helps in identifying the correct dome for your specific hearing aid and receiver.

PHONAK

Made in Switzerland



Standard



Power

Phonak Standard and Power (2 Pin)xReceivers



Image: A detailed chart illustrating the sizes (in mm and inches) for Open Smokey, Closed Smokey, Power Smokey, Open 4.0, Vented 4.0, Power 4.0, and Cap 4.0 domes, categorized for Old and New Hearing Aids.

COMPATIBILITY


These Phonak Power Smokey Domes are compatible with Phonak Power (2 Pin) xReceivers and a range of Phonak Audeo and Bolero hearing aid models manufactured from 2016 onwards. It is crucial to ensure your hearing aid model is listed below or uses a compatible receiver type.

Compatible Receivers

The domes are designed for use with Phonak Power (2 Pin) xReceivers. Ensure your receiver wire matches the 2-pin configuration.

PhonakPower (2 Pin)xRecievers

Phonak Power RIC xReceiver wire (2 Pin)
compatible with smokey power domes



Genuine OEM PHONAK DOMES

Image: Illustration showing various Phonak Power (2 Pin) xReceivers (0xP, 1xP, 2xP, 3xP) and their compatibility with Power Smokey Domes in sizes S, M, and L.

Audiologists' Recommend

Authentic Phonak Hearing aid Domes



Image: Depicts Phonak Standard (0xS) and Power (2xP) 2-pin xReceivers alongside examples of Closed Smokey, Open Smokey, and Power Smokey Domes.

Compatible Hearing Aid Models

These domes are compatible with the following Phonak Audeo and Bolero hearing aid models (manufactured from 2016 onwards) that use xReceiver Standard, xReceiver Power, SuperPower plus xReceiver Old Hearing Aids:

- Phonak Audeo B90-10, Phonak Audeo B90-312, Phonak Audeo B90-312T, Phonak Audeo B90-13, Phonak Audeo B70-10, Phonak Audeo B70-312, Phonak Audeo B70-312T, Phonak Audeo B70-13, Phonak Audeo B50-10, Phonak Audeo B50-312, Phonak Audeo B50-312T, Phonak Audeo B50-13
- Phonak Bolero B90-M, Phonak Bolero B90-P, Phonak Bolero B90-SP, Phonak Bolero B70-M, Phonak Bolero B70-P, Phonak Bolero B70-SP, Phonak Bolero B50-M, Phonak Bolero B50-P, Phonak Bolero B50-SP, Phonak Bolero B30-M, Phonak Bolero B30-P, Phonak Bolero B30-SP

PERFECT FIT - OEM Phonak Smokey Open, Closed and Power Domes are Compatible with following Hearing Aids

PLEASE DO BUY THIS SET FOR THE FOLLOWING HEARING AIDS ARE COMPATIBLE
COMPATIBLE WITH USE xReceiver Standard,xReceiver power, SuperPower plus
xReceiver Old Hearing Aids manufactured in 2016-

COMPATIBLE WITH Phonak Audeo B90-10, Phonak Audeo B90-312, Phonak Audeo B90-312T, Phonak Audeo B90-13, Phonak Audeo B70-10, Phonak Audeo B70-312, Phonak Audeo B70-312T, Phonak Audeo B70-13, Phonak Audeo B50-10, Phonak Audeo Phonak Audeo B90-10, Phonak Audeo B90-312, Phonak Audeo B90-312T, Phonak Audeo B90-13, Phonak Audeo B70-10, Phonak Audeo B70-312, Phonak Audeo B70-312T, Phonak Audeo B70-13 Phonak Audeo B50-10, Phonak Audeo B50-312, Phonak Audeo B50-312T, Phonak Audeo B50-13, Phonak Audeo B30-10, Phonak Audeo B30-312, Phonak Audeo B30-312T, Phonak Audeo B30-13 and Phonak Bolero B90-M, Phonak Bolero B90-P, Phonak Bolero B90-SP, Phonak Bolero B70-M, Phonak Bolero B70-P, Phonak Bolero B70-SP, Phonak Bolero B50-M, Phonak Bolero B50-P, Phonak Bolero B50-SP, Phonak Bolero B30-M, Phonak Bolero B30-P, Phonak Bolero B30-SP



Image: A graphic detailing the compatibility of Phonak Smokey Open, Closed, and Power Domes with various Phonak Audeo and Bolero hearing aid models.

INSTALLATION INSTRUCTIONS

Follow these steps to correctly attach a new Phonak Power Smokey Dome to your hearing aid receiver. Always handle domes with clean hands.

1. Gently remove the old dome from the receiver by pulling it straight off.
2. Ensure the receiver tip is clean and free from any debris.
3. Take a new Power Smokey Dome (Small, 10mm) from its packaging.
4. Align the opening of the dome with the tip of the receiver.
5. Push the dome firmly onto the receiver until it clicks securely into place. Ensure it is fully seated and not loose.
6. Verify that the dome is securely attached before inserting the hearing aid into your ear.

Video Guide: Attaching Hearing Aid Domes

Your browser does not support the video tag.

Video: This video demonstrates the different types of Phonak hearing aid domes (Open, Vented, Power) and provides a visual guide on how to correctly attach them to the hearing aid receiver. It shows the process for each dome type, ensuring a secure fit.

MAINTENANCE AND CARE

Regular cleaning and timely replacement of your hearing aid domes are crucial for optimal performance and hygiene.

- **Daily Cleaning:** Gently wipe the dome with a clean, dry cloth after each use to remove earwax and debris.
- **Avoid Liquids:** Do not use water or cleaning solutions directly on the domes unless specifically instructed by your audiologist or the manufacturer.
- **Regular Replacement:** Domes should be replaced regularly, typically every 2-4 weeks, or sooner if they appear discolored, stiff, or damaged. This prevents wax buildup and maintains sound quality.
- **Storage:** Store unused domes in a cool, dry place away from direct sunlight.

IMPORTANT INFORMATION

Hygiene Notice

Due to hygiene reasons, Phonak Domes are not refundable or exchangeable after opening the package. Only unopened packs can be returned. This policy is in place to ensure the health and safety of all users.

Made in Switzerland

PHONAK

GENUINE OEM
PHONAK DOMES



Image: Packaging of Phonak Power Smokey Domes displaying a prominent notice regarding the non-returnable policy for opened packages due to hygiene reasons.

Legal Disclaimer






Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a medical accessory and should be used as directed by a healthcare professional or audiologist.

SPECIFICATIONS

Feature	Detail
Product Name	Phonak Power Smokey Dome Small
Size	Small (10mm / 0.39 inches)
Quantity	10 Domes per pack
Manufacturer	Sonova
Origin	Switzerland

Product Dimensions	1.93 x 1.73 x 0.75 inches (packaging)
Weight	0.63 ounces (packaging)
ASIN	B0CK5SBMM2

Related Documents - Power Smokey Dome Small

	<p>Phonak RIC Hearing Aid Maintenance Guide for Marvel, Paradise, and Lumity</p> <p>A comprehensive guide from Phonak on how to clean, maintain, and troubleshoot your Marvel, Paradise, and Lumity RIC hearing aids, including filter replacement and moisture management. Learn about ear domes, filter types, drying procedures, and service options.</p>
	<p>Phonak RIC Coupling and Custom Earpiece Guide</p> <p>A comprehensive guide to Phonak RIC coupling options, including Infinio and Lumity receivers, domes, custom ear pieces, and CROS systems. Also covers custom product configurations for Phonak Virto models with shell and faceplate color options, fitting ranges, and wax protection.</p>
	<p>Phonak RemoteControl User Guide - Enhance Your Hearing Experience</p> <p>Comprehensive user guide for the Phonak RemoteControl, detailing setup, operation, troubleshooting, safety, and warranty information for managing hearing aid programs and volume.</p>
	<p>Phonak RightFit™ Custom Hearing Aid Ordering Guide</p> <p>A comprehensive guide for ordering Phonak RightFit™ custom hearing aids, detailing steps for 'as-small-as-possible' (HS) and 'dexterity' (FS) models. Learn about shell colors, faceplate options, and the benefits of precision fit and optimized sound quality for patients.</p>
	<p>Phonak Brio 3 R-C Hearing Aid - Features, Benefits, and Fitting</p> <p>Explore the Phonak Brio 3 R-C hearing aid, detailing its direct connectivity, hands-free calling, advanced features like HiFrequency Protect+ and AutoSense OS, and fitting ranges. This document provides an overview of its benefits and fitting updates.</p>

Phonak Charger Case RIC

Quick Guide | Kurzanleitung | Guide rapide | Guía rápida



Please refer to the User Guide for further instructions.
Bitte beachten Sie die Gebrauchsanweisung für weitere Anweisungen.
Veuillez vous référer au mode d'emploi pour plus d'informations.
Por favor consulte las instrucciones de uso para más información.

PHONAK
RIC RIC

[Phonak Charger Case RIC and Power Pack Quick Guide](#)

Quick guide for Phonak Charger Case RIC and Power Pack, detailing charging procedures and battery indicators. Includes instructions in English, German, French, and Spanish.