

ELENKER 3 Function Full Electric Hospital Bed

ELENKER 3 Function Full Electric Hospital Bed User Manual

Model: 3 Function Full Electric Hospital Bed

INTRODUCTION

This manual provides detailed instructions for the safe and effective use, assembly, and maintenance of your ELENKER 3 Function Full Electric Hospital Bed. Please read this manual thoroughly before operating the bed to ensure proper function and user safety. Keep this manual for future reference.

SAFETY INFORMATION

Important Safety Precautions:

- Always ensure the bed is on a stable, level surface.
- The maximum safe operating weight limit for this bed is **530 lbs**. Do not exceed this weight.
- Ensure all electrical connections are secure and free from damage. Do not operate the bed with damaged cords or plugs.
- Keep hands and feet clear of moving parts during adjustment.
- Lock the casters when the bed is in a stationary position to prevent unintended movement.
- Regularly inspect the bed for any signs of wear or damage. Discontinue use and contact support if issues are found.
- This bed is intended for medical and home care use. Consult a healthcare professional for specific patient positioning needs.

PRODUCT OVERVIEW

The ELENKER 3 Function Full Electric Hospital Bed is designed to provide adjustable positioning for enhanced comfort and care. Key components and features include:

- **Full Electric Adjustment:** Control bed height, backrest, and legrest positions electronically.
- **Remote Control:** Handheld remote for easy adjustment of bed functions.
- **High-Density Foam Mattress:** A 4.7-inch foam mattress with a removable and washable cover.

- **IV Pole:** Integrated adjustable IV pole for medical support.
- **Dining Table:** Convenient attachable dining table for meals or activities.
- **Foldable Guardrails:** Side guardrails for patient safety, which can be folded down.
- **Lockable Casters:** 5-inch casters for mobility, with locking mechanisms for stability.



Image: ELENKER 3 Function Full Electric Hospital Bed with a patient utilizing the dining table.



1
Handheld
Controller



2
Wire Storage
Hooks



3
4 Drainage
Bag Hooks



4
Dining Table
Board



5
Patient Information
Board



6
Corner
Bumper



Image: Overview of the ELENKER 3 Function Full Electric Hospital Bed, showcasing its premium features.




Image: Close-up view of the adjustable IV pole, designed for convenient medical support.



Image: Close-up view of the attachable bed tray, providing a surface for dining or activities.



Image: Close-up view of the foldable guardrail, showing its safety and convenience features.



Image: Close-up view of the 5-inch lockable casters, ensuring secure positioning and easy mobility.

SETUP

The ELENKER 3 Function Full Electric Hospital Bed requires minimal assembly. You will primarily need to install the headboards, footboards, and guardrails. The main bed frame comes largely pre-assembled.

1. **Unpacking:** Carefully remove all components from the packaging. It is recommended to have assistance due to the weight of the bed.
2. **Frame Placement:** Position the main bed frame in the desired location. Ensure adequate space around the bed for operation and access.
3. **Caster Installation (if separate):** Attach the 5-inch casters to the bed frame if they are not pre-installed. Ensure they are securely fastened.
4. **Headboard and Footboard Installation:** Attach the headboard and footboard to the designated slots on the bed frame. Secure them using the provided hardware.
5. **Guardrail Installation:** Install the foldable guardrails on both sides of the bed. Ensure they click securely into place.
6. **Mattress Placement:** Place the high-density foam mattress onto the bed frame. Ensure it is properly aligned.
7. **IV Pole and Dining Table:** Attach the IV pole and dining table to their respective mounting points as needed.
8. **Power Connection:** Connect the power cord to a suitable electrical outlet.

For a visual guide, please refer to the installation video below:

Video: Installation guide for the ELENKER 3 Function Hospital Bed. This video demonstrates the step-by-step assembly process for the bed components.


OPERATING INSTRUCTIONS

The ELENKER 3 Function Full Electric Hospital Bed is operated using a handheld remote control. The remote allows for independent adjustment of three key functions:

- **Backrest Adjustment:** Adjusts the angle of the backrest from 0° to 80°.
- **Legrest Adjustment:** Adjusts the angle of the leg position from 0° to 45°.
- **Full Bed Elevation:** Adjusts the overall height of the bed from 16.3 inches to 28.3 inches.

Using the Remote Control:

Hand controller



Large Buttons

Large buttons are easy to identify and operate

Image: Handheld remote controller with large, easy-to-identify buttons for various bed adjustments.

Press the corresponding buttons on the remote control to achieve the desired position. The remote typically features intuitive icons for each function (e.g., an upward arrow for raising, a downward arrow for lowering).

Backrest and Legrest Adjustment



Image: Visual representation of the backrest (0°-80°) and legrest (0°-45°) adjustment capabilities.

Height Adjustment

Height Range: 16.3-28.3"




Image: Diagram illustrating the bed's height adjustment range from 16.3 inches to 28.3 inches.

For a detailed demonstration of the bed's functions, please watch the video below:

Video: Detailed display and application of the ELENKER 3 Function Hospital Bed, demonstrating its various adjustable features.

MAINTENANCE

Proper maintenance ensures the longevity and safe operation of your hospital bed.

- **Cleaning the Bed Frame:** Wipe down the bed frame with a damp cloth and mild detergent. Avoid abrasive cleaners or solvents.
- **Mattress Care:** The mattress features a removable and washable cover. Follow the care instructions on the mattress label for cleaning. Ensure the mattress is completely dry before re-covering and use.
- **Caster Maintenance:** Periodically check the casters for debris and ensure they roll smoothly. Lubricate if necessary.
- **Electrical Components:** Inspect all electrical cords and connections regularly for fraying, cuts, or loose

connections. Do not attempt to repair damaged electrical components yourself.

- **Moving Parts:** Ensure all moving parts operate smoothly and are free from obstruction.

High Density Foam Mattress



Removeable & Washable
Mattress Cover



High elasticity
foam



Skin friendly and
comfortable



Healthy and
odorless

Image: High-density foam mattress, highlighting its removable and washable cover for easy maintenance.

TROUBLESHOOTING

If you encounter issues with your ELENKER 3 Function Full Electric Hospital Bed, refer to the following common troubleshooting steps:

- **Bed Not Responding to Remote:**

- Check if the bed is properly plugged into a working electrical outlet.
- Ensure the power switch on the bed (if present) is in the 'ON' position.
- Inspect the remote control cable for any damage.

- **Unusual Noises During Operation:**

- Check for any obstructions under or around the bed that might interfere with moving parts.
- Ensure all bolts and fasteners are tight.

- **Bed Not Moving Smoothly:**

- Verify that the casters are unlocked if you intend to move the bed.
- Check for any debris or obstacles hindering movement.

If these steps do not resolve the issue, please contact ELENKER customer support for further assistance.

SPECIFICATIONS

Feature	Detail
Model	3 Function Full Electric Hospital Bed
Overall Bed Dimensions	85.2 x 36.8 inches
Adjustable Height Range	16.3 - 28.3 inches
Backrest Adjustment	0° - 80°
Leg Position Adjustment	0° - 45°
Safe Operating Weight Limit	530 lbs
Item Weight	258 lbs
Mattress Thickness	4.7 inches
Caster Size	5 inches
Package Dimensions	86 x 40 x 18.5 inches

Dimensions and Weight Capacity




Image: Detailed diagram illustrating the dimensions and weight capacity of the hospital bed.

WARRANTY AND SUPPORT

For warranty information, product support, or to report any issues, please contact ELENKER customer service. Refer to your purchase documentation for specific warranty terms and contact details.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question