



Manuals.plus /

› DigiKuber /

› DigiKuber ECG Smart Watch TK10 User Manual

DigiKuber TK10

DigiKuber ECG Smart Watch TK10 User Manual

Model: TK10

INTRODUCTION

This manual provides comprehensive instructions for the setup, operation, and maintenance of your DigiKuber ECG Smart Watch TK10. Please read this manual carefully to ensure proper use and to maximize the functionality of your device. This smartwatch is designed to assist in monitoring various health parameters and activity levels.

PACKAGE CONTENTS

- DigiKuber ECG Smart Watch TK10
- USB Charging Cable
- User Manual (this document)

SETUP

1. Charging the Smart Watch

Before first use, fully charge your smartwatch. Connect the magnetic charging cable to the charging points on the back of the watch and plug the USB end into a standard USB power adapter (not included) or a computer USB port. A full charge typically takes approximately 3 to 4 hours.

2. Installing the Deepfit App

The DigiKuber Smart Watch requires the 'Deepfit' application for full functionality and data synchronization. You can download the Deepfit app from the Google Play Store (for Android devices) or the Apple App Store (for iOS devices).

1. Open your smartphone's app store.
2. Search for 'Deepfit'.
3. Download and install the app.

3. Pairing with Your Smartphone

Ensure Bluetooth is enabled on your smartphone and the watch is sufficiently charged.

1. Open the Deepfit app on your smartphone.
2. Follow the on-screen instructions to add a new device.
3. Select 'TK10' from the list of available devices.
4. Confirm the pairing request on both your watch and smartphone.
5. Once paired, the watch's date and time will automatically synchronize with your phone.

OPERATING INSTRUCTIONS

Basic Navigation

- **Touch Screen:** Swipe left/right, up/down to navigate through menus and functions.
- **Side Button:** Press to return to the home screen or wake the display.
- **Wrist Gesture:** Turn your wrist to illuminate the screen.

Health Monitoring Features

The watch offers various health monitoring capabilities. For accurate readings, ensure the watch is worn snugly on your wrist.



The DigiKuber TK10 Smart Watch interface showing various health metrics and its companion app displaying ECG details.

- **ECG Monitor:** Initiate an ECG measurement from the watch menu. The data, including AI diagnosis, can be viewed and exported from the Deepfit app.
- **Heart Rate (PPG):** Continuously monitors your heart rate.
- **Blood Pressure:** Provides estimated blood pressure readings.
- **Blood Oxygen (SpO2):** Measures your blood oxygen saturation levels.
- **Body Temperature:** Measures your body temperature for 24 hours, with data saved in the Deepfit app.
- **Sleep Monitor:** Automatically tracks your sleep quality and patterns when worn during sleep. Data syncs with the Deepfit app.
- **Stress Monitor:** Provides stress level assessments.

Note: The health data provided by this device is for reference only and should not be used for medical diagnosis or treatment. Consult a medical professional for any health concerns.

Activity Tracking

- **Pedometer:** Tracks your steps taken.
- **Calorie Counter:** Estimates calories burned based on activity.

- **Distance View:** Calculates distance covered.
- **Sports Modes:** Includes various modes such as running, cycling, rope jumping, basketball, table tennis, badminton, and general exercise. Select the appropriate mode before starting an activity to record specific data.

Smart Features

- **Call and Message Notifications:** Receive alerts for incoming calls and messages (SMS, social apps) directly on your watch. You can view message content on the watch.
- **Music Control:** Control music playback on your connected smartphone (play/pause, next/previous track).
- **Camera Control:** Remotely trigger your smartphone's camera shutter.
- **Stopwatch & Countdown:** Utilize built-in timing functions.
- **Sedentary Reminder:** Alerts you when you have been inactive for too long.
- **Find My Phone:** Helps locate your paired smartphone by making it ring.
- **Remote Care Function:** Share real-time health data with family members via the Deepfit app.

MAINTENANCE

Water Resistance (IP68)

The DigiKuber TK10 Smart Watch has an IP68 waterproof rating, meaning it is resistant to dust and can withstand immersion in water up to 1.5 meters for 30 minutes. It is suitable for daily use, such as hand washing, rain, and light splashes.

Important Precautions:

- Do **NOT** wear the watch when diving.
- Avoid wearing the watch during hot water baths, saunas, or in environments with high temperature and humidity, as steam can penetrate the seals and cause damage.
- Do not press any buttons underwater.
- After contact with water, especially saltwater, gently wipe the watch dry with a soft cloth.

Cleaning

Regularly clean your watch and strap to prevent skin irritation and maintain device performance.

- Wipe the watch screen and body with a soft, damp, lint-free cloth.
- For the strap, follow the manufacturer's cleaning recommendations for the specific material (silicone TPU or metal).
- Do not use harsh chemicals or abrasive cleaners.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Watch does not turn on	Low battery	Charge the watch for at least 30 minutes.
Cannot pair with smartphone	Bluetooth off, app not running, watch not discoverable, previous connection issues	Ensure Bluetooth is on. Restart the Deepfit app. Restart your smartphone. Restart the watch. Forget the device in your phone's Bluetooth settings and try pairing again.
Inaccurate health readings	Improper fit, movement during measurement, environmental factors	Ensure the watch is worn snugly on your wrist, about one finger's width from the wrist bone. Remain still during measurements. Note that these readings are for reference and not medical grade.
Notifications not received	App permissions, Bluetooth disconnected, notification settings in app	Check app permissions for notifications on your phone. Ensure the watch is connected via Bluetooth. Verify notification settings within the Deepfit app.
Short battery life	Frequent use of features, high screen brightness, continuous monitoring	Reduce screen brightness. Turn off continuous monitoring for features not constantly needed. Limit frequent use of power-intensive features.

SPECIFICATIONS

Feature	Detail
Model	TK10
Display	1.91 inch Color Touch Screen, IPS TFT 240 x 286
Compatible System	Android 6.0 and above, iOS 10 and above
Bluetooth	5.0 BLE (Low Consumption)
App Name	Deepfit
Battery Type	Lithium Polymer Battery
Battery Capacity	280mAh
Charging Type	Magnetic Absorption
Charging Time	Approx. 3 to 4 hours
Usage Time	5 to 7 days
Standby Time	20 to 30 days
Water Level	IP68
Dimensions	45 x 38 x 11.8 mm (Length, Width, Height)
Material	Case: Metal alloy; Lower Shell: PC+ABS; Bracelet: Silicone TPU or Metal

WARRANTY AND SUPPORT

Warranty Information

DigiKuber offers a one-year warranty for this product from the date of purchase. This warranty covers manufacturing defects and issues arising from normal use. It does not cover damage caused by misuse, accidents, unauthorized modifications, or failure to follow the instructions in this manual.

Customer Support

For technical support, warranty claims, or any questions regarding your DigiKuber ECG Smart Watch TK10, please contact DigiKuber customer service through the contact information provided on the official DigiKuber website or through your point of purchase.