

[Manuals.plus](#) /

› [Salter](#) /

› Salter 9127 SVSV3R Bathroom Analyser Scale User Manual

Salter 9127 SVSV3R

Salter 9127 SVSV3R Bathroom Analyser Scale User Manual

Comprehensive guide for your Salter Analyser Scale, covering setup, operation, and maintenance.

INTRODUCTION

The Salter 9127 SVSV3R Bathroom Analyser Scale is designed to help you monitor your health and fitness effortlessly. This scale provides comprehensive body composition analysis, including weight, body fat, body water, muscle mass, bone mass, Body Mass Index (BMI), and Basal Metabolic Rate (BMR). With its user-friendly features and sleek design, it's an ideal tool for tracking your progress towards your health goals.



Front view of the Salter Analyser Scale, showcasing its sleek glass design and the included AAA batteries.

SETUP

1. Battery Installation

The Salter Analyser Scale comes with 2 AAA batteries. Locate the battery compartment on the underside of the scale. Insert the batteries, ensuring correct polarity (+/-). Close the battery compartment securely.

2. Initial Placement

For accurate readings, place the scale on a firm, flat surface. Avoid placing it on carpets or uneven flooring. The scale's compact size allows it to fit away nicely in most bathrooms.



An angled perspective of the Salter Analyser Scale, highlighting its modern glass platform and compact design.

3. User Profile Setup

The scale features a 10-user memory, allowing multiple individuals to track their progress. To set up a user profile, use the control buttons on the display. You can preset your age and height. Once these details are entered, the scale will use them to calculate your Body Mass Index (BMI) and other body composition metrics. The scale saves this data for future reference.



The scale's display area, indicating its capacity for 10 user profiles and a dedicated athlete mode for specialized tracking.

OPERATING INSTRUCTIONS

1. Basic Weighing

To obtain an instant weight reading, simply step onto the scale. Ensure your feet are evenly placed on the toughened glass platform. The large, easy-to-read display will show your weight.



Image showing feet on the scale, demonstrating the simple 'step on' feature for immediate weight measurement.

2. Body Composition Analysis

For a full body composition analysis, step onto the scale barefoot. The scale will measure and display various metrics including weight, body fat, body water, muscle mass, bone mass, BMI, and BMR. These measurements provide a comprehensive overview of your body's composition.



An illustrative graphic detailing the comprehensive body composition analysis provided by the scale, including weight, body fat, body water, muscle mass, bone mass, BMI, and BMR.

3. Switching Units

The scale allows you to switch between metric (kilograms) and imperial (stones/pounds) units. Refer to the display for the current unit and use the appropriate button (usually on the underside or near the display) to toggle between them.



Easy-to-read backlit display



*Switch between
metric and imperial*



A detailed view of the scale's clear, backlit LCD display, highlighting the ability to switch between kilogram (KG), stone (ST), and pound (LB) units.

MAINTENANCE

To ensure the longevity and accuracy of your Salter Analyser Scale, follow these maintenance guidelines:

- **Cleaning:** Clean the scale with a damp cloth. Do not use abrasive cleaners or chemical solvents, as these can damage the glass surface or electronic components.
- **Storage:** Store the scale in a dry place at room temperature. Avoid extreme temperatures or high humidity.
- **Battery Replacement:** Replace batteries when the low battery indicator appears on the display. Always replace all batteries at the same time with new ones of the same type.

TROUBLESHOOTING

If you encounter issues with your Salter Analyser Scale, refer to the following common troubleshooting tips:

Problem	Possible Cause	Solution
"Lo" displayed	Low battery.	Replace all batteries with new AAA batteries.
"Err" displayed	Overload or measurement error.	Ensure weight is within the 180 kg capacity. Step off and then back on the scale.
Inaccurate readings	Scale on uneven surface or not calibrated.	Place scale on a hard, flat surface. Step on and off quickly to calibrate to zero before weighing.
No body composition readings	Not standing barefoot or profile not set up.	Ensure you are standing barefoot with clean, dry feet. Verify your user profile is correctly set up with age and height.

SPECIFICATIONS

- Brand:** Salter
- Model Number:** 9127 SVSV3R
- Product Dimensions:** 11.81 x 11.81 x 0.98 inches (30 x 30 x 2.5 cm)
- Weight Capacity:** Up to 180 kg (396 lbs / 28 st 4 lbs)
- User Memory:** 10 users
- Measurements:** Weight, Body Fat, Body Water, Muscle Mass, Bone Mass, BMI, BMR
- Display Type:** LCD
- Batteries:** 2 x AAA (included)
- Manufacturer:** Ultimate Products Ltd
- Color:** Silver



Diagram showing the compact dimensions of the Salter Analyser Scale, measuring 30cm in length and width, and 2.5cm in height.

WARRANTY AND SUPPORT

For detailed warranty information and customer support, please refer to the documentation included with your product or visit the official Salter website. Keep your purchase receipt as proof of purchase for any warranty claims.

© 2024 Salter. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question

and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0D14RW8CL>

More documents and user manuals: **<https://manuals.plus/>**