

[Manuals.plus](#) /

› [PMD Beauty](#) /

› PMD Beauty Clean Redvolution Facial Cleansing Brush User Manual

PMD Beauty Clean Redvolution

PMD Beauty Clean Redvolution Facial Cleansing Brush User Manual

Model: Clean Redvolution

[Visit PMD Beauty Official Website](#)

INTRODUCTION

The PMD Clean Redvolution combines the cleansing power of the PMD Clean with the rejuvenating treatment of red light therapy. This device features SonicGlow Technology for a gentle yet deep cleansing experience and 19 powerful red LED lights designed to visibly smooth, lift, and firm the skin.

KEY FEATURES

- **SonicGlow Technology:** Delivers 7,000 vibrations per minute for deep cleansing, effectively removing dirt and oil.
- **Red Light Therapy:** Features FDA-cleared medical-grade LEDs to promote collagen production, helping to smooth fine lines and improve skin radiance.
- **Touch Sensor:** LED lights activate only when the ring makes contact with the skin, ensuring targeted treatment.
- **Ultra-Hygienic Silicone:** Odor-resistant, hypoallergenic, and waterproof for safe and easy use.
- **Waterproof Design:** IPX7 waterproof protection makes it safe for use in the sink or shower.
- **Long-lasting & Rechargeable Battery:** Provides extended use between charges.



Image: The PMD Beauty Clean Redvolution device with callouts explaining its core features: Red Light Therapy, SonicGlow Technology, Touch Sensor, and Ultra-Hygienic Silicone.

ALL-IN-ONE BEAUTY



Waterproof &
Shower-Friendly
Design



Switch between
Cleansing &
Treatment



Deep Cleansing to
remove dirt, oil,
and blackheads



Red Light
Treatment targets
fine-lines and
supports skin
repair

Image: The PMD Clean Redvolution shown in use for deep cleansing with water and for red light treatment targeting fine lines.

SETUP

1. **Unboxing:** Carefully remove the PMD Clean Redvolution device and all accessories from the packaging.
2. **Initial Charge:** Before first use, fully charge the device using the provided charging cable. The device requires 1 Lithium Ion battery (included).
3. **Clean Device:** Rinse the silicone head of the device with warm water before its first application to your skin.

OPERATING INSTRUCTIONS

Cleansing Mode

1. Rinse your face with warm water.
2. Apply your preferred cleanser directly to the soft silicone bristles of the PMD Clean Redvolution.
3. Wet the device.
4. Turn on the device to mode 1 (gentle vibration) or mode 2 (intense vibration).
5. Using circular motions, gently cleanse your entire face.
6. Rinse your face and the device thoroughly after use.
7. Use 1-2 times daily for best results.

Red Light Therapy Mode

1. After cleansing, press the power button to activate the red light therapy mode.
2. Place the red light side of the device against the areas of your skin you wish to treat.
3. The device will vibrate upon contact with the skin, indicating the start of the treatment.
4. Hold the device to each desired area to help smooth fine lines and boost skin radiance.
5. This mode is safe for daily use.

HOW TO USE IN 3 SIMPLE STEPS



1



Apply cleanser

Apply


2



Cleanse in **circular motions** and rinse

Cleanse

3



Flip the device, turn on LED, and apply to skin

Treat

Image: Visual guide showing the three simple steps to use the PMD Clean Redvolution: apply cleanser, cleanse, and treat with red light.

Your browser does not support the video tag.

Video: An official PMD Beauty video demonstrating how the device helps fight breakouts and achieve cleaner, softer skin through its cleansing and red light features.

MAINTENANCE

To maintain your PMD Clean Redvolution device, simply wash the silicone head with warm water after each use. Ensure the device is dry before storing. No replacement brush heads are needed due to the ultra-hygienic silicone material.

TROUBLESHOOTING

- **Device Not Turning On:** Ensure the device is fully charged. Connect it to the charging cable and allow sufficient

time for charging.

- **Red Light Not Activating:** The red LED lights are designed to activate only when the ring makes contact with the skin. Ensure proper contact is made during use.
- **Reduced Vibration:** If the device's vibration seems weaker than usual, it may require recharging.
- **Cleaning Issues:** For stubborn residue, use a mild soap with warm water to clean the silicone head.

SPECIFICATIONS

Brand Name	PMD Beauty
Model	Clean Redvolution
ASIN	B0D24FZMDX
UPC	811485037074
Material Type	Silicone
Power Source	1 Lithium Ion battery (included)
Operation Mode	Automatic
Product Benefits	Pore Treatment, Fine Lines Treatment, Anti-Aging, Blemishes Treatment, Cleansing
Skin Type	All
Product Care	Wash with Warm Water

WARRANTY AND SUPPORT

PMD Beauty offers an unmatched 20-year warranty on the PMD Clean Redvolution device. For warranty claims, technical support, or any product-related inquiries, please visit the official PMD Beauty website or contact their customer service directly.

For further assistance, please refer to the contact information provided with your product packaging or visit the [PMD Beauty website](#).

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question