

Vimmk Y208

Vimmk Y208 Water Flosser Replacement Nozzles Instruction Manual

Brand: Vimmk | Model: Y208

1. INTRODUCTION

This manual provides essential information for the proper use and maintenance of your Vimmk Y208 Water Flosser Replacement Nozzles. These nozzles are designed to integrate seamlessly with your Vimmk Y208 portable professional oral irrigator, supporting your daily oral hygiene routine.

2. PACKAGE CONTENTS

The Vimmk Y208 replacement nozzle set includes 8 distinct nozzles, each designed for specific oral care needs. The different colored rings on the nozzles allow for easy identification among family members.

- 4 x Standard Nozzles
- 1 x Orthodontic Nozzle
- 1 x Plaque Brush Nozzle
- 1 x Periodontal Nozzle
- 1 x Tongue Cleaner Nozzle



Image 1: A set of eight Vimmk Y208 water flosser replacement nozzles. The top row shows four standard nozzles with colored rings (blue, green, black, grey) at their base. The bottom row displays specialized nozzles: an orthodontic nozzle, a plaque brush nozzle, another standard nozzle, and a tongue cleaner nozzle. All nozzles are transparent with a curved tip.

3. SETUP

To use a replacement nozzle, firmly connect it to the designated port on your Vimmk Y208 portable water flosser. Ensure the nozzle is securely seated to prevent leakage during operation.

4. OPERATING INSTRUCTIONS

Once the nozzle is attached to your water flosser, fill the water reservoir and select your desired pressure setting. Direct the water stream from the nozzle towards your gumline and between your teeth to remove food particles and plaque. Follow the specific instructions provided with your Vimmk Y208 water flosser for detailed operation.

- **Standard Nozzles:** For general daily cleaning and plaque removal.
- **Orthodontic Nozzle:** Designed for cleaning around braces and orthodontic appliances.
- **Plaque Brush Nozzle:** Features bristles for enhanced plaque removal in hard-to-reach areas.
- **Periodontal Nozzle:** Delivers a gentle stream for cleaning periodontal pockets.
- **Tongue Cleaner Nozzle:** For cleaning the tongue to help freshen breath.

5. MAINTENANCE

For optimal hygiene and performance, it is recommended to replace your water flosser nozzles every 3 months. After each use, rinse the nozzle thoroughly under running water to remove any debris. Store nozzles in a clean, dry place.

6. SAFETY WARNINGS

Attention: This oral irrigation device and its accessories are intended for use exclusively under adult supervision. Keep out of reach of children. Do not use if you have certain medical conditions without consulting a healthcare professional. Ensure you follow all usage and maintenance instructions to prevent damage or injury. If your water flosser is battery-powered, use only new, quality batteries as specified in its main manual.

7. TROUBLESHOOTING

These replacement nozzles are passive components. If you experience issues such as reduced water pressure, leakage, or improper fit, please refer to the troubleshooting section of your Vimmk Y208 water flosser's main instruction manual. Ensure the nozzle is correctly attached and not damaged. If problems persist, contact Vimmk customer support.

8. SPECIFICATIONS

- **Brand:** Vimmk
- **Model:** Y208 (Compatible with Vimmk Y208 portable water flosser)
- **Item Model Number:** 0767776761169
- **UPC:** 767776761169
- **Package Dimensions:** 13.5 x 7.1 x 3.9 cm
- **Package Weight:** 60 grams
- **Quantity:** 8 replacement nozzles per pack
- **Recommended Replacement:** Every 3 months

9. WARRANTY INFORMATION

Specific warranty details for replacement nozzles are typically covered under the warranty of the main Vimmk Y208 water flosser unit. Please refer to the warranty documentation provided with your original water flosser purchase or contact Vimmk customer service for detailed warranty terms and conditions.

10. SUPPORT INFORMATION

For further assistance, questions, or concerns regarding your Vimmk Y208 Water Flosser Replacement Nozzles, please contact Vimmk customer support through their official website or the retailer from whom you purchased the product.