

[Manuals.plus](#) /

› [upbright](#) /

› UpBright AC/DC Adapter Instruction Manual for Smith & Nephew Renasys Touch NPWT Devices

upbright 66801280 66801281

UpBright AC/DC Adapter Instruction Manual for Smith & Nephew Renasys Touch NPWT Devices

1. PRODUCT OVERVIEW

This UpBright AC/DC Adapter is designed to provide reliable power to compatible Smith & Nephew Renasys Touch Negative Pressure Wound Therapy (NPWT) devices. It functions as a power supply and battery charger, ensuring your medical device operates correctly.

Compatibility: This adapter is compatible with Smith & Nephew Renasys Touch NPWT devices, including models REF 66801280, 66801281, 66802134, REF 66801285, REF 66801286, and MENB1050A1900N02.



Image 1.1: Example of a compatible Smith & Nephew Renasys Touch NPWT device. The adapter is sold separately.

UpBright AC Adapter Power Supply Compatible with:

(Machine not Included)



Image 1.2: The UpBright AC/DC Adapter with its detachable power cord.

2. SAFETY INFORMATION

Please read and understand all safety instructions before using this product. Failure to follow these instructions may result in electric shock, fire, or damage to the product or connected device.


- **Indoor Use Only:** This adapter is designed for indoor use in dry locations. Do not expose to rain, moisture, or extreme temperatures.
- **Proper Ventilation:** Ensure the adapter has adequate ventilation during use. Do not cover it or place it in an enclosed space where heat cannot dissipate.
- **Voltage Compatibility:** Verify that the wall outlet voltage matches the input voltage specified on the adapter (100-240V AC).
- **Cable Integrity:** Regularly inspect the power cord and adapter for any signs of damage, such as frayed wires or cracked casing. Do not use if damaged.
- **Disconnect Safely:** Always grasp the plug, not the cord, when disconnecting from the wall outlet or the device.
- **Keep Away from Children:** Store the adapter and cords out of reach of children.
- **Medical Device Caution:** This adapter is intended for use with specific medical devices. Use only with the compatible Smith & Nephew Renasys Touch NPWT devices as specified.

3. PACKAGE CONTENTS

Upon opening the package, please ensure all items listed below are present and in good condition:

- 1 x UpBright AC/DC Adapter
- 1 x AC Power Cord (suitable for your region)

Come With A Suitable AC Power Cord (Only One)
to Fit Your Country/Area



THIS PRODUCT COMES WITH A UPBRIGHT®

LIMITED WARRANTY CERTIFICATE
+ UpBright® LIMITED WARRANTY



Image 3.1: The package includes one AC power cord, with the plug type varying by region (e.g., US, EU, UK, AU). A limited warranty certificate is also included.

4. SETUP

Follow these steps to properly set up your UpBright AC/DC Adapter:

1. **Connect AC Power Cord:** Insert the AC power cord firmly into the adapter's AC input port.
2. **Connect to Wall Outlet:** Plug the AC power cord into a standard wall electrical outlet. Ensure the outlet is easily accessible.
3. **Connect to Device:** Insert the DC output connector (barrel connector) from the adapter into the power input port of your Smith & Nephew Renasys Touch NPWT device.
4. **Verify Connection:** Ensure all connections are secure. The NPWT device should indicate that it is receiving power or charging.



Image 4.1: The AC/DC adapter connected to a wall outlet, ready for use.

5. OPERATING INSTRUCTIONS

Once properly connected, the UpBright AC/DC Adapter will automatically supply power to your Smith & Nephew Renasys Touch NPWT device. If the device has a rechargeable battery, the adapter will also charge the battery.

- Ensure the adapter is connected to a live power outlet and the NPWT device.
- Refer to your Smith & Nephew Renasys Touch NPWT device's user manual for specific operating instructions of the medical device itself.
- Disconnect the adapter from the wall outlet when not in use for extended periods.

6. MAINTENANCE

Proper maintenance ensures the longevity and safe operation of your AC/DC adapter.

- **Cleaning:** Disconnect the adapter from all power sources before cleaning. Use a soft, dry cloth to wipe the exterior.

Do not use liquid cleaners or abrasive materials.

- **Storage:** Store the adapter in a cool, dry place away from direct sunlight and moisture when not in use. Avoid tightly winding the cables to prevent damage.
- **Inspection:** Periodically inspect the adapter and power cords for any signs of physical damage, such as cuts, kinks, or discoloration. If damage is observed, discontinue use and replace the adapter.

7. TROUBLESHOOTING

If you encounter issues with your UpBright AC/DC Adapter, refer to the following common troubleshooting steps:

Problem	Possible Cause	Solution
Device not powering on or charging.	No power from wall outlet. Loose connection. Damaged adapter or cable. Incompatible device.	Test the wall outlet with another device. Ensure all connections (AC cord to adapter, adapter to wall, adapter to device) are secure. Inspect adapter and cables for damage. Replace if damaged. Confirm your device model is listed as compatible.
Adapter feels hot during use.	Normal operation (slight warmth). Poor ventilation. Overload or internal fault.	A slight warmth is normal. Ensure the adapter is in an open area with good airflow. Do not cover. If excessively hot or emitting smoke/smell, immediately disconnect from power and discontinue use.

If the problem persists after following these steps, please contact UpBright customer support.

8. SPECIFICATIONS

Feature	Detail
Brand	UpBright
Input Voltage	100-240V AC (Global Compatibility)
Power Source	AC/DC
Connector Type	Barrel Connector
Compatible Devices	Smith & Nephew Renasys Touch NPWT Devices (Models: 66801280, 66801281, 66802134, REF 66801285, REF 66801286, MENB1050A1900N02)
Product Dimensions	15.24 x 10.16 x 5.08 cm (approx. 6 x 4 x 2 inches)
Product Weight	255.15 g (approx. 0.56 lbs)

9. WARRANTY AND SUPPORT

This UpBright AC/DC Adapter comes with a Limited Warranty Certificate. Please refer to the certificate included in your package for specific terms and conditions regarding warranty coverage.

For technical support, warranty claims, or any questions regarding your UpBright product, please contact UpBright customer service through their official website or the retailer where the product was purchased.

Please retain your purchase receipt and warranty certificate for future reference.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0D3L67TN5>

More documents and user manuals: <https://manuals.plus/>