

[Manuals.plus](#) /

› [SuperHandy](#) /

› SuperHandy GoRise LT Portable Floor & Bath Lift User Manual

SuperHandy GBTS025

SuperHandy GoRise LT Portable Floor & Bath Lift User Manual

Model: GBTS025

INTRODUCTION

This manual provides essential information for the safe and effective use of your SuperHandy GoRise LT Portable Floor & Bath Lift. Please read this manual thoroughly before operating the device. This lift is designed to assist individuals with mobility challenges in transferring from floor level to a seated position, and for use in wet environments such as bathtubs, offering enhanced stability and ease of use.





Figure 1: SuperHandy GoRise LT Portable Floor & Bath Lift.

SAFETY GUIDELINES

- **Weight Capacity:** Do not exceed the maximum weight capacity of 330 lbs (150 kg). Overloading can cause damage to the lift and potential injury.
- **Stable Surface:** Always operate the lift on a flat, stable, and non-slip surface. Ensure the base is fully extended and stable before use.
- **Supervised Use:** The lift should always be used under the supervision of a trained caregiver. Do not allow unsupervised use by individuals who cannot operate it safely.
- **Waterproof Rating:** The device has an IP68 waterproof rating, making it suitable for wet environments. However, ensure the remote control's charging port is sealed when in use near water.
- **Clear Area:** Ensure the area around the lift is clear of obstructions during operation to prevent pinching or entanglement.
- **Battery Safety:** Use only the provided charger. Do not expose the battery or charger to extreme temperatures.

Disconnect the charger when not in use.

- **Regular Inspection:** Periodically inspect the lift for any signs of wear, damage, or loose components. Do not use if damaged.

PRODUCT FEATURES

- **Stable & Secure Lifting:** Lifts up to 330 lbs with a wider base measuring 23.6" x 16.1" (Depth x Width) for enhanced stability compared to standard models.
- **Compact Portable Design:** This full-body foldable structure is lightweight, weighing only 19 lbs. When folded, it measures 28.9" (Length) x 17.7" (Width) x 4.9" (Height), making it easy to store and transport.
- **IP68 Waterproof:** Suitable for wet areas like bathrooms, with a battery-powered design that ensures functionality is not compromised by water.
- **Ergonomic Comfort:** Offers a wider seat and intuitive control for height adjustment. Larger, strategically placed drain holes make our lift ideal for wet environments, ensuring efficient drainage and avoiding build-up.
- **Versatile with Removable Side Panels:** Includes removable side panels for easier lateral transfers, providing additional support and easy detachment for unrestricted access.



Figure 2: IP68 Waterproof design for use in wet environments.




Figure 3: Lift capacity and dimensions.

COMPONENTS

The SuperHandy GoRise LT lift consists of the main lift unit, a wired remote control, and removable side panels.

- **Main Lift Unit:** Comprises the seat, backrest, and the scissor-lift mechanism with a stable base.
- **Wired Remote Control:** Features power, up, and down buttons for intuitive operation. It houses the rechargeable lithium-ion battery.
- **Removable Side Panels:** Two panels that attach to the sides of the seat for additional support during transfers.

Portable Design

Our lift features a complete folding design, rubber latch system, and removable side panels, weighing just 19 lbs. Making it easy to transport and store.

Lightweight: 19 LBS

Folded Size: 28x 17.7 x 4.9 Inches



Figure 4: Wired Remote Control for the GoRise LT Lift.

SETUP & ASSEMBLY

The SuperHandy GoRise LT lift comes largely pre-assembled. Follow these steps for initial setup:

1. **Unpack:** Carefully remove the lift from its packaging.
2. **Connect Remote Control:** Locate the waterproof connector for the remote control cable and securely attach it to the corresponding port on the lift's main unit. Ensure a tight seal.
3. **Unfold the Lift:** Place the folded lift on a flat surface. Gently lift the backrest until the seat and base unfold into the operational position. Ensure all locking mechanisms are engaged.
4. **Attach Side Panels (Optional):** If desired, slide the removable side panels into the designated slots on either side of the seat until they click securely into place.
5. **Initial Charge:** Before first use, fully charge the lift's battery using the provided charger. The remote control will indicate charging status.

Your browser does not support the video tag.

Video 1: SuperHandy GUT170 floor lift Installation and Operation. This video demonstrates the connection of the remote, unfolding the lift, and attaching side panels.

OPERATING INSTRUCTIONS


The SuperHandy GoRise LT lift is operated using the wired remote control.

1. **Power On:** Press the yellow power button on the remote control to activate the lift. The indicator light will illuminate.
2. **Positioning:** Carefully position the lift near the individual requiring assistance. Ensure the individual is safely seated on the lift's seat.
3. **Lifting:** Press and hold the green "Up" arrow button on the remote control to raise the lift. Release the button when the desired height is reached.
4. **Lowering:** Press and hold the green "Down" arrow button on the remote control to lower the lift. Release the button when the desired height is reached, or when the lift is fully lowered to the floor.
5. **Transfer:** Once at a suitable height, assist the individual in safely transferring off the lift.
6. **Power Off:** After use, press the yellow power button to turn off the lift.



Figure 5: Demonstrating the lift in use for floor transfers.

One Step Folding



Compact Design

Figure 6: Step-by-step guide for using the lift.

CARE & MAINTENANCE

- **Cleaning:** Clean the lift regularly with a mild soap and water solution. Wipe dry with a soft cloth. Avoid abrasive cleaners or solvents.
- **Battery Charging:** Charge the lithium-ion battery when the low power indicator on the remote control illuminates. A full charge allows for approximately 30 complete up and down cycles.
- **Storage:** When not in use, fold the lift for compact storage. Ensure the remote control cable is neatly coiled and secured. Store in a dry place away from direct sunlight and extreme temperatures.





330 lbs

19 lbs

Aluminum

2 inches

30 Cycles

17.7" inches

16.1" inches

VS

Other Brands



300 lbs

30 lbs+

Plastic

3.4 inches+

10 Cycles

14.5" inches

12" inches

	Weight Capacity	
	Product Weight	
	Support Frame	
	Seat Height Clearance	
	Cycles Per Charge	
	Seat Width	
	Base Width	

Figure 7: The lift's compact and portable design for easy storage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Lift does not operate.	Battery is low or depleted. Remote control not properly connected. Power button not pressed.	Charge the battery. Ensure the remote control cable is securely connected. Press the yellow power button.
Lift operates slowly.	Low battery charge.	Recharge the battery.
Side panels are loose.	Not fully engaged.	Push firmly until they click into place.

If you encounter issues not listed here, or if solutions do not resolve the problem, please contact SuperHandy customer support.

TECHNICAL SPECIFICATIONS

Specification	Value
Model Number	GBTS025
Maximum Weight Recommendation	330 Pounds (150 kg)
Item Weight	18.92 pounds (8.6 kg)
Product Dimensions (Unfolded)	22"D x 29.5"W x 47.6"H
Folded Dimensions	28.9"L x 17.7"W x 4.9"H
Seat Height Range	2 inches (min) - 20 inches (max)
IP Rating	IP68 Waterproof
Battery Type	Lithium Ion (25.9V, 2000mAh, 51.8Wh)
Cycles Per Charge	Up to 30 cycles
Material	Aluminum, Plastic

WARRANTY & SUPPORT

The SuperHandy GoRise LT Portable Floor & Bath Lift comes with a manufacturer's warranty. Please refer to the warranty card included with your product for specific terms and conditions.

For technical support, replacement parts, or any inquiries, please contact SuperHandy customer service.

Brand: SuperHandy

Website: [Visit the SuperHandy Store on Amazon](#)

Phone: Refer to product packaging or manufacturer's website for current contact information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0D4B2RVPY>

More documents and user manuals: <https://manuals.plus/>