

[Manuals.plus](#) /

› [YOUSHARES](#) /

› YOUSHARES Insulin Cooler Travel Case with Temperature Display User Manual

YOUSHARES MPN_ISL_Case_H2F_BL_01

YOUSHARES Insulin Cooler Travel Case with Temperature Display

Model: MPN_ISL_Case_H2F_BL_01

Brand: YOUSHARES

INTRODUCTION

This manual provides essential instructions for the proper use and care of your YOUSHARES Insulin Cooler Travel Case. Designed for individuals requiring temperature-sensitive medication storage, this case offers portability and real-time temperature monitoring for insulin, GLP-1 medications, and other injectable peptides. Please read this manual thoroughly before first use to ensure optimal performance and longevity of your product.

PRODUCT FEATURES

- **Real-Time Temperature Monitoring:** Features a built-in Fahrenheit temperature display for continuous monitoring without opening the case.
- **Compact and Spacious Design:** Measures 8.5 × 3.5 × 2.3 inches and weighs 0.6 lb, designed to hold up to two insulin pens or GLP-1 injection pens along with ice packs.
- **Effective Insulation:** Constructed from durable PU, EVA, and aluminum foil materials to maintain a safe temperature range for medications.
- **TSA-Approved Ice Packs:** Includes two reusable ice packs, designed to keep medications cool for 6–8 hours after 8–12 hours of freezing. Cooling duration may vary with ambient temperature.
- **Durable Construction:** Waterproof PU cloth layer, shockproof EVA protective layer, and superior thermal aluminum insulation layer.



Image: The YOUSHARES Insulin Cooler Travel Case displaying its compact dimensions and portable design.



Image: A diagram illustrating the multi-layer insulation of the YOUSHARES Insulin Cooler Travel Case, including PU cloth, EVA protective layer, and aluminum insulation layer.

SETUP

1. **Unpack Contents:** Remove the travel case and two reusable ice packs from the packaging.
2. **Prepare Ice Packs:** Place both reusable ice packs in a freezer for 8 to 12 hours until completely frozen. Ensure they are flat to freeze evenly.
3. **Insert Ice Packs:** Once frozen, insert the ice packs into the designated compartments within the travel case. The case is designed to hold one ice pack on each side.
4. **Check Temperature Display:** The built-in temperature display will activate automatically. Verify that it shows a temperature within the safe range for your medication after the ice packs have been inserted for a short period.



Image: Close-up view of the YOUSHARES Insulin Cooler Travel Case's real-time temperature display, showing the internal temperature.

Your browser does not support the video tag.

Video: Demonstration of how to prepare and use the YOUSHARES Insulin Cooler Travel Case for keeping medication cool on the go.

OPERATING INSTRUCTIONS

1. **Load Medication:** After the ice packs are in place and the case is cool, carefully place your insulin pens, GLP-1 pens, or other temperature-sensitive medications into the elastic loops or mesh pockets inside the case. Ensure they are secure and not directly touching the frozen ice packs to prevent freezing.
2. **Close Securely:** Zip the case closed completely to maintain insulation.
3. **Monitor Temperature:** Periodically check the external temperature display to ensure your medication remains within its required temperature range.
4. **Travel Considerations:** The included ice packs are TSA-compliant for air travel. Always confirm current TSA regulations regarding medical supplies before flying.



Image: The interior of the YOUSHARES Insulin Cooler Travel Case, showing how insulin pens and other diabetic supplies can be organized.



Image: The YOUSHARES Insulin Cooler Travel Case positioned on an airplane tray table, illustrating its suitability for air travel.

Your browser does not support the video tag.

Video: A detailed review of the YOUSHARES Insulin Cooler Travel Case, highlighting its features and temperature alarm functionality.

MAINTENANCE

- **Cleaning the Case:** Wipe the exterior of the case with a damp cloth and mild soap as needed. Do not submerge the case in water. Allow it to air dry completely before storing or reuse.
- **Ice Pack Care:** The reusable ice packs can be washed with soap and water. Ensure they are completely dry before freezing again. Store them in a cool, dry place when not in use.
- **Temperature Display Battery:** The temperature display is powered by a small button cell battery. If the display becomes dim or stops working, replace the battery. Refer to the small battery compartment on the display unit for specific battery type and replacement instructions.
- **Storage:** Store the case in a cool, dry place away from direct sunlight when not in use.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Medication not staying cool for 6-8 hours	Ice packs not fully frozen; high ambient temperature; case not fully closed.	Ensure ice packs are frozen for 8-12 hours. Avoid prolonged exposure to high temperatures. Ensure zipper is fully closed.
Temperature display is blank or dim	Battery needs replacement.	Replace the button cell battery in the display unit.
Case appears damaged	Physical impact or improper handling.	Inspect for cracks or tears. If insulation is compromised, cooling performance may be affected. Contact customer support if damage impacts functionality.

SPECIFICATIONS

- **Product Dimensions:** 8.5 x 3.5 x 2.3 inches (21.5 x 9 x 5.8 cm)

- **Item Weight:** 9.52 ounces (270 g)
- **Materials:** Ethylene Vinyl Acetate (EVA), Polyurethane (PU), Aluminum Foil
- **Closure Type:** Zipper
- **Number of Compartments:** 1 main compartment with internal organization
- **Water Resistance Level:** Waterproof (exterior)
- **Cooling Duration:** Up to 6-8 hours (with frozen ice packs, ambient temperature 75°F–82°F / 24°C–28°C)
- **Included:** 1 x Insulin Pen Travel Case, 2 x Reusable Ice Packs

WARRANTY AND SUPPORT

YOUSHARES products are designed for reliability and performance. While specific warranty details are not provided in this manual, please retain your proof of purchase. For any product inquiries, technical support, or warranty claims, please contact YOUSHARES customer service through the retailer where the product was purchased or visit the official YOUSHARES brand store online.

For additional information and support, you may visit the [YOUSHARES Store on Amazon](#).

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0D6VPSDGL>

More documents and user manuals: <https://manuals.plus/>