

[Manuals.plus](#) /

› [XSOUL](#) /

› XSOUL IPL Hair Removal Device with Cooling - User Manual

XSOUL B0D8W5WH7M

XSOUL IPL Hair Removal Device with Cooling - User Manual

Model: B0D8W5WH7M

1. OVERVIEW

The XSOUL IPL Hair Removal Device with Cooling is designed for effective and comfortable at-home hair removal for both women and men. Utilizing advanced IPL (Intense Pulsed Light) technology, this device targets hair follicles to disrupt the hair growth cycle, leading to long-lasting hair reduction. Its integrated ice cooling system ensures a painless experience by maintaining a comfortable skin temperature during treatment.



Image: The XSoul IPL Hair Removal Device, showcasing its ergonomic design.

Key Features:

- **Advanced IPL Technology:** Features a red color wavelength (470 nm to 1200 nm) for effective hair follicle targeting.
- **Ice Cooling Technology:** Maintains skin temperature at approximately 50°F (10°C) to minimize discomfort and prevent redness.
- **9 Adjustable Energy Levels:** Allows customization for different skin types and hair removal areas, with higher levels offering greater efficacy.
- **Dual Operation Modes:** Includes Manual Mode for precise treatment on small areas (e.g., bikini line, armpits, face) and Auto Mode for efficient treatment on larger areas (e.g., arms, legs, back).
- **Long Lifespan:** Provides up to 999,999 flashes for extensive use.

2. IMPORTANT SAFETY INFORMATION

Read all instructions carefully before using the device. Failure to follow safety guidelines may result in injury or device damage.

General Safety Precautions:

- Always wear the provided protective glasses during operation to shield your eyes from the intense light.
- Do not use the device directly on tattoos, dark brown, black, or very dark skin tones, as this may cause skin discomfort or burns due to excessive light absorption.
- Avoid prolonged use or multiple light flashes on the same skin area during a single session to prevent skin irritation.
- Do not use on injured, irritated, or sunburnt skin.
- Consult a physician before use if you have any underlying skin conditions or are taking photosensitizing medications.
- Keep out of reach of children.

Skin and Hair Type Compatibility:

The effectiveness and safety of IPL hair removal depend on your skin tone and hair color. Refer to the chart below to determine if this device is suitable for you.



Image: Chart illustrating suitable skin tones and hair colors for IPL treatment.

- **Safe and Effective:** Light to medium skin tones with dark hair (black, dark brown, brown, light brown).
- **Unsafe/Ineffective:** Dark skin tones (dark brown, black) and very light hair colors (red, light blonde, white, grey).

3. PACKAGE CONTENTS

Verify that all items are present in the package upon unboxing:

- XSOUL IPL Hair Removal Device
- User Manual
- Protective Glasses
- Razor
- Power Adapter

Complete Accessories

Everything You Need Included



Image: All components included in the XSOUL IPL Hair Removal Device package.

4. DEVICE COMPONENTS

Familiarize yourself with the different parts of your XSOUL IPL device:

- **Flash Window:** Emits the IPL light.
- **Control Panel:** Displays energy level, mode, and remaining flashes.
- **Power Button:** Turns the device on/off and adjusts energy levels.
- **Flash Button:** Initiates a single flash in Manual Mode.
- **Ice Cooling Button:** Activates/deactivates the cooling function.
- **Air Vents:** For heat dissipation.



Image: Close-up of the device highlighting the 50°F Ice-Cooling feature.

5. SETUP

Follow these steps to prepare your device for use:

1. **Connect Power:** Plug the power adapter into the device and then into an electrical outlet. The device will power on automatically.
2. **Wear Protective Glasses:** Always put on the provided protective glasses before turning on the device or beginning treatment.

6. OPERATING INSTRUCTIONS

Preparation for Treatment:

1. **Shave the Area:** Ensure the skin area to be treated is clean, dry, and freshly shaved. Do not wax, pluck, or use depilatory creams, as these methods remove the hair follicle, which is necessary for IPL to be effective.
2. **Clean Skin:** Make sure there are no lotions, creams, or deodorants on the skin.

Device Operation:

1. **Power On:** Long press the Power Button to turn on the device. The display will light up.
2. **Activate Cooling (Optional):** Press the Ice Cooling Button to activate the cooling function. The cooling indicator will illuminate, and the cooling plate will begin to chill.
3. **Adjust Energy Level:** Short press the Power Button to cycle through the 9 energy levels. Start with level 1 for your first treatment and gradually increase as your skin adapts. Higher levels provide more effective results.
4. **Select Mode:**
 - **Manual Mode:** For small, sensitive areas like the bikini line, armpits, and face. Short press the Flash Button to emit a single flash.
 - **Auto Mode:** For larger areas such as arms, legs, and back. Long press the Flash Button to activate Auto Mode. The device will automatically flash when the flash window is in full contact with the skin.

Dual Flash Modes



Image: Illustration of Manual Flash for small areas and Auto Flash for large areas.

5. **Apply to Skin:** Place the flash window firmly against your skin, ensuring full contact. The device will not flash if there is insufficient contact.
6. **Perform Treatment:**
 - In Manual Mode, press the Flash Button for each flash, moving the device to the next area after each flash.
 - In Auto Mode, glide the device slowly and continuously over the skin area. The device will flash automatically.
7. **Treatment Schedule:** For optimal results, use the device 2-3 times a week for the first 1-4 weeks, then 1-2 times a week for weeks 5-8. After 8 weeks, use once every 2-4 weeks as needed for maintenance. Individual results may vary.
8. **Power Off:** Long press the Power Button to turn off the device after your session.

Your browser does not support the video tag.

Video: Official XSOUL video demonstrating the IPL Hair Removal Device with Cooling in use, showing how to operate the device and its features.

7. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your device.

Cleaning:

- Always unplug the device before cleaning.
- Use a soft, dry cloth to wipe the flash window and the device body. Do not use abrasive cleaners or harsh chemicals.
- Ensure no liquid enters the device.

Storage:

- Store the device in a cool, dry place, away from direct sunlight and moisture.
- Keep the device in its original packaging or a protective case to prevent damage.

8. TROUBLESHOOTING

If you encounter any issues with your device, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Power adapter not properly connected or power outlet issue.	Ensure the power adapter is securely plugged into both the device and a working electrical outlet.
Device does not flash.	Flash window not in full contact with skin, or incorrect mode selected.	Ensure the flash window is fully pressed against the skin. Check if the correct mode (Manual or Auto) is selected for the area.
Skin feels uncomfortable or hot.	Energy level is too high for your skin type, or cooling function is off.	Reduce the energy level. Activate the ice cooling function. Avoid treating the same area multiple times in one session.
No visible hair reduction after several weeks.	Inconsistent treatment schedule, incorrect energy level, or unsuitable skin/hair type.	Follow the recommended treatment schedule. Increase energy level if comfortable. Re-check skin and hair type compatibility.

9. SPECIFICATIONS

Technical details for the XSOUL IPL Hair Removal Device:

Feature	Detail
Model	B0D8W5WH7M
Brand	XSOUL
Technology	IPL (Intense Pulsed Light) with Ice Cooling
Wavelength	470 nm - 1200 nm
Energy Levels	9 adjustable levels (Max 19J)
Flashes	999,999
Operation Modes	Manual and Auto Flash

Feature	Detail
Cooling Temperature	Approximately 50°F (10°C)
Package Dimensions	8.7 x 6.97 x 3.78 inches
Weight	2.05 Pounds
UPC	850048244252
Manufacturer	XSOUL

10. WARRANTY & SUPPORT

XSOUL is committed to providing high-quality products and customer satisfaction.

Warranty Information:

- **36 Month Warranty:** Your device is covered by a 36-month manufacturer’s warranty from the date of purchase.
- **100 Day Money Back Guarantee:** If you are not completely satisfied with your purchase, you may be eligible for a full refund within 100 days.

Customer Support:

For any questions, concerns, or support needs regarding your XSOUL IPL Hair Removal Device, please contact our customer service team:

- **Email:** help@xsoulofficial.com

Our team is dedicated to resolving any issues and ensuring your satisfaction.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question