

Igria Callus Free

Igria Callus Free Foot Care Device

Model: Callus Free

1. INTRODUCTION

The Igria Callus Free device is designed for effective removal of dead skin and calluses from your feet. This rechargeable, ergonomic device features a built-in vacuum system for minimal mess and comes with three interchangeable grinding heads to suit various skin conditions. Please read this manual thoroughly before first use to ensure safe and optimal operation.

2. SAFETY INSTRUCTIONS

- Do not use the device on broken, irritated, or infected skin. Consult a medical professional if you have any skin conditions or concerns.
- This device is not suitable for individuals with diabetes or poor circulation. Consult a doctor before use.
- Keep the device out of reach of children.
- Ensure the device is fully dry before charging.
- Do not immerse the main body of the device in water beyond its IPX7 water-resistant rating. The grinding heads are washable.
- Avoid applying excessive pressure, as this may cause skin irritation or injury.
- Discontinue use immediately if you experience pain, discomfort, or adverse reactions.
- Use only the provided charging cable.

3. PRODUCT OVERVIEW

The Igria Callus Free device is designed for comfortable and efficient foot care. It features a rechargeable battery for cordless operation and a built-in vacuum to collect skin debris.

Package Contents:

Easy to Clean



Figure 1: Contents of the Igia Callus Free package, including the main device, three grinding heads, a charging cable, and a cleaning brush.

- 1 x Igia Callus Free Device with built-in vacuum
- 3 x Diamond Crystal Grinding Heads (Coarse, Regular, Fine)
- 1 x USB Charging Cable
- 1 x Cleaning Brush

Device Features:



Figure 2: The Igia Callus Free device in operation.

- **Ergonomic Design:** Comfortable to hold and use.
- **Rechargeable Battery:** Provides cordless convenience.
- **Built-in Vacuum:** Collects dead skin for a cleaner experience.
- **IPX7 Water-Resistant:** Allows for easy rinsing of the device head.
- **Interchangeable Grinding Heads:** Three types for customized exfoliation.



Coarse

For hard calluses and dead skin.



Regular

For newly formed calluses and normal skin care.



Fine

For hands or more delicate skin care.



Figure 3: Device dimensions for reference.

4. SETUP

4.1 Charging the Device

1. Before first use, fully charge the device.
2. Connect the USB charging cable to the charging port on the device and plug the other end into a suitable USB power adapter (not included) or a computer USB port.
3. The charging indicator light will illuminate during charging and change color or turn off when fully charged.
4. A full charge typically takes approximately 2-3 hours.

4.2 Attaching Grinding Heads

1. Ensure the device is turned off.
2. Select the desired grinding head (Coarse, Regular, or Fine).
3. Align the grinding head with the device's rotating mechanism and gently press it into place until it clicks securely.

4. To remove a grinding head, gently pull it away from the device.

5. OPERATING INSTRUCTIONS

Follow these steps for effective and safe use of your Igia Callus Free device:

1. Ensure your feet are clean and dry before use.
2. Select and attach the appropriate grinding head for your needs (refer to section 5.1).
3. Press the power button to turn on the device.
4. Gently glide the rotating head over the callused or rough areas of your feet. Do not apply excessive pressure.
5. Move the device in short, light strokes, avoiding staying in one area for too long to prevent irritation.
6. The built-in vacuum will collect the removed skin debris.
7. Continue until the desired smoothness is achieved.
8. After use, press the power button to turn off the device.
9. Moisturize your feet after treatment for best results.



Figure 4: Close-up of the device's spinning head in action.

5.1 Choosing the Right Grinding Head



Figure 5: Types of grinding heads and their applications.

- **Coarse Head:** Designed for removing stubborn dead skin and hard calluses.
- **Regular Head:** Suitable for tackling rough dead skin and newly formed calluses.
- **Fine Head:** Ideal for daily foot care and maintaining smooth skin, or for more delicate areas.



Figure 6: Illustrative results of using the device.

6. MAINTENANCE

6.1 Cleaning the Device

Regular cleaning ensures hygiene and extends the lifespan of your device.

1. Turn off the device and remove the grinding head.
2. Empty the removable vacuum bin by twisting and pulling it off. Dispose of the collected skin debris.
3. Use the provided cleaning brush to remove any remaining debris from the grinding head and the device's vacuum area.
4. Rinse the grinding head and the removable vacuum bin under running water.
5. The device itself is IPX7 water-resistant, meaning it can be rinsed under running water. Ensure the charging port cover is securely closed before rinsing. Do not submerge the entire device for extended periods.
6. Allow all parts to air dry completely before reassembling or storing.



Figure 7: Rinsing the grinding head under water.

6.2 Storage

Store the device in a cool, dry place, away from direct sunlight and extreme temperatures. Ensure it is fully charged before long-term storage to maintain battery health.

7. TROUBLESHOOTING

- **Device not turning on:** Ensure the device is fully charged. If the battery is low, connect it to the charger.
- **Reduced performance:** The grinding head may be worn out and needs replacement, or the vacuum bin might be full. Clean the device and replace the head if necessary.
- **Device stops during use:** This may indicate low battery. Recharge the device. It could also be a safety feature if excessive pressure is applied. Reduce pressure.
- **Skin irritation:** Ensure you are not applying too much pressure or using the device on sensitive or broken skin. Use the appropriate grinding head for your skin type.

8. SPECIFICATIONS

Brand	Igia
Model	Callus Free
Color	White
Material	Plastic
Item Weight	3 kg
Package Dimensions	20 x 20 x 15 cm
Battery	Rechargeable (type not specified)
Water Resistance	IPX7 (device head)

9. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation provided at the point of purchase or contact your retailer. Spare part availability information is not provided in this manual. For inquiries regarding spare parts or additional grinding heads, please contact Igia customer service or your authorized dealer.

