

## Aopvui BB

# Aopvui IPL Laser Hair Removal Device User Manual

**Brand:** Aopvui

**Model:** BB

## INTRODUCTION

The Aopvui IPL Laser Hair Removal Device is designed for permanent hair reduction at home, offering a safe and comfortable experience. Utilizing Intense Pulsed Light (IPL) technology with a Sapphire Ice-Cooling function, it aims to minimize discomfort during treatment by cooling the skin to approximately 5°C (41°F). This device features multiple energy levels and smart modes tailored for different body areas, providing an efficient and convenient alternative to salon treatments.

## IMPORTANT SAFETY INFORMATION

Please read all instructions carefully before using the device. Failure to follow safety guidelines may result in injury or device malfunction.

- **Skin Tone Suitability:** This device is NOT suitable for black or dark brown skin tones. Darker skin contains more melanin, which can absorb too much IPL energy, potentially causing discoloration or burns. Refer to the skin tone chart below to determine suitability.
- **Hair Color Suitability:** Not suitable for light blonde, red, or white/gray hair. IPL targets melanin in the hair follicle, and these hair colors lack sufficient melanin for effective treatment.
- **Pre-Treatment:** Always shave the hair on the skin surface before use to prevent burning.
- **Eye Protection:** Always wear the provided safety glasses during operation to protect your eyes from the intense light.
- **Device Placement:** Ensure the device is placed vertically at 90 degrees to the skin surface for proper function and safety. The device will only flash when it detects full skin contact.
- **Corded Operation:** This is a corded device and must be plugged in during use. It is not wireless.
- **Not Suitable For:** Individuals with tattoos, skin diseases, sunburn, or pregnant women.

# Three Smart Modes

give you precise and painless hair removal, safe for the whole body



## Face Mode

Suitable for Fine Soft Hair



## Body Mode

Suitable for Regular to Dense Hair



## Bikini Mode

Suitable for Coarse Stubborn Hair



Figure 1: Skin Tone and Hair Color Suitability Chart. Green check indicates applicable, red cross indicates not applicable.



Figure 2: Conditions and areas where the IPL device should not be used, including irritated skin, sunburn, pregnancy, moles, freckles, and tattoos.

## WHAT'S INCLUDED

Your Aopvui IPL Laser Hair Removal Device package includes the following items:

- Aopvui IPL Laser Hair Removal Device

- Power Adapter
- Safety Glasses
- Cleaning Cloth (for light window)
- Two Razors (for pre-treatment shaving)
- User Manual

## PRODUCT OVERVIEW

---

The Aopvui IPL device is designed for ease of use and effectiveness. Key features include:

- **Sapphire Ice-Cooling Function:** The light head cools to 5°C (41°F) upon skin contact, providing a comfortable and painless experience by protecting the epidermis from heat.
- **9 Energy Levels:** Adjustable intensity levels to suit different skin sensitivities and hair types. Higher levels generally yield better results.
- **3 Smart Modes:** Optimized settings for specific body areas:
  - **Face Mode:** Suitable for fine, soft hair.
  - **Body Mode:** Suitable for regular to dense hair.
  - **Bikini Mode:** Suitable for coarse, stubborn hair.
- **2 Flash Modes:**
  - **Automatic Mode:** For large areas like arms, legs, and back, allowing for faster treatment.
  - **Manual Mode:** For small, precise areas like lips, bikini line, and underarms.
- **Unlimited Flashes:** Designed for long-term use, sufficient for whole family treatment.
- **LCD Display:** Clear display showing adjustable modes and energy levels.



Figure 3: Front view of the Aopvui IPL Laser Hair Removal Device with its display screen.

# Safe and Effective **Sapphire Ice Cooling** Function



Speedy Cooling



Lowers to 41°F



Long-lasting Cooling



Unlimited Flashes



Figure 4: Illustration of the Speedy Cooling, 41°F temperature, Long-lasting Cooling, and Unlimited Flashes features.

# 360° Full Touch Sapphire Cooling System

The **12 cm<sup>2</sup>** sapphire contact surface enables comfortable and painless hair removal treatments on the whole body and sensitive areas

Aopvui



**41°F**

**12 cm<sup>2</sup>** full-touch sapphire

Ordinary



**67°F**

**3cm<sup>2</sup>** aluminum




Figure 5: Comparison of Aopvui's 12 cm<sup>2</sup> sapphire contact surface cooling to 41°F versus an ordinary 3 cm<sup>2</sup> aluminum surface cooling to 67°F.

## The Salon Experience Without the Salon Cost

|                     | Our Sapphire Ice Light Window      | Ordinary Ice                        | Other equipment              |
|---------------------|------------------------------------|-------------------------------------|------------------------------|
| Sapphire Ice Effect | 5°C/41°F With Ice cooling effect   | 25°C/77°F Ice cooling is incomplete | X                            |
| Convenience         | Use at Home Anytime                | Needed appointment                  | Inconvenient operation       |
| Cost                | Unlimited Flashes                  | X                                   | Expensive and short lifetime |
| Body Parts          | From head to toe                   | Pay per body parts                  | √                            |
| User Friendly       | 24 hours/7 day online team support | X                                   | X                            |




Figure 6: Visual representation of Face Mode (1-6 energy levels), Body Mode (1-9 energy levels), and Bikini Mode (1-3 energy levels).

## SETUP

Before your first use, ensure your skin is clean and dry. Follow these steps for initial setup:

1. **Shave the Area:** Use one of the provided razors or your preferred method to shave the hair from the area you intend to treat. This ensures the IPL energy is directed to the hair follicle, not the surface hair.
2. **Connect Power:** Plug the power adapter into the device and then into a wall outlet. The device is corded and requires continuous power during use.
3. **Wear Safety Glasses:** Put on the included safety glasses to protect your eyes from the bright flashes of light.
4. **Power On:** Press and hold the power button (located on the side of the device) until the LCD screen illuminates.

# Easy professional hair removal in just 3 steps

### 1. Shave



### 2. Choose mode



### 3. Start treatment



☒ Only 2-3 session per week

☒ Full body treatment in just 15 mins

☒ Nearly Painless even on sensitive areas

Figure 7: Visual guide for the three main steps: Shave, Choose Mode, and Start Treatment.

## OPERATING INSTRUCTIONS

Once the device is set up, follow these steps for effective hair removal:

1. **Select Energy Level:** The device defaults to Level 1. Press the power button briefly to cycle through the 9 energy levels. Start with Level 1 for your first treatment and gradually increase as your skin tolerates it.
2. **Choose Mode:** Use the mode selection button (often integrated with the main button or a separate one) to switch between Face, Body, and Bikini modes. The LCD screen will indicate the selected mode and recommended energy levels.
3. **Apply to Skin:** Place the light window of the device firmly and vertically (90 degrees) against the skin area you wish to treat. The device has a skin contact sensor and will only flash when full contact is detected.
4. **Initiate Flash:** Press the main flash button to emit a pulse of light. You will feel a cool sensation from the Sapphire Ice-Cooling function.
5. **Move to Next Area:** After each flash, lift the device slightly and move it to the next adjacent area. For larger areas, use the Automatic Mode by holding down the flash button and gliding the device over the skin. For smaller areas, use Manual Mode by pressing the button for each flash.
6. **Repeat:** Continue treating the desired area until you have covered it completely. Avoid overlapping flashes too much on the same spot.



Figure 8: Detailed steps for device usage, emphasizing shaving, eye protection, power connection, and proper skin contact.



Figure 9: Visual guide to suitable treatment areas including bikini line, legs, arms, armpits, back, and face.

## USAGE SCHEDULE AND EXPECTED RESULTS

---

Consistency is key for optimal results. Follow the recommended usage schedule:

- **Frequency:** Use the device 2-3 sessions per week for the initial phase.
- **Visible Results:** Many users experience visible hair reduction in just 4-6 weeks.
- **Long-Term:** Continued use will lead to permanent hair reduction as the IPL energy destroys hair follicles.

Individual results may vary based on hair type, skin tone, and adherence to the treatment schedule.

# Is Aopvui Right for my Skin?

## Step 1 Identify Your Skin Tone



## Step 2 Identify the Color of Your Body Hair



✓ Applicable

✗ Not Applicable

**Attention:** Not suitable for tattoos , skin disease,sunburn, Pregnant women

Figure 10: Visual timeline showing hair growth before treatment, after 4 weeks, and after 8 weeks, demonstrating visible hair reduction.

## MAINTENANCE

To ensure the longevity and optimal performance of your Aopvui IPL device, follow these maintenance guidelines:

- **Cleaning the Light Window:** After each use, gently wipe the light window with the provided cleaning cloth to remove any residue or hair. Ensure the device is powered off and unplugged before cleaning.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and moisture. Keep it in its original packaging or a protective case when not in use.
- **Avoid Liquids:** Do not immerse the device in water or any other liquid.

## TROUBLESHOOTING

If you encounter any issues with your Aopvui IPL device, refer to the common problems and solutions below:

| Problem                               | Possible Cause                                                                 | Solution                                                                                                                                                                                                                              |
|---------------------------------------|--------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Device does not power on.             | Not properly plugged in; power button not held long enough.                    | Ensure the power adapter is securely connected to the device and a working wall outlet. Press and hold the power button for a few seconds until the screen lights up.                                                                 |
| Device does not flash.                | Insufficient skin contact; incorrect skin tone/hair color; device malfunction. | Ensure the light window is in full, flat contact with the skin (90-degree angle). Verify your skin tone and hair color are suitable for IPL treatment (refer to Safety Information). If the issue persists, contact customer support. |
| Skin feels too warm or uncomfortable. | Energy level is too high; prolonged contact with one area.                     | Reduce the energy level. Ensure you are moving the device consistently and not holding it on one spot for too long. The Sapphire Ice-Cooling function should help mitigate discomfort.                                                |

## SPECIFICATIONS

- **Brand:** Aopvui
- **Model Name:** Aopvui (BB)
- **Power Source:** Corded Electric
- **Special Feature:** 2 Modes and Unlimited Light Flashes for Whole Family Use, High Efficiency with 9 Levels 3 Smart Body Part Choices, Speedy Sapphire Ice Cooling, Ergonomic Design, Enjoy Painless During Hair Removal Process
- **Recommended Uses For Product:** Armpit, Arm, Bikini Line, Legs, Back, Chest
- **Package Dimensions:** 8.7 x 6.26 x 4.21 inches
- **Item Weight:** 1.94 Pounds
- **Manufacturer:** Aopvui

## WARRANTY AND SUPPORT

For any product inquiries, technical assistance, or warranty claims, please contact Aopvui customer support. The product is supported by a 24 hours/7 day online team support, ensuring assistance is available whenever needed. Please refer to your purchase documentation for specific warranty details and contact information.

### ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

**Details**

Model number, symptoms, what you tried, or the section of the manual you are using.

**Submit question**