

[Manuals.plus](#) /

› [foreverLily](#) /

› foreverLily Rechargeable LED Facial Mask User Manual

foreverLily B0DDBZLYT7

foreverLily Rechargeable LED Facial Mask User Manual

Model: B0DDBZLYT7

INTRODUCTION

Thank you for choosing the foreverLily Rechargeable LED Facial Mask. This device utilizes advanced LED light technology to support various skin care needs, including firming, anti-aging, and overall skin appearance improvement. Designed for convenient home use, this manual provides essential information for safe and effective operation.

LED BEAUTY MASK

Restore Your Beauty



USB Charge



Cordless use



100%



Visual eye and
mouth window



270 LEDs



Image: The foreverLily LED Facial Mask highlighting its key features: USB charging, cordless operation, 100% UV-free light, clear eye and mouth openings, and 270 LED lights.

SAFETY INFORMATION

Please read all safety instructions before using the device. Failure to follow these instructions may result in injury or damage to the product.

- Do not use if you have a history of light-induced seizures or epilepsy.
- Consult a physician before use if you are pregnant, breastfeeding, or have any serious medical conditions.
- Avoid direct eye contact with the LED lights. The mask is designed with eye protection, but prolonged direct exposure should be avoided.
- Keep out of reach of children.
- Do not use on broken, irritated, or infected skin.
- Do not immerse the device in water or other liquids.
- Use only the provided charging cable.
- Discontinue use if you experience any discomfort or adverse reactions.

PRODUCT COMPONENTS

The foreverLily LED Facial Mask package includes the following items:

- LED Facial Mask
- USB Charging Cable
- Adjustable Head Strap
- User Manual (this document)



Image: Detailed view of the mask's interior, highlighting the 270 LED lamp beads, soft silicone pads for comfort, the VPS material construction, and the location of the charging port.



Image: The LED Facial Mask featuring rose gold metallic accents around the eye and nose areas, designed for full-face coverage with precise openings for eyes and mouth.

SETUP

1. **Charge the Device:** Before first use, fully charge the LED facial mask. Connect the USB charging cable to the mask's charging port and a suitable USB power adapter. The charging indicator light will show red during charging and turn green when fully charged. The device cannot be operated while charging.



Image: The LED mask shown from the side, with a charging cable connected to its port, illustrating the charging setup.

2. **Adjust the Head Strap:** Attach the adjustable head strap to secure the mask comfortably on your face. Ensure the mask fits snugly but without excessive pressure.



Image: Hands are shown adjusting the black fabric strap of the mask, illustrating how to achieve a secure and comfortable fit.

OPERATING INSTRUCTIONS

Follow these steps for optimal use of your LED Facial Mask:

1. **Prepare Your Face:** Cleanse your face thoroughly to remove any makeup, dirt, or oils. You may choose to apply a serum suitable for your skin type before using the mask.

4 Easy Steps



1. Wash your face, you can choose to apply serum



2. Put on the mask and choose the color



3. Lie flat with your eyes closed or playing phone, just enjoy 15 mins



4. Apply lotion

Image: A visual guide showing the four easy steps for using the mask: 1. Wash your face and apply serum. 2. Put on the mask and select a color. 3. Lie flat for 15 minutes. 4. Apply lotion after use.

2. **Power On and Select Light Mode:** Press and hold the power button (often a touch switch on the top of the mask) for 1.5 seconds to turn on the device. Briefly touch the switch to cycle through the available light colors.

Wireless Touch Screen

Touch the switch to change the color of you want



Image: A detailed view of the mask's top, indicating the "Touch Switch" for changing light colors and the "Boot/key" for power on/off. It also notes that a red light indicates charging, green for full charge, and the device cannot be opened while charging.




Image: A finger is shown pressing the touch switch on the top of the mask, with a small green indicator light illuminated, demonstrating how to operate the controls.

- Treatment Session:** Lie down comfortably with your eyes closed and wear the mask for approximately 15 minutes per session. The mask will automatically shut down after 15 minutes.



Image: A person is shown lying down comfortably while wearing the LED facial mask, illustrating the recommended posture during a treatment session.

4. **Post-Treatment:** After the session, remove the mask and apply your preferred moisturizer or lotion.



Image: A woman is applying a white moisturizing cream to her cheek, demonstrating the recommended skincare step after using the LED mask.

Understanding the LED Light Modes

The mask features multiple LED light colors, each designed to address specific skin concerns:

- **Red Light:** Promotes collagen production and helps reduce the appearance of fine lines and wrinkles.
- **Green Light:** Helps reduce blackheads and acne, and can improve overall skin tone.
- **Blue Light:** Soothes skin and targets acne-causing bacteria.
- **Purple Light:** A combination of red and blue light, beneficial for skin repair and reducing acne marks.
- **Orange Light:** Brightens skin tone and can help reduce dullness.
- **Bright Green Light:** Boosts skin metabolism.
- **White Light:** Helps firm the skin and improve overall texture.

7 Colors LED Lights

Different Lights Have Different Effects



Red light:



Green Light:



Blue Light:



Purple Light



Orange light:



Bright Green:



White Light:

Image: A visual representation of the seven LED light colors available on the mask, detailing the specific benefits of each color for skin care.




Image: A comprehensive chart outlining the benefits of each light color in phototherapy, such as red light for collagen and wrinkle reduction, green for blackheads and acne, blue for soothing, orange for brightening, bright green for metabolism, and white for skin firming.




Image: A side-by-side comparison of the skin around an eye, illustrating a reduction in fine lines and wrinkles in the "after" image compared to the "before" image.




Image: A before-and-after comparison of facial skin, with the "before" image showing redness and blemishes, and the "after" image displaying clearer, smoother skin.

MAINTENANCE

- **Cleaning:** After each use, gently wipe the mask with a soft, dry cloth. For deeper cleaning, use a slightly damp cloth with mild soap, then wipe dry immediately. Do not use abrasive cleaners or immerse the mask in water.
- **Storage:** Store the mask in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Mask does not turn on.	Low battery or not charged.	Ensure the mask is fully charged before use. The charging indicator should be green.
Lights are dim or flickering.	Low battery.	Recharge the mask.

Problem	Possible Cause	Solution
Mask feels uncomfortable.	Strap too tight or incorrect positioning.	Adjust the head strap for a comfortable fit. Ensure the mask is centered on your face.
Charging indicator not working.	Faulty cable or power source.	Try a different USB cable or power adapter. Ensure connections are secure.

SPECIFICATIONS

- **Brand:** foreverLily
- **Model:** B0DDBZLYT7
- **Product Dimensions:** 24 x 21.7 x 11.2 cm
- **Weight:** 650 g
- **Power Source:** Rechargeable Battery (1 x 9V type needed, likely internal)
- **Special Features:** Rechargeable, 7 LED Light Colors
- **Skin Type Compatibility:** All skin types
- **Benefits:** Anti-aging, Skin Firming, Acne Reduction, Skin Tone Improvement

WARRANTY AND SUPPORT

This foreverLily product is covered by a standard manufacturer's warranty. For specific warranty details, please refer to the documentation included with your purchase or contact the seller directly. For technical support, troubleshooting assistance, or any questions regarding your LED Facial Mask, please contact the seller through your purchase platform or visit the official foreverLily brand store on Amazon: [foreverLily Store](#).