

[Manuals.plus](#) /

› [Nicefeel](#) /

› Nicefeel Water Dental Flosser Replacement Classic Heads User Manual

Nicefeel Tip 2

User Manual

NICEFEEL WATER DENTAL FLOSSER REPLACEMENT CLASSIC HEADS

Model: Tip 2

Product Overview

These Nicefeel water dental flosser replacement tips are designed to maintain optimal oral hygiene. They deliver a powerful yet gentle stream of water to effectively remove plaque and debris from hard-to-reach areas between teeth and below the gumline. Made from high-quality materials, these tips ensure consistent performance and longevity.



Image: Two Nicefeel Water Dental Flosser Replacement Classic Heads and their packaging. The packaging is orange with the "Nicefeel" logo and text "2x Standard Tips Jet Tips for Water Flosser".

Setup and Compatibility

These replacement heads are compatible with most Nicefeel water dental flossers. They have been thoroughly tested to ensure a seamless fit with various models.

- **Compatible Models:** NF A1, FC159, FC2561, and FC5090.

Installation and Operation

Replacing the nozzle on your Nicefeel water dental flosser is a simple process. These classic nozzle tips are designed for deep cleaning, providing a thorough oral care experience.

1. **Press the Nozzle Replacement Button:** Locate and press the nozzle replacement button, typically found on the top

of your Nicefeel water dental flosser device.

- 2. Remove Old Nozzle:** Carefully pull out the old nozzle from its slot.
- 3. Insert New Nozzle:** Insert the new replacement tip into the slot. Ensure it is securely in place before operating the oral irrigator.

Once installed, the water flosser is ready for use. The design of these tips facilitates effective plaque removal and deep cleaning.

Maintenance

Regular replacement of your water dental flosser heads is crucial for maintaining peak performance and hygiene. This ensures your oral care routine remains effective.

- Replacement Frequency:** It is recommended to swap out the water dental flosser heads every 3-6 months. The exact frequency may vary depending on your usage habits.

Troubleshooting

If you encounter any issues with the installation or performance of your Nicefeel replacement classic heads, please refer to your water flosser's main instruction manual or contact Nicefeel customer support for assistance.

Specifications

Feature	Detail
Product Dimensions	2.15 x 1.02 x 4.72 inches
Item Weight	0.32 ounces
Model Number	Tip 2
Manufacturer	Nicefeel
Special Feature	Portable
Product Benefits	Gum Health, Flossing
Unit Count	2 Count

Support

For any questions or assistance regarding your Nicefeel Water Dental Flosser Replacement Classic Heads, please contact the Nicefeel support team. They are available to help with product-related inquiries.

Note: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.