



i02, i03, i04, i05, i06, i07, i08, i09, i010

Replacement Brush Heads User Manual

Compatible with Oral-B iO Series Electric Toothbrushes

INTRODUCTION

This manual provides essential information for the proper use, installation, and maintenance of your new replacement brush heads. These brush heads are designed to be compatible with specific Oral-B iO series electric toothbrushes, ensuring effective and hygienic oral care.



Figure 1: Overview of the Replacement Brush Head. This image displays the brush head's design, featuring a white body and a circular arrangement of blue and white bristles at the top, indicating its compatibility and function for electric toothbrushes.

SAFETY INFORMATION

Please read and understand all safety information before using the product.

- Keep out of reach of children.
- Do not use on sensitive areas of the mouth if irritation occurs.
- Dispose of properly after use, following local regulations for plastic waste.
- May contain materials that could cause an allergic reaction in sensitive individuals. Discontinue use if any adverse reaction is observed.

COMPATIBILITY

These replacement brush heads are specifically designed for compatibility with the following Oral-B iO series electric toothbrushes:

- Oral-B iO2
- Oral-B iO3

- Oral-B iO4
- Oral-B iO5
- Oral-B iO6
- Oral-B iO7
- Oral-B iO8
- Oral-B iO9
- Oral-B iO10

They are also compatible with Oral-B iO toothbrush types: **3753, 3758, 3769, 3776, 3794.**

SETUP AND INSTALLATION

Follow these steps to properly install a new brush head onto your Oral-B iO electric toothbrush handle:

1. Ensure your toothbrush handle is turned off.
2. Align the replacement brush head with the metal shaft on the toothbrush handle. The brush head should slide down easily onto the shaft.
3. Push the brush head firmly down until it clicks into place. Do not force it. A secure fit is indicated by a slight click and no wobble.
4. To remove a brush head, simply pull it straight up and off the handle.

OPERATING INSTRUCTIONS

Once the brush head is installed, you can use your electric toothbrush as usual. For optimal cleaning performance, follow these general guidelines:

- Wet the brush head and apply a small amount of toothpaste.
- Place the brush head against your teeth at a slight angle, covering one tooth at a time.
- Turn on the toothbrush. Allow the brush head to do the work, gently moving it from tooth to tooth. Do not scrub aggressively.
- Brush each quadrant of your mouth for approximately 30 seconds, or as recommended by your dentist.
- After brushing, rinse your mouth and the brush head thoroughly under running water.

MAINTENANCE AND CARE

Proper maintenance extends the life of your brush heads and ensures hygienic use.

- **Cleaning:** After each use, rinse the brush head thoroughly under running water with the toothbrush handle turned on briefly to remove any remaining toothpaste or debris. Then, turn off the handle and remove the brush head to rinse it more thoroughly.
- **Drying:** Allow the brush head to air dry in an upright position. Avoid covering it or storing it in a closed container, which can promote bacterial growth.
- **Replacement:** Dentists recommend replacing your toothbrush head every 3 months, or sooner if the bristles become frayed or worn. Worn bristles are less effective at cleaning and can harbor more bacteria.

TROUBLESHOOTING

If you encounter any issues with your replacement brush heads, consider the following:

- **Brush head does not fit:** Ensure you have an Oral-B iO series toothbrush handle. These brush heads are specifically designed for iO models and are not compatible with other Oral-B series (e.g., Vitality, Pro, Genius).
- **Bristles fraying quickly:** This may indicate excessive pressure during brushing. Electric toothbrushes are designed to clean effectively with minimal pressure.
- **Poor cleaning performance:** Check if the bristles are worn out. If so, it's time to replace the brush head. Ensure you are using proper brushing technique as described in the "Operating Instructions" section.
- **Unusual noise or vibration:** Ensure the brush head is securely attached to the handle. If the issue persists, consult your toothbrush handle's user manual or contact the manufacturer of your toothbrush handle.

SPECIFICATIONS

Feature	Detail
Brand	Generic
Compatibility	Oral-B iO2, iO3, iO4, iO5, iO6, iO7, iO8, iO9, iO10 series toothbrushes (Types 3753, 3758, 3769, 3776, 3794)
Quantity	4 Pack
Color	White
Material	Plastic
Item Weight	3 ounces (total for pack)
Manufacturer	Aftermarket

WARRANTY AND SUPPORT

For specific warranty information or product support, please refer to the packaging or contact the seller directly. As these are aftermarket replacement parts, warranty terms may differ from those of the original toothbrush manufacturer.

© 2024 Generic. All rights reserved.

This manual is for informational purposes only. Oral-B and iO are registered trademarks of their respective owners.