

- Manuals.plus /
- › LifePro /
- › Lifepro Red Light Therapy Belt User Manual

LifePro LP/ALVRD/BLK

Lifepro Red Light Therapy Belt User Manual

Model: LP/ALVRD/BLK | Brand: LifePro

INTRODUCTION

Welcome to your Lifepro Red Light Therapy Belt manual. This device provides professional-grade infrared and red light therapy for various body areas, including back, shoulders, and joints. It features dual LED settings for comprehensive wellness and relaxation.



The Lifepro Red Light Therapy Belt in use, demonstrating its flexible design and the visible red light LEDs.

SETUP GUIDE

The Lifepro Red Light Therapy Belt is designed for ease of use and portability.

ON-THE-GO RELIEF PORTABLE & ADJUSTABLE

Velcro and light weight materials provide
a perfect fit and a portable design.




Dimensions and weight of the portable and adjustable therapy belt.

The belt is lightweight and portable, weighing approximately 1.2 lbs (0.54 kg). It can be powered via a USB port from your computer, laptop, or a 5V-1A charging adapter.

UNPLUG & GO LIFEPRO POWERBANK

Pair with the Lifepro Powerbank for hassle-free portability and fewer wires for your wellness routine.

Powerbank sold separately.



Connecting the therapy belt to a power source for on-the-go relief.

Ensure all jewelry is removed from the area of application before use to prevent interference or discomfort.



FITS JUST RIGHT



EXTENSION BELT

The velcro extension belt stretches from 8 to 43 inches, providing easy comfort and targeted wellness.




The included extension belt provides a customizable fit for various body types.

An extension belt is included to ensure a comfortable and secure fit for various body sizes, stretching from 8 to 43 inches.

OPERATING YOUR THERAPY BELT

The belt features two distinct LED settings: 660nm Red Light and 850nm Near Infrared Light, offering dual-wavelength therapy for deep healing.

A man with grey hair and a beard, wearing a grey long-sleeved shirt and blue shorts, is holding a black and blue light therapy device. The device has a grid of red lights on the back. A red glow is visible on his lower back where the device is positioned. The background shows a living room with a grey sofa and a green plant.

DUAL-WAVELENGTH THERAPY FOR DEEP HEALING

ALL BENEFITS OF INFRARED & RED LIGHT

660 NM

Red Light



COMBINED MODE
for More Benefits




850 NM

Invisible Near
Infrared Light

Understanding the benefits of Red Light and Near Infrared Light wavelengths.

The 660nm Red Light is known for improving skin issues, accelerating wound healing, promoting blood circulation, stimulating collagen regeneration, and reducing skin wrinkles. The 850nm Near Infrared Light penetrates deeper tissues, aiding in muscle repair, joint stiffness relief, reducing inflammation, and boosting immunity.



CUSTOMIZE YOUR WELLNESS EASY-TO-USE REMOTE

Adjust time, intensity and
mode in one click.



The easy-to-use remote allows for simple adjustments of therapy settings.

Use the easy-to-use remote control to adjust the time, intensity, and mode of your therapy session.

ALLEVARED BELT EASE TENSION & DISCOMFORT

Let red light therapy alleviate tension
and help relax your muscles.



The hands-free design allows for comfortable use in various positions.

The hands-free design allows you to use the belt comfortably while sitting, lying down, or moving around, providing targeted relief where you need it most.

CARE AND MAINTENANCE

To ensure the longevity and optimal performance of your Lifepro Red Light Therapy Belt, proper care and maintenance are essential.

- Keep the device clean and dry.
- Wipe the surface with a soft, damp cloth. Do not use harsh chemicals or abrasive cleaners.
- Do not immerse the device in water or any other liquid.
- Store the belt in a cool, dry place, away from direct sunlight and extreme temperatures.

- When not in use, store the belt in its original packaging or a protective bag to prevent damage.

TROUBLESHOOTING COMMON ISSUES

If you encounter any issues with your Lifepro Red Light Therapy Belt, please refer to the following common troubleshooting steps:

- **Device Does Not Power On:** Ensure the belt is properly connected to a functional power source (USB port or 5V-1A adapter). Check the power indicator on the remote control.
- **Lights Not Illuminating:** Verify that the correct mode and intensity settings are selected on the remote control. Ensure the belt is securely plugged into the power source.
- **Discomfort During Use:** Ensure all jewelry is removed from the area of application. Adjust the belt for a comfortable fit. If discomfort persists, discontinue use and consult a healthcare professional.
- **Remote Control Not Responding:** Check the connection between the remote and the belt. Ensure there is no damage to the cable.

For any persistent issues not resolved by these steps, please contact Lifepro customer support or refer to the official user manual PDF.

PRODUCT SPECIFICATIONS

Feature	Detail
Item Package Dimensions L x W x H	11.14 x 8.7 x 3.62 inches
Package Weight	1.35 Kilograms
Item Dimensions LxWxH	49.6 x 7.08 x 1 inches
Brand Name	LifePro
Warranty Description	Lifetime Customer Support
Model Name	Red Light Therapy Belt
Color	Blue
Suggested Users	unisex-adult
Number of Items	1
Manufacturer	LifePro
Part Number	LP/ALVRD/BLK
Style	LifePro

Feature	Detail
Size	One Size
ASIN	B09CLM44TM
Customer Reviews	4.4 out of 5 stars (2,135 ratings)
Best Sellers Rank	#29,971 in Health & Household, #41 in Light Therapy Products
Date First Available	August 13, 2021

WARRANTY AND CUSTOMER SUPPORT

The Lifepro Red Light Therapy Belt comes with a lifetime warranty, ensuring long-term support for your product and peace of mind.

For any questions, concerns, or assistance with your device, please contact Lifepro customer service. Their dedicated team is available to provide support and guidance.

For detailed information regarding product usage, safety guidelines, and comprehensive troubleshooting, please refer to the [Official User Manual \(PDF\)](#).

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Answered questions about this manual

Why does the control panel on my Lifepro AllevaRed Light Therapy Belt (LP/ALVRD/BLK) turn on, but the LEDs on the belt do not illuminate?

If the control panel on your Lifepro AllevaRed Light Therapy Belt (LP/ALVRD/BLK) is functioning but the LEDs on the belt are not illuminating, please follow these troubleshooting steps: Check the Belt Connection: Ensure the...

[Read full answer](#)

The panel on my Lifepro AllevaRed Light Therapy Belt lights up, but why won't the belt LEDs illuminate?

If your Lifepro AllevaRed Light Therapy Belt (LP/ALVRD/BLK) has a functioning control panel but the LEDs on the belt are not illuminating, please follow these troubleshooting steps: Check the Belt Connection: Ensure the male...

[Read full answer](#)

Open the live page

<https://manuals.plus/asin/B0DNG85YT8>

More documents and user manuals: <https://manuals.plus/>