

CULT MAS0015

CULT Impact Deep Tissue Massage Gun User Manual

Model: MAS0015 | Brand: CULT

1. INTRODUCTION

Thank you for choosing the CULT Impact Deep Tissue Massage Gun. This device is designed to provide effective muscle relief and recovery through percussion massage. Featuring a powerful motor, long-lasting battery, and interchangeable massage heads, it offers a professional-grade massage experience at home or on the go. Please read this manual thoroughly before use to ensure safe and optimal operation.

2. SAFETY INFORMATION

To reduce the risk of electric shock, fire, and injury, please observe the following precautions:

- Use the device only as described in this manual.
- Keep the device clean and dry. Do not submerge it in water or expose it to moisture.
- Do not operate continuously for more than one hour. Allow 30 minutes of rest before reusing.
- Do not use on injured skin, open wounds, or areas with swelling, bruising, or inflammation.
- Consult a physician before use if you have any medical conditions, are pregnant, or have a pacemaker.
- Keep out of reach of children and pets.
- Do not block the ventilation ports.
- Only use the provided charging cable.

3. PACKAGE CONTENTS

Your CULT Impact Deep Tissue Massage Gun package includes:

- CULT Impact Deep Tissue Massage Gun
- 4 Interchangeable Massage Heads (Silicone Round, Flat, Cone, Spade)
- Type-C Charging Cable
- User Manual

4. PRODUCT OVERVIEW

The CULT Impact Massage Gun is designed for ergonomic handling and powerful performance. It features a durable body, an intuitive control panel, and a robust motor for deep tissue percussion.



Figure 1: Key Features of the CULT Impact Massage Gun

SOFT SILICONE ROUND HEAD



Anti-Slip Grip



Allergy Friendly



Durable Design



Body Adaptive



Figure 2: Soft Silicone Round Head Features

POWERFUL & DURABLE MOTOR

Powered by a **30W motor** for smooth, quiet,
and long-lasting performance.



Figure 3: Powerful 30W Copper Brush Motor

3 HOURS LONG BATTERY LIFE

Powered by a **3000 mAh** battery for
uninterrupted relaxation



Figure 4: 3000mAh Long-Lasting Battery

5. SETUP

5.1 Charging the Device

Before first use, fully charge the massage gun. Connect the Type-C charging cable to the charging port on the device and plug the other end into a compatible USB power adapter (not included). The LED indicators will show the charging status. A full charge takes approximately 2 hours and provides up to 4 hours of continuous use.

5.2 Attaching Massage Heads

To attach a massage head, align the head with the opening on the front of the gun and push it firmly into place. To remove, gently pull the head out. Ensure the device is turned off before changing heads.

6. OPERATING INSTRUCTIONS

Watch this video for a visual guide on operating your CULT Impact Deep Tissue Massage Gun:

Your browser does not support the video tag.

Video 1: Unboxing and Basic Operation of the CULT Impact Massage Gun

This video demonstrates the unboxing, components, and basic operation of the CULT Impact Deep Tissue Massage Gun, including how to attach different massage heads and adjust speed settings.

6.1 Power On/Off

To turn on the device, long-press the power button located on the top of the gun. The LED indicators will light up. To turn off, long-press the power button again.


6.2 Speed Adjustment

The CULT Impact Massage Gun offers 6 intensity modes. Once powered on, short-press the power button to cycle through the speed settings. The LED indicators will show the current speed level.

6.3 Using Different Massage Heads

Each massage head is designed for specific muscle groups and massage types:

- **Silicone Round Head:** Ideal for large muscle groups like arms, back, hips, thighs, and calves. It provides a comfortable and gentle massage, soothing tension.




Silicone Round Head

Soothes tension in large muscle groups.

Usage: Silicone Round Head

- **Flat Head:** Suitable for relieving stiffness in chest, abs, and leg muscles. Its flat surface provides a broader contact area for general relaxation.



Flat Head

Relieves stiffness in chest, abs, and legs.

Usage: Flat Head

- **Cone Head:** Designed to target soft tissues and specific trigger points, promoting deeper relaxation. It mimics a finger massage for precise application.

Cone Head

Targets soft tissues to promote relaxation.



Usage: Cone Head

- **Spade Head:** Best for soothing neck and shoulder muscle discomfort. Its unique shape allows for effective massage along the spine and shoulder blades.

Spade Head

Soothes neck and shoulder muscle discomfort.



Usage: Spade Head




Figure 5: Targeted Application on Various Muscle Groups

7. USAGE GUIDELINES

Gently glide the massage gun over dry, clothed soft tissue. Massage any one muscle area for up to 60 seconds. Avoid applying excessive pressure or massaging bony areas. Listen to your body and adjust the intensity as needed.

TARGETED RELIEF

AT YOUR FINGERTIPS



Figure 6: Targeted Relief Application

8. MAINTENANCE

8.1 Cleaning

Wipe the device with a slightly damp cloth after each use. Do not use abrasive cleaners or solvents. Ensure the device is completely dry before storing.

8.2 Storage

Store the massage gun in a cool, dry place away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.

9. TROUBLESHOOTING

- **Device not turning on:** Ensure the battery is charged. Long-press the power button to activate.
- **Weak percussion:** Check the battery level. If low, recharge the device.
- **Unusual noise:** Ensure the massage head is securely attached. If the issue persists, discontinue use and contact customer support.

- **Device stops during use:** The device may have overheated or the battery is low. Allow it to cool down or recharge.

10. SPECIFICATIONS

Feature	Detail
Model Number	MAS0015
Dimensions (LxWxH)	21 x 6.3 x 20 cm
Item Weight	390 g
Power Source	Battery Powered
Battery Capacity	3000mAh Lithium-ion
Runtime	Up to 4 hours
Charging Time	Under 2 hours (Type-C)
Motor Power	30W
Force	Up to 6.5kg
Amplitude	4mm
Frequency (RPM)	Up to 3200 RPM
Speed Settings	6
Material	Acrylonitrile Butadiene Styrene
Color	Grey
Country of Origin	China

11. WARRANTY AND SUPPORT

Your CULT Impact Deep Tissue Massage Gun is backed by a 1-year warranty. For warranty claims, technical support, or any inquiries, please refer to the contact information provided on your warranty card or visit the official CULT website. Keep your purchase receipt for warranty validation.