

[Manuals.plus](#) /

› [RENPHO](#) /

› RENPHO ELIS 3X Smart Body Composition Scale User Manual

RENPHO ELIS 3X

RENPHO ELIS 3X Smart Body Composition Scale User Manual

Model: ELIS 3X

1. INTRODUCTION

This manual provides comprehensive instructions for the safe and effective use of your RENPHO ELIS 3X Smart Body Composition Scale. This device is designed to measure body weight and analyze various body composition metrics. Please read this manual thoroughly before first use and retain it for future reference.



Image: RENPHO ELIS 3X Smart Scale highlighting its user-friendly dial operation and display metrics.

2. PACKAGE CONTENTS

- RENPHO ELIS 3X Smart Body Composition Scale
- User Manual

3. SETUP

3.1 Initial Charging

The RENPHO ELIS 3X scale is equipped with a rechargeable battery. Before first use, fully charge the device using the provided USB-C cable.

- A full charge takes approximately 3 hours.
- A single charge can provide up to 180 days of use.


USB-C Rechargeable – Never Buy Batteries Again



3-Hour
Charging



180-Day
Battery Life



USB-C
Rechargeable




Image: The RENPHO ELIS 3X Smart Scale showing its USB-C charging port and overall dimensions.

3.2 Placement

For accurate measurements, place the scale on a hard, flat surface. Avoid carpets or uneven flooring, as this can affect accuracy. The scale features anti-slip foot pads for stability.

3.3 Unit Conversion

The scale supports various units for weight (kg, lb, st) and height (cm, inch). You can typically change these settings via the dial on the scale or within the RENPHO Health App.



Image: Close-up of the RENPHO ELIS 3X Smart Scale illustrating the dial for profile editing and a button for unit conversion.

4. OPERATING INSTRUCTIONS

4.1 Creating a User Profile (On-Scale)

The ELIS 3X scale allows you to create and manage user profiles directly on the device using the integrated dial.

1. Turn the dial to select "USER 1-8" or "GUEST" and press the dial to enter profile settings.
2. Turn the dial to select gender and press to confirm.
3. Turn the dial to select age and press to confirm.
4. Turn the dial to select height and press to confirm.

- ① Dial
- ② ITO-Coated Full Coverage Sensor
- ③ Tempered Glass
- ④ BIA Chip
- ⑤ Seamless Wrapped Base
- ⑥ 4 Anti-Slip Foot Pads



Image: Visual guide demonstrating the four steps to set up a user profile directly on the scale's display using the control dial.

4.2 Taking a Measurement

Step onto the scale with bare, dry feet. Ensure your feet are positioned correctly on the ITO-coated full-coverage sensor area for accurate body composition analysis. The scale will display your metrics.

ITO-Coated Full-Coverage Sensor for 99.9% Data

Maximum Capacity: 400 lbs




Image: The RENPHO ELIS 3X Smart Scale highlighting its ITO-coated full-coverage sensor for accurate data collection.

4.3 Displayed Metrics

The large VA display on the scale provides instant insights into 6 key health metrics:

- Weight
- Body Fat Percentage

- Muscle Percentage
- Visceral Fat Index
- BMI (Body Mass Index)
- 20-day Weight Trend

All-in-One Health Monitor

See 6 metrics and 20-day weight trend directly on the display



Image: The RENPHO ELIS 3X Smart Scale's display showing weight, body fat, muscle, visceral fat, BMI, and a 20-day weight trend graph.


4.4 Optional App Connection (RENPHO Health App)

For a more detailed analysis of 18 body metrics and long-term tracking, connect your scale to the RENPHO Health App.

1. Download the RENPHO Health App from your device's app store.
2. Register and add users within the app.
3. Step on the scale and connect it to your smartphone via Bluetooth.
4. Match your profile in the app and set your goals.
5. The scale will display your user ID and progress towards goals.

Connect to RENPHO Health App and Track 18 Body Metrics

Optional App Connection



5 Body Goal Metrics



Body Type



Target to optimal fat mass



Optimal Weight



Target to optimal weight



Target to optimal muscle mass

Work with

RENPHO



Image: A user interacting with the RENPHO Health App on a smartphone, displaying various body metrics tracked by the scale.

Auto-Identify User by Weight




Image: A visual representation of the four steps required to sync the smart scale with the RENPHO Health App for data tracking.

4.5 Multi-User Functionality

The scale can store up to 8 user profiles. It can auto-identify users by weight or allow manual selection via the dial.

- **Auto-Identify:** The scale attempts to match the current measurement to an existing profile based on weight.
- **Manual Selection:** If auto-identification is not desired or if weights are similar, use the dial to select your profile before weighing.

One Scale for Whole Family

Store up to 8 user profiles on the scale





Image: A family using the RENPHO ELIS 3X Smart Scale, illustrating its capability to store and track data for multiple users.

How to Set Up a User Profile with Dial?



Step 1

Turn the dial to select "USER 1-8" or "GUEST" and press the dial to enter the profile settings.



Step 2

Turn the dial to select gender and press to confirm.



Step 3

Turn the dial to select age and press to confirm.



Step 4

Turn the dial to select height and press to confirm.

Image: The RENPHO ELIS 3X Smart Scale demonstrating its feature to automatically identify different users based on their weight.

Sync with App




Step 1

Download RENPHO Health app, register, and add users.




Step 2

Step on the scale and connect to smartphone.



Step 3

Match your profile and set your goals.



Step 4

Now, see your name & goals on scale.

Image: The RENPHO ELIS 3X Smart Scale illustrating the process of manually selecting a user profile using the integrated dial.

4.6 Pregnancy Mode

The scale includes a Pregnancy Mode, designed for expectant mothers and individuals with implanted medical devices. This mode provides weight tracking without bio-electrical impedance analysis (BIA), ensuring safety.



Real-Time Guidance for Perfect Measurements

Image: A pregnant woman near the RENPHO ELIS 3X Smart Scale, with an inset showing the Pregnancy Mode setting enabled in the app.

5. MAINTENANCE

5.1 Cleaning

Wipe the scale's surface with a damp cloth and mild detergent. Avoid abrasive cleaners or submerging the scale in water. Ensure the scale is dry before storage or next use.

5.2 Storage

Store the scale in a cool, dry place away from direct sunlight and extreme temperatures.

6. TROUBLESHOOTING

- **Scale not turning on:** Ensure the scale is fully charged. Connect it to a USB-C power source for at least 3 hours.
- **Inaccurate readings:**
 - Ensure the scale is placed on a hard, flat surface.
 - Step onto the scale with bare, dry feet.
 - Avoid moving during measurement.
 - Recalibrate the scale by stepping on it briefly, then stepping off, and waiting for it to display '0.0'.
- **App not syncing:**
 - Ensure Bluetooth is enabled on your smartphone.
 - Make sure the RENPHO Health App is open and running in the foreground.
 - Check if the scale is within Bluetooth range of your smartphone.
 - Restart both the scale and your smartphone.
- **Body composition metrics not displayed:** Ensure you are stepping on the scale with bare, dry feet, making full contact with the ITO-coated sensor area.

7. SPECIFICATIONS

Feature	Specification
Brand	RENPHO
Model	ELIS 3X
Display Type	Digital VA Display
Weight Limit	400 Pounds (180 kg)
Power Source	Rechargeable (USB-C)
Battery Life	Up to 180 days on a single charge
Body Composition Metrics	6 on-scale, 18 via app
User Profiles	Up to 8
Dimensions	11.02 x 11.02 x 1.26 inches (approx.)
Color	Dark Blue



Image: An exploded diagram showing the internal components and layers of the RENPHO ELIS 3X Smart Scale, including the dial, ITO sensor, tempered glass, BIA chip, base, and anti-slip foot pads.

8. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the official RENPHO website or contact their customer service directly. Details are typically available on the product packaging or within the RENPHO Health App.
Manufacturer: RENPHO

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0DQLDXGJ7>

More documents and user manuals: <https://manuals.plus/>