

[Manuals.plus](#) /

› [CINCOM](#) /

› CINCOM Leg Massager with Heat and Air Compression User Manual

CINCOM CM-121A

CINCOM Leg Massager with Heat and Air Compression User Manual

Model: CM-121A

1. PRODUCT OVERVIEW

The CINCOM Cordless Leg Massager is designed to provide soothing heat therapy combined with air compression pulses for comprehensive leg massage. It is ideal for improving circulation, relieving pain, and reducing edema in the calves. Its portable and rechargeable design allows for use anywhere, from post-workout recovery to relaxation at home or in the office.



Image: The CINCOM Leg Massager, showing the two main units.

2. WHAT'S INCLUDED

Your CINCOM Leg Massager package includes the following components:

- Two Massager Wraps
- User Manual
- Type-C Charging Cable and Adapter
- Storage Bag

CINCOM Provide



Two Massager Warps



User Manual



Type-c



Storage Bag

Image: All items included in the CINCOM Leg Massager package.

3. SAFETY INFORMATION

Please read all safety warnings and instructions before using this product to prevent injury or damage.

- Avoid using for more than 20 minutes at a time.
- Consult a healthcare professional if you have certain medical conditions before use.
- Exercise caution when using heat therapy, especially if you have sensitive skin or circulatory issues.

- Take precautions to prevent electrical hazards; do not use with damaged cords or in wet environments.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

4. SETUP

4.1 Wearing the Massager

The massager is designed to fit calves up to 24.5 inches. Adjust the size using the Velcro fasteners to ensure a comfortable and secure fit. It can be used on both upper and lower legs, as well as arms and forearms for versatile relief.

Large Adjustable Size Up to **24.5 in**



Image: Woman demonstrating the adjustable size of the leg massager.

4.2 Charging the Device

The CINCOM Leg Massager features a 2500mAh battery, providing extended use on a single charge. To charge, connect the included Type-C charging cable to the device and a power source. The indicator lights on the control panel will show the charging status. A full charge can last for approximately a week of typical use.

Cordless & Portable Design



2500mAh
Reliable Performance



3 Heating
104°F / 113°F / 122°F



Intensity
Low / Mid / High



Power Button
20 Minutes Shut-Off



3 Modes
Sequence / Circulation / Combination



Image: Close-up of the control panel showing power, heat, intensity, and mode indicators.

5. OPERATING INSTRUCTIONS

5.1 Powering On/Off

Press and hold the power button on the control panel to turn the massager on or off.

5.2 Adjusting Settings

The massager offers 3 heat levels, 3 massage modes, and 3 intensities, allowing for 27 possible combinations to customize your massage experience.

- **Heat Levels:** Select from 104°F, 113°F, or 122°F to apply soothing warmth.
- **Massage Modes:** Choose between Sequence, Circulation, or Combination modes for varied massage patterns.
- **Intensities:** Adjust the compression strength to Low, Mid, or High.

Heated Massager with Air, Anywhere



Wireless



Air Compression




3 Modes for Relief




3 Heating Levels

Image: Woman relaxing with the massager, highlighting wireless, air compression, 3 modes, and 3 heating levels.

Legs × Heat = ∞ Comfort



Heating area



Quick heating




Image: Diagram illustrating the heating area and quick heating temperature levels.

5.3 Usage Recommendations

The massager is equipped with a 20-minute smart timer that automatically shuts off the device after each session, ensuring optimal muscle relaxation and preventing overuse. It can be used for post-exercise recovery, managing swelling, or simply for daily relaxation.

Post-exercise Recovery





Image: Man demonstrating the massager's use on arms, calves, thighs, and forearms.

Anytime, Anywhere

Enjoy CINCOM Wireless Calf Massager




Image: People using the massager at home, while working, and outdoors, emphasizing its portability.

6. MAINTENANCE

To ensure the longevity and optimal performance of your CINCOM Leg Massager, follow these maintenance guidelines:

- Wipe the surface of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- Store the massager in its included storage bag in a cool, dry place away from direct sunlight and extreme

temperatures.

- Ensure the device is fully charged before long-term storage to maintain battery health.

7. TROUBLESHOOTING

If you encounter any issues with your CINCOM Leg Massager, please refer to the following common solutions:

- **Device not turning on:** Ensure the battery is fully charged. Connect the massager to the charging cable and allow it to charge for at least 2-3 hours.
- **Weak compression:** Check that the massager wraps are securely fastened around your legs. Adjust the intensity setting to a higher level.
- **No heat:** Verify that the heat function is activated and set to a desired temperature. Allow a few minutes for the heat to build up.
- **Uncomfortable fit:** Readjust the Velcro straps to achieve a snug yet comfortable fit. Ensure the massager is not too tight.

If the problem persists, please contact CINCOM customer support for assistance.

8. SPECIFICATIONS

Feature	Detail
Product Dimensions	8.2 x 3.4 x 14.8 inches
Weight	3.31 Pounds
Model Number	CM-121A
Power Source	Battery Powered
Material	Nylon
Color	Gray
Battery Capacity	2500mAh
Heat Levels	3 (104°F, 113°F, 122°F)
Massage Modes	3 (Sequence, Circulation, Combination)
Intensity Levels	3 (Low, Mid, High)
Auto Shut-off	20 minutes

9. WARRANTY AND SUPPORT

CINCOM is a reputable seller of health massagers. For any inquiries, technical support, or warranty claims, please refer to the contact information provided in your product packaging or visit the official CINCOM website. Our dedicated customer service team is ready to assist you.