

Vive MOB1117

Vive Mobility Lightweight Folding Walker with Wheels User Manual

Model: MOB1117

INTRODUCTION

The Vive Mobility Lightweight Folding Walker with Wheels is designed to provide stable and effortless mobility for individuals, particularly seniors, navigating various environments. Its compact, narrow frame and ultra-lightweight design make it ideal for use in small spaces, crowded areas, and for convenient transport and storage. This manual provides essential information for the safe and effective use, setup, and maintenance of your walker.

KEY FEATURES

- **Compact, Foldable, & Ultra-Lightweight:** Effortlessly folds in one step for easy transport and storage. Weighs just 9.5 pounds yet supports up to 300 pounds. Features a patent-pending locking system to keep it securely folded.
- **Compact, Narrow Frame:** Ideal for small hallways, crowded areas, and sharp turns. The open-frame design eliminates obstructive rear bars, allowing for a natural walking stride.
- **Adjustable Height & Comfort Grip Brakes:** Customize fit with height-adjustable handles (36" - 42") for optimal posture. Ergonomic hand grips provide a secure, slip-resistant hold, while hand brakes ensure smooth stopping and control.
- **Durable All-Terrain Wheels:** Durable 6-inch wheels glide smoothly over various surfaces including carpets, tile, pavement, and uneven terrain.
- **Pre-Assembled:** Arrives fully assembled; simply unfold and use.

SETUP

Unfolding the Walker

The Vive Mobility Lightweight Folding Walker arrives pre-assembled. To prepare it for use:

1. Carefully remove the walker from its packaging.
2. Place the folded walker on a flat, stable surface.
3. Grasp both handles and pull the frame apart until it fully extends and locks into place. You should hear a click indicating it is securely open.



COMPACT FOLDING DESIGN

Folds flat for storage and travel convenience



BUILT FOR ALL DAY USE

Made from durable aluminium



ULTRA-LIGHT CONSTRUCTION

Weight only 9.5 lbs, easy to lift and carry



DURABLE SUPPORT

Holds up to 300 lbs



PRE-ASSEMBLED AND READY TO USE

Enhanced Walking Space

Figure 1: The Vive Mobility Lightweight Folding Walker in its fully unfolded and ready-to-use state.

Adjusting Handle Height

For optimal posture and comfort, adjust the handle height to fit your body:

1. Locate the height adjustment lock on each front upright.
2. Loosen the lock by unscrewing the star knob.
3. Move the handle to your desired height. The handles should be at wrist level when your arms are hanging loosely at your sides.
4. Once the desired height is achieved, securely tighten the star knob to lock the handle in place. Ensure both handles are at the same height.

Ultra Lightweight

Weighs only 9.5 lbs



Figure 2: The walker features adjustable handles with a height range of 33 to 39 inches for personalized comfort.

Video: About the Vive 4-Wheel AirFrame Walker

Your browser does not support the video tag.

Video 1: This video provides an overview of the Vive 4-Wheel AirFrame Walker, demonstrating its features and ease of use.

OPERATING INSTRUCTIONS

Walking with the Walker

To use the walker for mobility:

1. Ensure the walker is fully unfolded and the handles are adjusted to your comfortable height.
2. Grip the ergonomic handles with a secure, slip-resistant hold.
3. Push the walker forward a short distance.
4. Step forward, maintaining your balance and allowing the open-frame design to provide unobstructed walking space.
5. Repeat the process, moving the walker and then stepping forward.

Open-Frame Design

Clear path for smooth, comfortable walking

OTHER BRANDS



RISK OF KNEE BUMPS

VIVE WALKER



NO KNEE-BUMPING BAR

Figure 3: The open-frame design allows for a natural walking motion without obstruction.

Using the Brakes

The walker is equipped with hand brakes for control and stopping:


- To slow down or stop, gently squeeze the brake levers located on the handles.
- To engage the parking brake (lock the wheels), push the brake levers down until they click into the locked position. This is useful when standing or resting.
- To release the parking brake, pull the brake levers up.

Designed for All-Day Comfort

Soft-grip handles reduce hand fatigue



**SOFT PADDING FOR
EXTENDED-USE COMFORT**



**ARTHRITIS & WEAK GRIP
STRENGTH FRIENDLY**



Figure 4: Ergonomic hand grips and easy-to-use brake levers provide comfort and control.

Folding for Storage and Transport

The walker features a quick and easy one-step folding mechanism:

1. Ensure the parking brakes are disengaged.
2. Locate the folding lever, typically near the center of the frame.
3. Pull the lever upwards to release the locking mechanism.
4. Gently push the sides of the walker together until it is fully folded. The patent-pending locking system will secure it in the folded position.

Foldable, Compact, & Portable

Travel-friendly design with easy storage



Figure 5: The walker folds compactly for easy storage and transport in vehicles or small spaces.

Video: Top Features | Vive Lightweight Walker

Your browser does not support the video tag.

Video 2: This video highlights the top features of the Vive Lightweight Walker, including its folding mechanism and maneuverability.

MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your walker:

- **Cleaning:** Wipe down the frame and handles with a damp cloth and mild soap. Avoid using abrasive cleaners or solvents that could damage the finish.
- **Wheel Inspection:** Periodically check the 6-inch wheels for wear and tear. Ensure they rotate freely and are free from debris.
- **Brake Check:** Regularly test the hand brakes to ensure they engage and disengage properly, providing reliable stopping power.
- **Frame Integrity:** Inspect the aluminum frame for any signs of damage, cracks, or loose components. Do not use the walker

if any part is compromised.

- **Storage:** Store the walker in a dry place when not in use to prevent rust or damage from moisture.

TROUBLESHOOTING

If you encounter any issues with your Vive Mobility Lightweight Folding Walker, refer to the following common solutions:

- **Walker does not unfold/fold smoothly:** Ensure all locking mechanisms are fully disengaged before attempting to fold or unfold. Check for any obstructions in the folding mechanism.
- **Wheels are not rolling smoothly:** Inspect the wheels for any trapped debris (hair, lint, etc.) and clean them. Ensure the wheel axles are not bent or damaged.
- **Brakes are not engaging/disengaging:** Check the brake cables for any kinks or damage. Ensure the brake levers are properly adjusted and not stuck. If issues persist, contact customer support.
- **Handles are loose after adjustment:** Re-tighten the star knobs on the height adjustment locks firmly.

For issues not covered here, please contact Vive customer support.

SPECIFICATIONS

Feature	Specification
Product Dimensions (D x W x H)	21"D x 26"W x 42"H
Item Weight	9.5 Pounds
Maximum Height	39 Inches
Minimum Height	33 Inches
Weight Capacity Maximum	300 Pounds
Color	Midnight Blue
Material Type	Aluminum
Is Foldable	Yes
Walker Type	Standard Walker
UPC	810113579771

WARRANTY AND SUPPORT

The Vive Mobility Lightweight Folding Walker comes with a **LIFETIME GUARANTEE** on its frame, providing lasting confidence in its durability. For any questions, concerns, or warranty claims, please contact Vive Health Direct customer support. This product is also eligible for purchase with Health Savings Account (HSA) or Flexible Spending Account (FSA) funds.

Vive Lifetime Guarantee



Figure 6: The Vive Health Lifetime Guarantee ensures product quality and customer satisfaction.

© 2026 Vive Health. All rights reserved.
For further assistance, visit vivehealth.com

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question