

[Manuals.plus](#) /

› [Vive](#) /

› Vive Mobility Lightweight Walker User Manual

Vive B0DT7NZDSZ

Vive Mobility Lightweight Walker User Manual

Model: B0DT7NZDSZ

1. INTRODUCTION AND OVERVIEW

Thank you for choosing the Vive Mobility Lightweight Walker. This manual provides essential information for the safe and effective use of your new walker. Designed for effortless handling, this walker is compact, foldable, and ultra-lightweight, making it ideal for seniors and individuals seeking enhanced mobility.

Key features include a patent-pending locking system for secure folding, a narrow frame perfect for tight spaces, an open-frame design for a natural walking motion, adjustable height handles, and comfort grip brakes. Constructed from high-strength aluminum, it supports up to 300 pounds while weighing just 9.5 pounds.



Image: The Vive Mobility Lightweight Walker, showcasing its sleek design and lightweight nature (9.5 lbs).



COMPACT FOLDING DESIGN

Folds flat for storage and travel convenience



BUILT FOR ALL DAY USE

Made from durable aluminium



ULTRA-LIGHT CONSTRUCTION

Weight only 9.5 lbs, easy to lift and carry



DURABLE SUPPORT

Holds up to 300 lbs



PRE-ASSEMBLED AND READY TO USE

Enhanced Walking Space

Image: Diagram highlighting key features: compact folding, durable aluminum, ultra-light construction (9.5 lbs), 300 lbs support, and pre-assembled design.

Sturdy & Long Lasting

Supports up to 300 lbs



Image: A man using the walker outdoors, emphasizing its sturdy and long-lasting construction, supported by a lifetime guarantee and patent-pending locking system.

2. SETUP AND ASSEMBLY

The Vive Mobility Lightweight Walker arrives fully assembled, requiring minimal setup before use. Simply unfold the walker and ensure all locking mechanisms are securely engaged.

2.1 Unfolding the Walker

To unfold the walker, gently pull the sides apart until the frame clicks into its fully open and locked position. Verify that the patent-pending locking system is engaged to ensure stability.

2.2 Adjusting Handle Height

For optimal posture and comfort, adjust the height of the handles to fit your body. Locate the height adjustment locks on each handle. Loosen the lock, slide the handle to your desired height (between 33 and 39 inches), and then firmly relock it. Ensure both handles are set to the same height for balanced support.

Foldable, Compact, & Portable

Travel-friendly design with easy storage



Image: A woman demonstrating the compact folding of the walker, easily fitting it into a car trunk for travel and storage.

Ultra Lightweight

Weighs only 9.5 lbs



Image: Visual representation of the walker's dimensions (21"D x 26"W x 42"H) and adjustable handle height range (33"-39"), highlighting its ultra-lightweight design.

3. OPERATING INSTRUCTIONS

The Vive Mobility Lightweight Walker is designed for smooth and natural movement across various surfaces.

3.1 Walking with the Walker

Grip the ergonomic handles with a secure, slip-resistant hold. The open-frame design eliminates obstructive rear bars, allowing you to walk naturally without interference. Move effortlessly with the durable 6-inch wheels that glide smoothly over carpets, tile, pavement, and uneven surfaces, making it suitable for both indoor and outdoor use.

3.2 Using the Brakes

The walker features hand brakes for controlled stopping and speed management. Squeeze the hand brakes to slow down or stop. For parking, push the brake levers down until they lock into place, securing the walker. To release the parking brake, pull the levers back up.

3.3 Maneuvering in Tight Spaces

The narrow, space-saving frame and 6.5-inch rotating wheels make this walker ideal for navigating small hallways, crowded areas, and sharp turns with ease.

Open-Frame Design

Clear path for smooth, comfortable walking

OTHER BRANDS



RISK OF KNEE BUMPS

VIVE WALKER



NO KNEE-BUMPING BAR

Image: A visual comparison demonstrating the open-frame design of the Vive Walker, which provides a clear path for walking without knee-bumping bars, unlike other traditional walkers.

Easy to Maneuver

Navigate tight spaces and turns



6.5" Rotating Wheels

Maneuver through sharp turns and rough terrains with ease and comfort



Narrow Design

Move effortlessly through narrow halls and doorways

Image: Demonstrating the walker's ease of maneuverability, showing it navigating sharp turns outdoors and moving effortlessly through narrow doorways indoors.

Designed for All-Day Comfort

Soft-grip handles reduce hand fatigue




Image: Detailed view of the soft-grip handles and hand brakes, designed for extended-use comfort and friendly for users with arthritis or weak grip strength.

Your browser does not support the video tag.

Video: An official product video from Vive Health Direct, titled "About the 4 Wheel AirFrame Walker," demonstrating the walker's features, folding mechanism, and ease of use in various environments.

4. MAINTENANCE

To ensure the longevity and optimal performance of your Vive Mobility Lightweight Walker, regular maintenance is recommended:

- **Cleaning:** Wipe down the frame and handles with a damp cloth and mild soap as needed. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the wheels for any debris, wear, or damage. Ensure they rotate freely.
- **Brake Check:** Test the hand brakes regularly to ensure they engage and disengage properly, providing reliable stopping power.
- **Frame Integrity:** Inspect the aluminum frame for any signs of bending, cracks, or loose connections. Ensure all screws and fasteners are tight.

5. TROUBLESHOOTING

If you encounter any issues with your Vive Mobility Lightweight Walker, please refer to the following common concerns:

- **Walker does not fold/unfold smoothly:** Ensure no debris is obstructing the folding mechanism. Check that the locking system is fully disengaged before folding and fully engaged after unfolding.
- **Wheels are not rolling smoothly:** Inspect wheels for hair, lint, or other obstructions. Clean the wheel axles if necessary.
- **Brakes are not engaging/disengaging:** Check the brake cables for any kinks or damage. Ensure the brake levers are not stuck. If issues persist, contact customer support.
- **Walker feels unstable:** Verify that the handle heights are evenly adjusted and securely locked. Ensure the frame is fully unfolded and locked.

For any issues not covered here, or if problems persist, please contact Vive Health customer support.

6. SPECIFICATIONS

Brand	Vive
Model	B0DT7NZDSZ
Color	Midnight Blue
Product Dimensions	21"D x 26"W x 42"H (21 x 26 x 42 inches)
Item Weight	9.5 Pounds
Maximum Weight Capacity	300 lbs
Adjustable Handle Height Range	33" - 39" (Maximum Height: 39 Inches)
Wheel Size	6 inches
Manufacturer	Vive Health
Country of Origin	China

7. WARRANTY AND SUPPORT

The Vive Mobility Lightweight Walker comes with a **LIFETIME GUARANTEE** on its frame, providing lasting confidence in its durability. This product is also eligible for purchase with Health Savings Account (HSA) or Flexible Spending Account (FSA) funds.

For product support, warranty claims, or any questions, please contact Vive Health Direct, the seller and manufacturer of this product. Refer to your purchase documentation for specific contact details or visit the Vive Health official website.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question

and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

What is the purpose of the red lever on the brake handle of the Vive Mobility Lightweight Walker?

On the Vive Mobility Lightweight Walker, the levers on the brake handles serve two primary purposes depending on how they are moved:Speed Control and Stopping: Squeezing the levers upward toward the handle allows you...

[Read full answer](#)