

BEBIRD R37R/D39R/B31R/R1Ultra

Bebird Earsight Plus Ear Cleaner Replacement Tips Instruction Manual

Official Guide for Models R37R, D39R, B31R, and R1 Ultra

1. PRODUCT OVERVIEW

The Bebird Original Earsight Plus Ear Cleaner Replacement Tips are genuine and safe accessories designed for use with specific Bebird ear cleaning devices. These tips are crafted from medical-grade materials to ensure comfort, durability, and safe ear cleaning. The kit provides a comprehensive set of tools for various earwax removal needs.

This accessory kit is designed exclusively for the following Bebird Earsight Plus models: **R37R**, **D39R**, and **B31R**. It is also compatible with the **BebirdR1 Ultra**. Please note that these tips are *not* compatible with Bebird R1 or other Bebird devices not listed.



Image 1.1: Compatible Bebird Earsight Plus devices, including models R37R, D39R, B31R, and R1 Ultra, for which these replacement tips are designed.

2. KIT COMPONENTS

This 14-piece accessory kit includes a variety of specialized cleaning tips to address different types of earwax and cleaning requirements. The kit is suitable for use by the entire family, offering gentle options for children and effective solutions for adults.

The kit typically includes:

- **Tweezers:** For removing stubborn earwax.
- **Scooper:** For general earwax removal.
- **Scraper:** For extracting debris.
- **Sponge:** For removing wet earwax.
- Additional assorted tips and accessories for comprehensive ear care.



Image 2.1: Overview of the 14-piece Bebird replacement tips kit, displaying the various types of tips and their storage container.

3. SETUP AND INSTALLATION

To replace or install a cleaning tip on your Bebird Earsight Plus device, follow these steps:

1. **Remove the Current Tip:** Gently unscrew and remove the existing ear spoon head or tip from the device.
2. **Attach the New Tip:** Select the desired replacement tip. Align the new tip with the lens of your Bebird device and carefully screw it into place until it is securely fastened. Ensure it is snug but do not overtighten.



Image 3.1: Visual guide for attaching replacement tips to the Bebird device, illustrating the unscrewing and screwing process.

4. OPERATING INSTRUCTIONS

The Bebird replacement tips are designed for precision cleaning, catering to various earwax types. Always use the Bebird device as per its primary instructions and select the appropriate tip for the task.

4.1. Matching Tips to Earwax Types

Different tips are optimized for specific earwax conditions:

- **Large Earwax:** May require initial loosening or specific tools for removal.
- **Stubborn Earwax:** Best addressed with the tweezer tip for secure gripping and extraction.
- **Earwax Debris:** The scraper tip is effective for gently dislodging and removing small particles.
- **Dry Earwax / Wet Earwax:** The sponge tip is suitable for absorbing and removing wet earwax, while other tips can manage dry earwax.

4 Types of Cleaning Tips

Precision Cleaning for Every Earwax Type



Large earwax



Stubborn earwax



Earwax debris



Dry ear wet ear



Image 4.1: Visual representation of different earwax types, illustrating the need for specialized cleaning tips.

4.2. Tip-Specific Usage

Refer to the following guide for the intended use of each primary tip:

- **Tweezers:** Use for grasping and removing stubborn or larger pieces of earwax. Position the tweezers around the earwax and gently pull it out.
- **Scooper:** Ideal for general earwax removal. Carefully scoop out earwax along the ear canal walls.
- **Scraper:** Designed to extract earwax debris. Gently scrape along the ear canal to dislodge and remove small particles.
- **Sponge:** Best for wet earwax. The sponge material absorbs and helps remove moist earwax.

Match Tips to Your Earwax




Image 4.2: Detailed illustrations demonstrating the application of each cleaning tip (tweezers, scooper, scraper, sponge) within the ear canal.

4.3. Safety and Family Use

The tips are designed to be soft, safe, and sanitary for all ages, including seniors, adults, and children. Always exercise caution and gentle movements when cleaning ears. Consult a medical professional if you experience discomfort or have pre-existing ear conditions.

Soft. Safe & Sanitary Perfect for All Ages



Seniors



Adults



Kids



Image 4.3: Depiction of the tips' suitability for all age groups, emphasizing their soft, safe, and sanitary design.

5. MAINTENANCE AND CARE

To ensure hygiene and prolong the life of your Bebird replacement tips, it is recommended to clean them after each use. The tips are designed to withstand multiple sterilizations.

- **Cleaning:** After each use, remove the tip from the device. Clean the tip thoroughly with an appropriate disinfectant or sterilizing solution.
- **Drying:** Allow the tips to air dry completely before storing them.
- **Storage:** Store the cleaned and dried tips in their original container or a clean, dry place to prevent contamination.

6. TROUBLESHOOTING

If you encounter any issues with the Bebird Earsight Plus Ear Cleaner Replacement Tips, consider the following:

- **Tip Not Attaching Securely:** Ensure the tip is correctly aligned with the device's lens and screwed in straight. Do not force it.
- **Ineffective Cleaning:** Verify that you are using the correct tip for the type of earwax you are attempting to remove. Ensure the tip is clean and free of previous debris.
- **Discomfort During Use:** Immediately cease use if you experience any pain or discomfort. Ensure gentle movements and do not insert the device too deeply into the ear canal.

For persistent issues or questions, please refer to the support information in Section 8.

7. SPECIFICATIONS

Package Dimensions	3.15 x 2.36 x 0.55 inches
Item Weight	0.16 ounces
Manufacturer	Heifeng Zhizao (Shenzhen) Technology Co., Ltd
ASIN	B0DX2BJWF2
Date First Available	February 13, 2025




8. WARRANTY AND SUPPORT

Bebird stands by the quality of its products. These replacement tips come with a quality guarantee. If you encounter any issues or have questions regarding the product, please do not hesitate to contact Bebird customer support. They are available to provide hassle-free assistance.

For support, please visit the official Bebird store or contact their customer service channels as provided on the product packaging or official website.

© 2025 BEBIRD. All rights reserved.

Related Documents

	<p>Bebird Smart Visual Ear Cleaner SE & Plus Series User Manual</p> <p>Comprehensive guide for Bebird Smart Visual Ear Cleaner SE and Plus series, covering device overview, app connection, accessory installation, troubleshooting, product specifications, and safety warnings.</p>
	<p>Bebird R1+ Series Smart Visual Ear Cleaner User Manual Setup, Specs, Troubleshooting</p> <p>Comprehensive user manual for the Bebird R1+ Series Smart S Visual Ear Cleaner. This guide covers product setup, accessory lists for models R1+, R1+ S, and R1+ Plus, accessory replacement instructions, troubleshooting common issues, detailed product specifications, and important safety precautions. Learn how to connect, use, and maintain your Bebird ear cleaning device.</p>
	<p>bebird R3 Smart Visual Ear-Clean Rod User Manual</p> <p>Comprehensive user manual for the bebird R3 Smart Visual Ear-Clean Rod, detailing setup, usage, troubleshooting, and specifications. Learn how to connect, replace accessories, and safely use this innovative ear cleaning device.</p>



[bebird Home 30S Smart Ear Cleaning Device User Manual](#)

Comprehensive user manual for the bebird Home 30S Smart Ear Cleaning Device, covering setup, app connection, accessory usage, troubleshooting, and safety guidelines.