

[Manuals.plus](#) /

› [ALLJOY](#) /

› ALLJOY Hand Massager Model AJHDM2402 Instruction Manual

ALLJOY AJHDM2402

ALLJOY Hand Massager Instruction Manual

Model: AJHDM2402 | Brand: ALLJOY

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your ALLJOY Hand Massager. Please read all instructions carefully before operating the device and retain this manual for future reference.



Image 1: ALLJOY Hand Massager unit with its retail packaging.

2. IMPORTANT SAFETY INFORMATION

- Avoid use on broken skin.
- Do not overuse the device.
- Ensure the device is used on a stable surface.
- Consult a physician before use if you have pre-existing medical conditions.
- Keep out of reach of children.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This device is intended for targeted treatment, mainly of the hands, and can be used by adults over 18 years of age.

3. PACKAGE CONTENTS

- ALLJOY Hand Massager Unit
- USB Charging Cable
- Finger Massager Tool
- Storage Bag

Air Compression & Kneading



Image 2: The hand massager unit shown with its charging cable and finger massager tool.

Who Should Not Use Hand Massagers?

Not recommended for the following users



1. With Pacemaker



2. Recent Hand Injury



3. Skin Conditions



4. Pregnant Women

Image 3: A closer view of the included finger massager tool.

4. PRODUCT FEATURES

- **Upgraded Kneading Massage:** Targets acupressure points in the hands, improving circulation and promoting relaxation.
- **Multiple Functions:** Offers 2 levels of shiatsu intensity, 3 massage modes, 3 massage intensities, and 2 heating levels for a customizable experience.
- **Soothing Heat Therapy:** Features 3 heat settings (104°F, 113°F, 122°F) to enhance blood flow and loosen muscles, soothing soreness and stiffness.
- **Portable & Chargeable:** Equipped with a built-in 2500 mAh lithium battery, providing up to 170 minutes of continuous use on the lowest setting. Features a 15-minute auto-off timer.

5. SETUP

5.1 Charging the Device

Before first use, fully charge the hand massager. Connect the provided USB charging cable to the device and a suitable power source (e.g., laptop, power bank, wall adapter). A full charge typically takes approximately 2.5 hours. The device cannot be used while charging.

Full Coverage Hand Massager



Image 4: The hand massager unit being charged via a USB cable connected to a laptop.

5.2 Inserting Your Hand

Place your hand into the massager opening, ensuring your palm is aimed at the massage head for optimal effect. The open design accommodates most hand sizes.

6. OPERATING INSTRUCTIONS

The ALLJOY Hand Massager features an intuitive control panel for easy operation.

3 Levels of Heat Adjustable

Constant hot compress on the hand,
quickly relax the hand muscles



Image 5: The control panel of the hand massager, showing various settings.

6.1 Power On/Off

- To turn on: Press and hold the Power button for 2 seconds.
- To turn off: Press and hold the Power button for 2 seconds.

6.2 Adjusting Settings

- **Shiatsu:** Select between 2 levels of shiatsu intensity or turn off.
- **Heat:** Choose from 3 heat levels: Level 1 (104°F), Level 2 (113°F), Level 3 (122°F).
- **Intensity:** Adjust the massage intensity through 3 levels.
- **Mode:** Select from 3 different massage modes.
- **Timer:** The device features a 15-minute auto-off timer.

Personalized Massage at Your Fingertips

One-touch physical button start



Shiatsu

— low — high — off



Heating

— low — medium
— high — off



Intensity

— low — medium — high



Mode

3 modes to choose



Power On/Off

Long press for 2s to turn on or off

Image 6: Illustration of the heat therapy function of the hand massager.

6.3 Product Demonstration Video

Your browser does not support the video tag.

Video 1: Official product video demonstrating the features and operation of the ALLJOY Hand Massager.

Your browser does not support the video tag.

Video 2: A short preview video showcasing the CINCOM Hand Massager in various settings.

7. MAINTENANCE

7.1 Cleaning

Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water. Ensure the device is unplugged before cleaning.

7.2 Storage

Store the massager in a cool, dry place away from direct sunlight and extreme temperatures. Use the provided storage bag to protect the device when not in use.

8. TROUBLESHOOTING

- **Device not turning on:** Ensure the device is fully charged. If the issue persists, contact customer support.
- **No heat:** Verify that the heat function is activated on the control panel.
- **Weak massage intensity:** Adjust the intensity level using the control panel.

9. SPECIFICATIONS

Attribute	Value
Product Dimensions	3 x 3 x 3 inches
Item Model Number	AJHDM2402
Date First Available	September 14, 2024
Use for	Hands
Power Source	Battery Powered
Material	Polyethylene (PE)
Item Weight	0.8 Kilograms (1.76 Pounds)
Battery Capacity	2500 mAh



Image 7: Dimensions of the ALLJOY Hand Massager.

10. WARRANTY AND SUPPORT

The ALLJOY Hand Massager comes with a 365-day worry-free warranty. For any issues or inquiries, please contact our friendly and timely customer support team.