

LEIPUT S3

LEIPUT S3 2-in-1 Visual Otoscope Ear Wax Removal Tool User Manual

Model: S3 | Brand: LEIPUT

1. INTRODUCTION

Thank you for choosing the LEIPUT S3 2-in-1 Visual Otoscope Ear Wax Removal Tool. This device combines a high-definition camera with versatile ear cleaning tools, including omnidirectional tweezers and silicone ear spoons, to provide a clear and effective solution for ear wax removal and ear canal inspection. Please read this manual thoroughly before use to ensure safe and optimal operation.

2. SAFETY INFORMATION

- Do not insert the device too deeply into the ear canal to avoid injury.
- Always use the silicone ear spoon covers when operating in Ear Spoon Mode to prevent direct contact of metal with the ear canal.
- Consult a medical professional if you experience pain, discomfort, or have a history of ear infections or perforated eardrums.
- Keep out of reach of children. This device is not a toy.
- Ensure the device is fully charged before use to prevent unexpected power loss during operation.
- Do not attempt to disassemble or modify the device.
- The device is IP67 waterproof for the otoscope, ear spoon, and lens only. Do not immerse the entire device in water.

3. PACKAGE CONTENTS

Verify that all items are present in your package:

- LEIPUT S3 Ear Camera with Omnidirectional Tweezers (1)
- Silicone Earplug Cover (6)
- Earplug Cover Storage Bottle (1)
- Charging Cable (1)
- Instruction Manual (1)



Image 1: The LEIPUT S3 2-in-1 Visual Oscope Ear Wax Removal Tool shown with its included accessories: the main device with tweezers, silicone ear spoon covers, a storage bottle for the covers, and a USB charging cable. The device is blue.

4. PRODUCT OVERVIEW

The LEIPUT S3 integrates a high-resolution camera with a versatile ear cleaning mechanism. Key components include:

- **Visual Oscope Camera:** Provides a clear view of the ear canal.
- **Omnidirectional Tweezers:** Designed for precise removal of larger earwax pieces.
- **Silicone Ear Spoons:** Attach to the tweezers for gentle cleaning of remaining earwax.
- **Integrated LED Lights:** Illuminate the ear canal for optimal visibility.
- **6-Axis Gyroscope:** Ensures stable and smooth image transmission.
- **Wi-Fi Module:** Connects to your smartphone for live viewing.

We have upgraded ear wax removal, 6-Axis Gyroscope and 20fps transmission speed make the transmission more stable and smooth, visual effect is clearer.



4Led
Lights



3.5mm
mini lens



2K
Camera



6-Axis
gyroscope



Image 2: An exploded diagram illustrating the internal components of the LEIPUT S3 device, highlighting the 4 LED lights, 3.5mm mini lens, 2K camera, and 6-axis gyroscope, which contribute to its stable and clear visual output.

5. SETUP GUIDE

5.1 Charging the Device

1. Connect the provided charging cable to the device's charging port.
2. Plug the other end of the cable into a standard USB power adapter (not included) or a computer's USB port.
3. The indicator light on the device will show charging status (refer to device for specific light behavior). Charge fully before first use.

5.2 App Installation and Wi-Fi Connection

1. **Download the App:** Search for the "LEIPUT Ear Cleaner" app (or similar name as indicated in the packaging) on the Apple App Store (for iOS) or Google Play Store (for Android). Install the app on your smartphone.
2. **Power On the Device:** Press and hold the power button on the LEIPUT S3 device until the indicator light turns on.
3. **Connect via Wi-Fi:**
 - Open your smartphone's Wi-Fi settings.
 - Locate and select the Wi-Fi network named "LEIPUT-XXXX" (where XXXX is a unique identifier).
 - No password is typically required for this direct device connection.
4. **Launch the App:** Open the installed LEIPUT app. The live view from the otoscope camera should now be displayed on your screen.

6. OPERATING INSTRUCTIONS

6.1 General Usage Tips

- Always ensure a stable Wi-Fi connection for a clear and uninterrupted view.

- Move the device slowly and gently within the ear canal.
- Observe the image on your phone screen carefully to guide your movements.
- Avoid touching the eardrum.

6.2 Tweezers Mode

This mode is suitable for removing larger pieces of earwax or foreign objects visible in the ear canal.

1. Ensure no silicone ear spoon cover is attached to the tweezers.
2. Carefully insert the otoscope into the ear canal while observing the live feed on your phone.
3. Position the tweezers around the earwax or object you wish to remove.
4. Gently close the tweezers to grasp the item.
5. Slowly and carefully withdraw the device from the ear canal.

6.3 Ear Spoon Mode

This mode is for cleaning residual earwax and is generally safer for general ear cleaning.

1. Select a silicone ear spoon cover from the storage bottle.
2. Gently press the 360-degree switch on the tweezers to open them slightly.
3. Carefully place the silicone ear spoon cover onto the tweezers. The tweezers will remain open to secure the cover and increase friction, preventing it from dropping.
4. Insert the otoscope with the attached ear spoon into the ear canal.
5. Using the live view, gently scoop or scrape earwax from the ear canal walls.
6. Withdraw the device and clean the ear spoon as needed.

Immersion cleaning cerumen

Safer to be visible

Put on the silicone
cover and you can
clean your ears

Open the tweezers
and remove the
large pieces of earwax

Tweezers Mode



Grab huge
chunks of earwax



Grabbing foreign
objects in the ear canal



Suitable for people prone
to ear canal infections

Ear Spoon Mode



Easy cleaning
of earwax

Image 3: This image demonstrates the two primary operating modes of the LEIPUT S3. On the left, "Tweezers Mode" shows the device's metal tweezers designed for grabbing larger earwax chunks or foreign objects. On the right, "Ear Spoon Mode" illustrates the silicone cover attached to the tweezers, used for gentle scraping and cleaning of earwax.

7. MAINTENANCE AND CLEANING

Proper cleaning ensures hygiene and extends the lifespan of your device.

- After each use, remove the silicone ear spoon cover and dispose of it if it's single-use, or clean it thoroughly if reusable.
- The otoscope tip, ear spoon, and lens are IP67 waterproof. They can be rinsed under running water or wiped with an alcohol swab.
- **Important:** Do not immerse the entire device body in water. Only the tip is waterproof.
- Use a soft, dry cloth to wipe the main body of the device.
- Store the device and accessories in a clean, dry place.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not power on.	Low battery.	Charge the device fully using the provided USB cable.
Cannot connect to Wi-Fi.	Device not powered on; incorrect Wi-Fi network selected; Wi-Fi interference.	Ensure device is on. Select the correct "LEIPUT-XXXX" network. Try restarting the device and your phone. Move closer to the device.
No image or blurry image in app.	Poor Wi-Fi connection; dirty lens; app malfunction.	Check Wi-Fi connection. Clean the camera lens carefully with an alcohol swab. Restart the app.
Tweezers/Ear Spoon not functioning correctly.	Obstruction; mechanical issue.	Ensure no debris is blocking the mechanism. If the issue persists, contact customer support.

9. SPECIFICATIONS

- **Model:** S3
- **Dimensions:** 7.24 x 3.11 x 0.67 inches
- **Weight:** 3.28 ounces
- **Battery:** 1 Lithium Polymer battery (included)
- **Camera Resolution:** 2K Quad HD (as per image, implied)
- **Lens Diameter:** 3.5mm (as per image)
- **Gyroscope:** 6-Axis
- **Waterproof Rating:** IP67 (for otoscope tip, ear spoon, and lens)
- **Connectivity:** Wi-Fi (iOS & Android compatible)

10. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the contact details provided on the product packaging or visit the official LEIPUT website. Please have your model number (S3) and purchase date available when contacting support.

11. LEGAL DISCLAIMER

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0F3J59DYG>

More documents and user manuals: **<https://manuals.plus/>**