

Aopvui FZ-608

Aopvui IPL Hair Removal Device FZ-608 User Manual

Brand: Aopvui | Model: FZ-608

1. INTRODUCTION

The Aopvui IPL Hair Removal Device FZ-608 utilizes Intense Pulsed Light (IPL) technology to effectively reduce unwanted hair growth for long-lasting results. Designed for both women and men, this device features an advanced ice-cooling system for a nearly painless experience. With 9 energy levels and 2 flash modes, it offers customizable treatment for various body areas, including armpits, legs, arms, and the bikini line. This manual provides essential information for safe and effective use of your device.



Figure 1: Aopvui IPL Hair Removal Device FZ-608 and accessories.

2. SAFETY INFORMATION

Before using the Aopvui IPL Hair Removal Device FZ-608, please read all instructions carefully. This device is intended for hair removal and should not be used for any other purpose. Always wear the provided protective goggles during operation to shield your eyes from the intense pulsed light. Do not use the device on irritated, sunburned, or broken skin. Consult a physician if you have any pre-existing skin conditions or concerns.

Skin and Hair Type Compatibility: IPL technology is most effective on lighter skin tones with darker hair. It is generally unsuitable for very dark skin tones (black, dark brown) and very light hair colors (red, light blond, white/grey) as it may not be effective or could cause adverse reactions. Refer to the compatibility chart in Section 4 for detailed guidance.

Important Note: The device should not be used more than three times on the same spot during a single treatment session to avoid skin irritation.

3. PACKAGE CONTENTS

Your Aopvui IPL Hair Removal Device FZ-608 package includes the following items:

- 1 x Aopvui IPL Hair Removal Device FZ-608
- 1 x Shaver
- 1 x Protective Goggles
- 1 x User Manual
- 1 x Power Adapter

4. SKIN AND HAIR TYPE COMPATIBILITY

For optimal results and safety, it is crucial to determine if your skin tone and hair color are suitable for IPL hair removal. The device targets melanin in the hair follicle, so a contrast between hair color and skin tone is necessary.



Figure 2: Skin Tone and Hair Color Compatibility Chart. Green checkmarks indicate suitable combinations, while red crosses indicate unsuitable combinations.

Unsuitable Combinations:

- **Dark Skin Colors:** Black, Dark Brown. The high melanin content in very dark skin can absorb too much light energy, leading to burns or discoloration.
- **Light Hair Colors:** Red, Light Blond, White, Grey. These hair colors contain insufficient melanin for the IPL device to effectively target and disable the hair follicle.

5. SETUP

1. **Prepare the Skin:** Before each treatment, ensure the skin area is clean, dry, and free of any lotions, oils, or deodorants. Shave the area to be treated to remove surface hair. This allows the IPL light to effectively reach the hair follicles beneath the skin.
2. **Wear Protective Goggles:** Always put on the provided protective goggles before plugging in or operating the device to protect your eyes from the intense light flashes.
3. **Connect Power:** Plug the power adapter into the device and then into a power outlet. The device display will illuminate.

6. OPERATING INSTRUCTIONS

6.1 Power On/Off and Mode Selection

1. **Power On:** Press and hold the power button (usually marked with a circle and vertical line) for a few seconds until the display turns on.
2. **Mode Selection:** Short press the power button to switch between HR (Hair Removal), SC (Skin Care), and RA (Rejuvenation) modes. The default mode is typically Hair Removal (HR).
3. **Energy Level Adjustment:** Use the '+' and '-' buttons to adjust the energy level from 1 (lowest) to 9 (highest). It is recommended to start with Level 1 for your first treatment and gradually increase as your skin tolerates it.
4. **Flash Modes:** The device offers two flash modes:
 - **Manual Mode:** Short press the flash button to emit a single flash. Ideal for small or sensitive areas like the upper lip or bikini line.
 - **Auto Mode:** Press and hold the flash button for 3 seconds to activate Auto Mode. The device will automatically flash when it detects proper skin contact. This mode is suitable for larger areas like legs, arms, and back.



Figure 3: Auto Mode for large areas (arms, back, legs) and Manual Mode for small areas (lip, armpit, bikini line).

6.2 Application Technique

1. **Skin Contact:** Place the device window firmly against your skin, ensuring full contact. The device will only flash when proper skin contact is detected.
2. **Treatment:** For Manual Mode, press the flash button. For Auto Mode, glide the device slowly over the skin. Do not treat the same area more than three times per session.
3. **Ice-Cooling System:** The integrated ice-cooling system maintains a comfortable skin temperature (around 8°C) during treatment, minimizing discomfort and making the process nearly painless.

AFTER 3 WEEKS WITH VISIBLE RESULT



Figure 4: The ice-cooling function ensures a comfortable treatment experience.

6.3 Official Product Videos

Your browser does not support the video tag.

Video 1: Product Introduction Video for Aopvui IPL Hair Removal Device FZ-608. This video provides a general overview of the device and its features.

Your browser does not support the video tag.

Video 2: Introduction Video for FZ608G. This video demonstrates the basic functions and usage of the Aopvui IPL Hair Removal Device.

7. TREATMENT SCHEDULE AND EXPECTED RESULTS

Consistency is key for effective IPL hair removal. Follow the recommended schedule for optimal results:

- **Weeks 1-2:** Use 3-4 times a week. Hair growth may slow down.
- **Weeks 3-5:** Use 2-3 times a week. You should observe visibly reduced hair.

- **Weeks 6-8:** Use 1-2 times a week. Significant hair reduction and thinner hair growth are expected.
- **After 8 Weeks:** Use once every 2-4 weeks for maintenance, or as needed to maintain hair-free skin.

The Benefits of At-Home Hair Removal VS Other Hair Removal Methods



Hair removal creams



Waxing



Razors



Salon Laser Experience

Figure 5: Illustrative results of hair reduction after consistent IPL treatment.

8. MAINTENANCE

To ensure the longevity and optimal performance of your Aopvui IPL Hair Removal Device FZ-608:

- **Cleaning:** After each use, unplug the device and gently wipe the treatment window with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and moisture. Keep it out of reach of children.

9. TROUBLESHOOTING

If you encounter any issues with your device, please refer to the following common troubleshooting tips:

- **Device Not Turning On:** Ensure the power adapter is securely plugged into both the device and a working power outlet.
- **No Flash:** Check that the device window is in full contact with your skin. The device has a built-in skin contact sensor and will not flash if contact is insufficient. Also, ensure the energy level is set and the flash button is pressed correctly (or held for Auto Mode).
- **Skin Discomfort/Redness:** If you experience discomfort, reduce the energy level. If redness persists, discontinue use and consult a healthcare professional.
- **Ineffective Results:** Verify that your skin tone and hair color are compatible with IPL technology (refer to Section 4). Ensure you are following the recommended treatment schedule consistently.

10. SPECIFICATIONS

- **Model Number:** FZ-608
- **UPC:** 783821620021
- **ASIN:** B0F4MFZ5ZL
- **Package Dimensions:** 11.69 x 7.2 x 2.64 inches
- **Item Weight:** 2.23 Pounds
- **Energy Levels:** 9
- **Flash Modes:** 2 (Manual, Auto)
- **Flashes:** 999,999
- **Manufacturer:** Shenzhen IONKA Medical Technology Co., Ltd.

11. WARRANTY AND SUPPORT

For warranty information or technical support regarding your Aopvui IPL Hair Removal Device FZ-608, please refer to the contact details provided in the user manual included in your package or visit the official Aopvui website. Keep your purchase receipt as proof of purchase for any warranty claims.

© 2025 Aopvui. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question