

BTFO G28014

BTFO USB Charging Cable for Momcozy Breast Pumps - User Manual

1. INTRODUCTION

Thank you for choosing the BTFO USB Charging Cable. This manual provides essential information for the safe and effective use of your new charging cable. Please read these instructions carefully before use and retain them for future reference.

This product is designed to provide reliable power delivery for compatible Momcozy hands-free breast pumps.

2. PRODUCT OVERVIEW

The BTFO USB Charging Cable is a durable and convenient solution for charging your Momcozy breast pumps. It features a Y-split design, allowing two breast pumps to be charged simultaneously from a single USB power source.

Key Features:

- **Dual Charging Capability:** Charge two compatible breast pumps at once.
- **Optimal Length:** 22.44 inches (57 cm) for flexible use.
- **Durable Construction:** Made from strong PVC material for longevity.
- **Wide Compatibility:** Designed for Momcozy M5, M6, M9, S9 Pro, and S12 Pro models.

Compatibility:

This charging cable is compatible with the following Momcozy hands-free breast pump models:

- Momcozy M5
- Momcozy M6
- Momcozy M9
- Momcozy S9 Pro
- Momcozy S12 Pro

COMPATIBILITY



Figure 2.1: The BTFO USB charging cable with visual indicators of compatible Momcozy breast pump models (M5, S12 Pro, M9, S9 Pro, M6).



Figure 2.2: Detailed view of the BTFO USB charging cable, highlighting the USB-A input and dual micro-USB outputs.

3. SETUP

Follow these steps to set up your BTFO USB Charging Cable:

1. **Connect to Power Source:** Insert the USB-A end of the charging cable into a compatible USB power adapter (not included), a computer USB port, or a power bank. Ensure the power source provides adequate power for charging.
2. **Connect to Breast Pump(s):** Insert the micro-USB ends of the cable into the charging ports of your Momcozy breast pump(s). You can connect one or two pumps simultaneously.
3. **Verify Connection:** Ensure both ends of the cable are securely connected. The breast pump's charging indicator should activate, typically showing a light or screen icon indicating it is receiving power.

CHARGER FOR MOMCOZY M5 BREAST PUMP



Figure 3.1: Example of two Momcozy breast pumps connected to a wall adapter for charging using the BTFO cable.

4. OPERATING INSTRUCTIONS

Once connected, the BTFO USB Charging Cable will begin to charge your Momcozy breast pump(s). The charging process is automatic. Monitor the breast pump's battery indicator to determine charging status and when the battery is full.

Charging Tips:

- For optimal charging speed, use a USB power adapter with sufficient output (e.g., 5V/2A).
- The cable allows for convenient charging in various locations, including at home, in the office, or during travel.

CONVENIENT OUTDOOR CHARGING



Figure 4.1: The BTFO charging cable facilitating convenient charging of breast pumps in a vehicle setting.

5. MAINTENANCE

Proper care and maintenance will extend the lifespan of your BTFO USB Charging Cable.

Cleaning:

- Wipe the cable with a soft, dry cloth to remove dust or dirt.
- Do not use harsh chemicals or abrasive cleaners.

Storage:

- Store the cable in a cool, dry place away from direct sunlight and extreme temperatures.
- Avoid tightly coiling or bending the cable at sharp angles to prevent internal wire damage.

Material Construction:

The cable features a durable PVC insulation sheath and premium copper wire for efficient and safe power transmission.



SAFE AND RELIABLE

Premium copper wire




Figure 5.1: Internal structure of the cable, highlighting the copper wire and PVC insulation for safety and reliability.

6. TROUBLESHOOTING

If you encounter issues with your charging cable, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Breast pump not charging	Loose connection Faulty power source Incompatible breast pump model	Ensure all connections are secure. Try a different USB power adapter or port. Verify your breast pump model is listed as compatible (Momcozy M5, M6, M9, S9 Pro, S12 Pro).
Slow charging	Low power output from source Using a non-optimal power adapter	Use a USB power adapter with higher output (e.g., 5V/2A). Ensure the power source is not overloaded.
Cable appears damaged	Physical stress or wear	Discontinue use immediately. Replace the cable to prevent potential electrical hazards.

7. SPECIFICATIONS

Feature	Detail
Brand	BTFO
Model Name	G28014
Cable Type	USB (USB-A to Dual Micro-USB)
Connector Type	Auxiliary
Cable Length	22.44 inches (57 cm)
Color	White
Material	PVC
Compatible Devices	Momcozy M5, M6, M9, S9 Pro, S12 Pro Hands-Free Breast Pumps
Special Feature	Lightweight
Unit Count	2.0 Count (2 cables per package)
Item Weight	1.48 ounces (42 grams)

8. WARRANTY INFORMATION

Specific warranty details for this product are not provided in the available information. Please refer to the retailer's return policy or contact the seller directly for warranty inquiries.

9. SUPPORT

For further assistance or questions regarding your BTFO USB Charging Cable, please contact the seller through the platform where the purchase was made. Ensure you have your order details available for faster service.