

[Manuals.plus](#) /

› [Jaioyios](#) /

› Jaioyios AF1332-002 Movable Electric Floor Lift Chair User Manual

## Jaioyios AF1332-002

# Jaioyios AF1332-002 Movable Electric Floor Lift Chair User Manual

Model: AF1332-002

## 1. INTRODUCTION

This manual provides essential instructions for the safe and effective use of your Jaioyios Movable Electric Floor Lift Chair. This device is designed to assist individuals in rising from the floor and transferring to a bed or sofa. Please read this manual thoroughly before initial use and retain it for future reference.

### Key Features:

- **Movable Design:** Equipped with four wheels for easy repositioning.
- **Electric Lift Mechanism:** Smooth and stable lifting from the floor.
- **Foldable Structure:** Enhanced design for swift setup and compact storage.
- **Durable Construction:** Made from high-grade metal with a PU seat, supporting up to 400 lbs.
- **High-Capacity Battery:** Provides over 60 lift cycles per charge.
- **Swivel Wheels with Lock:** Two front swivel wheels for 360° rotation and locking mechanism for stability.

## 2. IMPORTANT SAFETY GUIDELINES


Adherence to these safety guidelines is crucial for preventing injury and ensuring proper operation of the lift chair.

- **Indoor Use Only:** This lift chair is designed exclusively for indoor environments. Do not use it outdoors.
- **Not a Wheelchair:** This device is not intended to be used as a wheelchair for transportation.
- **Weight Capacity:** Do not exceed the maximum weight capacity of 400 lbs (181 kg).
- **User Capability:** Suitable for conscious users who can actively participate in the transfer process and use hand strength independently.
- **Inappropriate Use:** Do not use for individuals with paralysis, dementia, or those unable to self-support or follow instructions.
- **Surface Compatibility:** Avoid using on deep-pile or long-haired carpets, as this may impede wheel movement and stability.
- **Wheel Locks:** Always engage the wheel locks before a user sits on or gets off the chair, and before initiating a lift.
- **Supervision:** Always ensure a trained caregiver is present to assist and supervise during operation.

- **Clear Area:** Ensure the area around the lift chair is clear of obstructions before operation.
- **Power Cord:** Keep the power cord away from heat sources, sharp edges, and wet areas.
- **Battery Care:** Use only the provided charger. Do not attempt to open or modify the battery pack.



Image: The lift chair is intended for indoor use only and should not be used as a wheelchair.



## A Thoughtful Gift for Parents and Loved Ones

Image: This diagram illustrates that the lift chair is suitable for conscious users who can assist themselves, and not for individuals with paralysis, dementia, or inability to self-support.

### 3. PACKAGE CONTENTS

Upon unpacking, verify that all components are present and undamaged.

- Jaioyios Movable Electric Floor Lift Chair (main unit)
- Armrests (2)
- Handheld Remote Control
- Rechargeable Lithium-Ion Battery
- Power Adapter/Charger
- User Manual (this document)

### 4. ASSEMBLY

The Jaioyios lift chair features an enhanced foldable design for quick and easy setup.

1. **Unpack:** Carefully remove all components from the packaging.
2. **Unfold the Chair:** Place the main unit on a flat, stable surface. Gently unfold the chair until it locks into its upright position.
3. **Attach Armrests:** Slide the armrests into the designated slots on each side of the chair until they click securely into place.
4. **Connect Remote and Battery:** Connect the handheld remote control and the battery to their respective ports on the chair. Ensure connections are firm.



Image: The Jaioyos lift chair in its unfolded, ready-to-use configuration, showing the attached armrests and remote control.

# High-Capacity Battery


Charging time: 2-3 hours



60 lift operations  
per full charge



Screen showing  
battery level



Rechargeable

Image: This illustration shows the lift chair's dimensions when unfolded and its compact form when folded for storage.

## 5. OPERATION

### 5.1. Charging the Battery

The lift chair is powered by a rechargeable Lithium-Ion battery.

1. Locate the battery pack and the charging port.
2. Connect the power adapter to the battery pack, then plug the adapter into a standard electrical outlet.
3. Charging time is approximately 2-3 hours for a full charge. The battery display will indicate the charging status.
4. A fully charged battery provides over 60 lift cycles.

# LIFT CHAIR DIMENSIONS

Assembly-free folding design: unfold for easy use



PRODUCT WEIGHT

**44.2 LBS**

**Packing Illustration  
(Assembly-Free)**



Image: The 2200mAh high-capacity battery, indicating a charging time of 2-3 hours and providing over 60 lift operations per full charge.

## 5.2. Moving the Lift Chair

The lift chair is equipped with four wheels for easy movement. The two front wheels are swivel wheels with a locking mechanism.

1. **Unlock Wheels:** To move the chair, ensure the wheel brakes on the front swivel wheels are disengaged. Press the lever up to release the brake.
2. **Positioning:** Gently push or pull the chair to the desired location. The swivel wheels allow for 360° rotation.
3. **Lock Wheels:** Once the chair is in position, press the lever down on the front swivel wheels to engage the brakes and secure the chair.



## Other Lift Chair

**Stuck in Place –  
Heavy & Hard to Move**



## Our Lift Chair

**Effortless Mobility –  
no straining, no hassle!**



Image: This image highlights the effortless mobility of the lift chair due to its wheels, contrasting it with a stationary chair.

# Electric Lift Chair

Equipped with 2 front swivel wheels for 360° rotation and locking.



Press the brake to lock the wheels



Release the brake before moving


Image: A detailed view of the lift chair's wheel, showing how to press the brake to lock the wheels and release it before moving.

## 5.3. Assisting a User from the Floor

Always ensure the user is conscious and able to follow instructions and assist with their own transfer.

1. **Prepare the Chair:** Move the lift chair close to the user on the floor. Ensure the chair is in its lowest, reclined position. Lock the front swivel wheels.
2. **Position the User:** Guide the user to position themselves onto the seat of the lift chair. Ensure they are seated securely and comfortably.
3. **Initiate Lift:** Using the handheld remote control, press the "UP" button to slowly raise the chair. Maintain visual contact with the user and communicate throughout the process.
4. **Monitor Stability:** Observe the user's stability and comfort during the lift. If any discomfort or instability occurs, pause the lift.

5. **Transfer:** Once the user is at a suitable height, unlock the wheels and carefully move the chair to the target location (e.g., bed, sofa). Lock the wheels again before the user attempts to transfer off the chair.
6. **Lowering:** To lower the chair, press the "DOWN" button on the remote control.



## How to Use the Lift Chair




Image: This visual guide demonstrates the process of unfolding the chair, locking wheels, assisting a person onto the chair, unlocking wheels for movement, and transferring the person to a bed or sofa.

## 6. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your lift chair.

- **Cleaning:** Wipe down the chair's surfaces, including the PU seat, with a damp cloth and mild detergent. Avoid abrasive cleaners.
- **Inspection:** Periodically inspect all moving parts, connections, and the power cord for any signs of wear, damage, or loose components.
- **Battery Storage:** If storing the chair for an extended period, ensure the battery is partially charged (around 50%) and store it in a cool, dry place. Recharge every few months to maintain battery health.
- **Professional Service:** Do not attempt to repair the motor or electrical components yourself. Contact customer support for professional service.

## 7. TROUBLESHOOTING

If you encounter issues with your lift chair, refer to the following common problems and solutions.

| Problem                      | Possible Cause                                                    | Solution                                                                                         |
|------------------------------|-------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| Chair does not lift or move. | Battery is low or not connected. Remote control is not connected. | Ensure the battery is fully charged and securely connected. Check the remote control connection. |

| Problem                                | Possible Cause                                   | Solution                                                                                                             |
|----------------------------------------|--------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| Chair movement is jerky or unstable.   | Uneven surface. Obstruction in wheels. Overload. | Ensure the chair is on a flat, stable surface. Check wheels for obstructions. Verify user weight is within capacity. |
| Wheels do not lock or unlock properly. | Brake mechanism is dirty or damaged.             | Clean the brake mechanism. If damaged, contact customer support.                                                     |

If the problem persists after attempting these solutions, please contact Jaioyios customer support.

## 8. SPECIFICATIONS

| Feature                       | Detail                       |
|-------------------------------|------------------------------|
| Model Number                  | AF1332-002                   |
| Brand                         | Jaioyios                     |
| Maximum Weight Capacity       | 400 lbs (181 kg)             |
| Product Dimensions (Unfolded) | Approx. 27"D x 19.7"W x 28"H |
| Product Dimensions (Folded)   | 19"D x 17.7"W x 4"H          |
| Product Weight                | 44.2 lbs (20 kg)             |
| Battery Type                  | Lithium Ion (2200mAh)        |
| Charging Time                 | 2-3 hours                    |
| Lift Cycles per Charge        | Over 60 cycles               |
| Color                         | Grey                         |

## 9. WARRANTY AND CUSTOMER SUPPORT

For warranty information, technical support, or any questions regarding your Jaioyios Movable Electric Floor Lift Chair, please contact the manufacturer or your point of purchase. Keep your purchase receipt as proof of purchase.  
Manufacturer: Jaioyios

