

FlyingJoy All Terrain Walkers

FlyingJoy 300lb Foldable 2-in-1 Rollator Walker and Transport Wheelchair User Manual

Model: All Terrain Walkers

1. SAFETY INFORMATION

For your safety and the safety of others, please read and follow all instructions carefully before using this product.

- Always lock the brakes after sitting and ensure the walker is stable before attempting to walk or transfer.
- When adjusting height, ensure all height holes are aligned and the adjustment knob is securely locked.
- When using as a transport chair, utilize the safety fastening straps on the waist and footrest for improved stability and user security.
- Regularly inspect all components for wear and tear. Do not use if any part is damaged or missing.
- Ensure the product is fully unfolded and locked into position before use.

2. SETUP AND ASSEMBLY

This rollator walker arrives 95% pre-assembled for quick and easy setup. No tools or screws are required for initial assembly.

2.1 Unfolding the Rollator

Carefully remove the rollator from its packaging. Stand the unit upright and gently pull the handles apart until the frame locks into its open position.

2.2 Attaching Wheels and Handles

The wheels and handles attach using a push-button mechanism. Insert each wheel and handle into its designated slot until the button clicks into place, indicating a secure connection.



Image: All components of the rollator walker are laid out, demonstrating the tool-free assembly process.

2.3 Converting to Transport Chair Mode

To convert from rollator to transport chair mode, flip the backrest to the opposite side. Attach the detachable footrests to the front frame. Ensure they are securely locked in place.

ROLLATOR WALKER & TRANSPORT CHAIR COMBO

2 IN 1



ROLLATOR WALKER



TRANSPORT CHAIR




Image: The 2-in-1 functionality is shown, with the product configured as a rollator walker on the left and a transport chair on the right.

2.4 Assembly Video Tutorial

Video: This video demonstrates the assembly process and conversion between walker and wheelchair modes, highlighting ease of setup.

3. OPERATING INSTRUCTIONS

3.1 Adjusting Height

The handles are adjustable from 33.5 inches to 40 inches to accommodate users from 5'1" to 6'4". The seat height can be adjusted from 20.8 inches to 22.8 inches for proper leg support. To adjust, locate the adjustment pins on the handles and

seat frame, press them in, and slide the component to the desired height. Ensure the pins fully engage in the adjustment holes before use.

Adjustable & Detachable



Image: A diagram illustrating the adjustable handlebar height (33.5" to 40") and seat height (20.8" to 22.8").

3.2 Braking System

The rollator features a 3-mode braking system: Free-roll, Resistance, and Lock. Squeeze the handbrakes to engage the brakes for slowing down or stopping. Push the brake levers down to lock the wheels for secure sitting. Pull the levers up to release the lock.

BRAKING SAFETY

Easy to operate, brake immediately by gripping the handle, triangular design bracket, safer



**Brake-Provides Resistance
Safer for Going Downhill**



**Lock-Use when Sitting or to
Prevent From Rolling**



Image: A person is being pushed in the transport chair, with insets showing how to operate the brake levers for resistance and locking.

3.3 All-Terrain Wheels

Equipped with 10-inch front and 8-inch rear airless tires, this rollator is designed to navigate various surfaces smoothly, including grass, gravel paths, cracked sidewalks, indoor marble, and wood floors. The non-marking rubber protects indoor surfaces.




Image: Close-up of the 10-inch front wheel and 8-inch rear wheel, with illustrations of various terrains like grass, uneven ground, mud, marble, and wood.

3.4 Storage and Accessories


The rollator includes a large removable storage basket under the seat with dividers for organizing personal items. A hidden Velcro-secured pouch is also available for valuables. An anti-slip cup holder and a cane/umbrella holder are integrated into the design.

FREE YOUR HANDS NOT FOR STANDARD DOOR

Mesure before purchase



Cup Holder



Cane/
UmbrellaHolder



Image: Close-up views of the cup holder and the cane/umbrella holder attached to the rollator frame.

4. MAINTENANCE

The FlyingJoy rollator is designed for easy maintenance. Regular cleaning with a damp cloth is recommended. The modular design allows for individual replacement of wheels, brakes, cables, and seat cushion using basic tools. For detailed replacement instructions, refer to the video tutorial or contact customer support.

Video: This video provides an overview of the FlyingJoy rollator, including features relevant to its construction and potential maintenance.

5. TROUBLESHOOTING

If you encounter any issues with your FlyingJoy rollator, please review the following common solutions:

- **Brakes not engaging:** Check if the brake cables are properly tensioned and connected. Ensure no debris is obstructing the brake mechanism.
- **Wheels not rolling smoothly:** Inspect wheels for debris or damage. Ensure all wheel attachment pins are fully engaged.
- **Unstable frame:** Verify that the rollator is fully unfolded and all locking mechanisms are engaged. Check all adjustable components (handles, seat) are securely locked at the desired height.

For further assistance, please contact FlyingJoy customer service.

6. SPECIFICATIONS

Feature	Detail
Max Weight Capacity	300 pounds (136 kg)
Product Net Weight	24.9 pounds (11.3 kg)
Folding Size	33.5"L x 14.6"W x 35.3"H (850x370x896 mm)
Overall Width	25.2 inches (639 mm)
Overall Length	28.5 - 32.4 inches (720-780 mm)
Handle Height Adjustment	33.5 - 40 inches (850-1015 mm)
Seat Cushion Size	13.9" x 13.9" (353 x 353 mm)
Seat Height Adjustment	20.8 - 22.8 inches
Wheels	Front: 10 inches, Rear: 8 inches
Material of Frame	Aircraft-grade Aluminum
Recommended Body Height	5'1" to 6'4" (1494-1707 mm)
Color	Red

7. WARRANTY AND SUPPORT

FlyingJoy is committed to providing quality products and customer satisfaction. While specific warranty details are not provided in this manual, our friendly customer service team is available 7 days a week to assist you with any questions or concerns regarding your rollator walker.

Please refer to your purchase documentation or the FlyingJoy website for the most up-to-date warranty information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question