

[Manuals.plus](#) /

> [Kolrry](#) /

> Kolrry Ear Wax Removal Tool Camera with 1080P HD and 6 Spoons Instruction Manual

Kolrry B0F9FK5NV3

Kolrry Ear Wax Removal Tool Camera with 1080P HD and 6 Spoons Instruction Manual

Model: B0F9FK5NV3

INTRODUCTION

The Kolrry Ear Wax Removal Tool Camera is designed to provide a clear, real-time view of the ear canal for precise inspection and safe earwax removal. This device combines a high-definition camera with specialized tools to assist in personal ear hygiene.



Image: The Kolrry Ear Wax Removal Tool Camera and its components, highlighting its wireless connectivity and 1080P HD camera capabilities.

SAFETY INFORMATION

- Do not use the product while charging.
- Always supervise children when using this product.
- Be cautious when inserting and removing ear spoons to avoid injury.
- Keep the product out of reach of children when not in use.
- Before each use, ensure the silicone sleeve on the ear spoon is tightened to prevent scratching the ear canal.
- Get accustomed to how the otoscope feels in the ear canal during initial use, and move it gently and slowly.

KIT CONTENTS

The Kolrry Ear Wax Removal Tool Camera kit includes the following items:

- Otoscope with light (1)
- Type-C charging cable (1)
- Silicone ear spoons (6) - featuring 3 kinds of curved soft spoons (Low arc, Middle arc, High arc)
- Manual (1)
- 7-piece ear wax removal tool kit (1)



Image: A comprehensive view of the Kolrry Ear Wax Removal Tool Camera kit, showcasing the main device, interchangeable silicone tips, and additional metal tools.

PRODUCT FEATURES

- **1080P HD Camera:** Provides a clear, real-time picture of the ear canal.
- **Bright LED Lights:** Illuminates the dark environment of the ear canal for enhanced visibility.
- **Upgraded Rigid Integrated Design:** The camera and ear spoon are firmly fixed with a non-flexible connection to prevent wobbling or detachment, ensuring maximum safety.
- **Interchangeable Silicone Ear Spoons:** Five highly flexible silicone ear spoons with 3 distinct curved shapes (Low arc, Middle arc, High arc) are included for comfortable and effective cleaning.
- **Wide Compatibility:** Supports both Android and iOS devices.
- **Type-C Charging:** Convenient and fast charging.
- **360° Wide Angle Lens:** Offers a broad view of the ear canal.
- **IP67 Waterproof:** Ensures durability and easy cleaning of the tip.




Image: Visual representation of the device's key features, such as the 1080P camera, integrated design, and compatibility with iOS and Android phones.

SETUP

1. App Download

1. Scan the QR code provided in the manual or search for the "Suear" app in your device's app store (Google Play Store for Android, App Store for iOS).
2. Download and install the "Suear" app.

2. Device Connection

1. Press and hold the ear wax cleaner On/Off button for 1 to 2 seconds to turn it on. A blue light will illuminate.
2. Go to your device's Wi-Fi settings and connect to the Wi-Fi network named "Suear-xxxx" (where xxxx is a unique identifier).
3. Return to the "Suear" app and click "Connect the device" or "Begin to Clean Ears" to start using the device.

Note: It is normal for the prompt "Network unavailable" to appear when connecting, as the device's Wi-Fi does not support internet networking.

3. Attaching Silicone Ear Spoons

1. Select the desired silicone ear spoon based on the type of earwax or ear canal shape.
2. Carefully twist and tighten the silicone sleeve onto the tip of the otoscope. Ensure it is securely fastened to prevent detachment during use.

Wide Compatibility

Work on IOS and Android system



Android



iPhone



iPad



Image: The Kolrry device demonstrating compatibility with various mobile devices, including Android phones, iPhones, and iPads.

OPERATING INSTRUCTIONS

1. Powering On/Off

- To power on, press and hold the power button located on the device for 1-2 seconds until the blue indicator light illuminates.
- To power off, press and hold the power button again for 1-2 seconds until the light turns off.

2. Using the Camera and LED Lights



- Once connected to the app, the 1080P HD camera will provide a real-time view of your ear canal on your connected device.
- The bright LED lights will automatically illuminate the ear canal, ensuring a clear and detailed image.
- Utilize the wide-angle lens to navigate and observe different areas of the ear canal.

3. Gentle Wax Removal

- Carefully insert the otoscope tip into the ear canal while observing the real-time view on your device.
- Using the selected silicone ear spoon, gently scoop or dislodge earwax. Avoid applying excessive force or pushing wax deeper into the ear.
- Move the device slowly and deliberately to prevent accidental injury to the ear canal or eardrum.
- The different curved spoons are designed to adapt to various earwax types and ear canal anatomies for effective removal.

1080P HD Clear No delay Video Image

Macro camera 1080 HD view of the microscopic world inside the ear



WIFI
connection



6 Bright
Led Lights



360° wide
angle lens



70
minutes

Image: A close-up view of a smartphone screen showing the 1080P HD live feed from the ear camera, illustrating the clarity and detail provided during use.

1. Cleaning the Device

- After each use, carefully remove the silicone ear spoon.
- Clean the silicone ear spoon and the camera lens with a soft, damp cloth or the provided cleaning brush.
- The camera tip is IP67 waterproof, allowing for easy cleaning with water. Ensure the main body of the device is not submerged.
- Allow all parts to dry completely before reassembling or storing.

2. Charging

- Use the provided Type-C charging cable to charge the device.
- Do not use the product while it is charging.
- The device uses 1 Lithium Ion battery (included).

Product Details Display




Image: A diagram illustrating the various parts of the Kolrry Ear Wax Removal Tool Camera, including the charging port, power button, camera, and LED indicator.

TROUBLESHOOTING

- **"Network unavailable" message:** This is normal. The device's Wi-Fi is for direct connection to your phone/tablet only and

does not provide internet access. Proceed with using the app.

- **No image in app:** Ensure the device is powered on (blue light is active) and your phone/tablet is connected to the "Suear-xxxx" Wi-Fi network. Restart both the device and the app if necessary.
- **Unstable connection:** Ensure your device is close to the ear wax removal tool. Avoid interference from other Wi-Fi networks if possible.
- **Ear spoon detachment:** Always ensure the silicone ear spoon is securely tightened onto the otoscope tip before use.

SPECIFICATIONS

- **Brand:** Kolrry
- **Model:** B0F9FK5NV3
- **Material:** Plastic
- **Power Source:** Battery Powered (1 Lithium Ion battery included)
- **Package Dimensions:** 6.73 x 2.6 x 0.75 inches
- **Item Weight:** 2.4 ounces
- **Manufacturer:** Kolrry
- **UPC:** 739768029267

CONTACT AND SUPPORT

For further assistance, product inquiries, or support, please refer to the contact information provided in your product packaging or visit the official Kolrry website.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0F9FK5NV3>

More documents and user manuals: <https://manuals.plus/>