

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

- > [Milavan](#) /
- > [Milavan Smart Health Ring 2.0 User Manual](#)

## Milavan TK20

### Milavan Smart Health Ring 2.0: Your Clear Guide to Wellness



**Smart & Stylish Design**

#### SECTION 1: GETTING STARTED

**STEP 1: FIND YOUR FIT**  
  
Measure Finger Circumference for Perfect Size.

**STEP 2: CHARGE UP**  
  
Place in Case, Charge 1-2 Hours. Up to 75 Days with Case.

**STEP 3: DOWNLOAD & PAIR**  
  
Install 'H RING' App & Follow On-Screen Instructions.

#### SECTION 2: HEALTH MONITORING

**HEART RATE & HRV**  
Automatic Measures Every 30 Mins.

**BLOOD OXYGEN**  
Continuous SpO2 Levels Tracking.

**BLOOD PRESSURE**  
Regular Monitoring Intervals.

**STRESS LEVELS**  
Real-Time Stress Tracking.

#### SECTION 3: ACTIVITY & LIFESTYLE TRACKING & CONVENIENCE

**187 SPORT MODES**  
Select Workout in App, Track Activity & Calories.

**INTELLIGENT SLEEP**  
Detailed Analysis: Deep, Light, REM Sleep & Wake Times.

**LIFESTYLE FEATURES**  
Remote Camera Control, Device Finding, Health Data Sharing, Women's Health Management.

#### SECTION 4: EVERYDAY USE SPECS & CARE

**WATERPROOF**  
IP68 Rated - Showering, Swimming (up to 50m) OK.

**BATTERY LIFE**  
Ring Lasts 3-5 Days, 75 Days with 500mAh Case.

**CLEANING**  
Wipe Gently with Soft, Damp Cloth. Avoid Harsh Chemicals.

Discover the Milavan Smart Health Ring 2.0 with this comprehensive guide covering setup, health monitoring (heart rate, SpO2, blood pressure, stress), activity tracking (187 sport modes, sleep), and everyday use features.

# Milavan Smart Health Ring 2.0 User Manual

Model: TK20



## 1. INTRODUCTION

The Milavan Smart Health Ring 2.0 is a sophisticated fitness tracker designed to monitor various health metrics. This manual provides essential information for setting up, operating, and maintaining your device to ensure optimal performance and accurate health insights.



Image: The Milavan Smart Health Ring 2.0, its charging case, and the mobile application interface.

## 2. WHAT'S IN THE BOX

Upon unboxing your Milavan Smart Health Ring 2.0, verify that all the following items are included:

- Smart Ring x 1
- Charging Case x 1
- USB Type-C charging cable x 1
- Magnetic USB charging cable x 1
- User Manual x 1

## 3. SIZING GUIDE


Accurate sizing is crucial for comfortable wear and precise data collection. It is recommended to use a sizing kit for the best fit. If a sizing kit is not available, follow these steps to measure your finger circumference:

1. Use a thin, inelastic string and wrap it around the finger where you intend to wear the ring.

2. Mark the point where the string overlaps.
3. Measure the length of the string up to the mark to determine your finger's circumference.


Refer to the size chart below to find the corresponding ring size.

## Size Comparison



**Step1:**  
When Measuring, Please Use A Thin, Inelastic String And Wrap It Around Your Finger, And Mark It By A Pen.

**Step2:**  
When Measuring, Please Use A Thin, Inelastic String And Wrap It Around Your Finger, And Mark It By A Pen.



Model	Inner Diameter	Circumference
<b>6#</b>	<b>16.5mm</b>	<b>49-51mm</b>
<b>7#</b>	<b>17.3mm</b>	<b>52-54mm</b>
<b>8#</b>	<b>18.1mm</b>	<b>55-57mm</b>
<b>9#</b>	<b>18.9mm</b>	<b>58-60mm</b>
<b>10#</b>	<b>19.8mm</b>	<b>61-62mm</b>
<b>11#</b>	<b>20.6mm</b>	<b>63-65mm</b>
<b>12#</b>	<b>21.4mm</b>	<b>66-67mm</b>
<b>13#</b>	<b>22.2mm</b>	<b>68-70mm</b>
<b>14#</b>	<b>23.4mm</b>	<b>71-74mm</b>

Image: Size chart for Milavan Smart Health Ring 2.0, detailing inner diameter and circumference measurements for various ring sizes.

## 4. SETUP

### 4.1 Charging the Ring

Before initial use, fully charge your Milavan Smart Health Ring 2.0. Place the ring into its charging case. Connect the charging case to a power source using the provided USB Type-C or magnetic USB charging cable. The ring typically charges in 1-2 hours.

# CHARGING CASE



MILAVAN 1.0

**Upgraded**



MILAVAN 2.0

## **Charging Stability:**

Upgraded by 50%

## **Larger Battery:**

Upgraded from  
400mAh to 500mAh

## **Battery Health:**

+50% Durability

## **Data cable:**

Upgraded with  
Dual Data Cables

Image: The Milavan Smart Health Ring 2.0 charging case, illustrating how the ring is placed for charging.

## **4.2 Downloading the App**

The Milavan Smart Health Ring 2.0 requires the companion "H RING" app for full functionality. The app is compatible with both iOS and Android devices.

- For iOS devices: Search for "H RING" in the Apple App Store.
- For Android devices: Search for "H RING" in the Google Play Store.

Download and install the app on your smartphone.

-  Ultra-light for all-day comfort
-  Crafted with premium stainless steel and resin
-  Gentle on sensitive skin
-  Inner layer specially treated for a smooth, skin-friendly feel

## ALL-DAY COMFORT EXPERIENCE

Milavan 2.0 features a lightweight build and ergonomic design, perfectly blending style with functionality. Whether at work, during workouts, or while relaxing, it ensures natural comfort and adapts effortlessly to all-day wear.

**Upgrade Points**

  
 Upgraded Chip

  
 Enhanced Charging Case

  
 Scratch-Resistant Design

  
 All-Day Comfort

Image: Instructions for downloading the Milavan 2.0 companion app from the App Store and Google Play Store.

### 4.3 Pairing the Ring

After installing the app, open it and follow the on-screen instructions to pair your Milavan Smart Health Ring 2.0 via Bluetooth. Ensure your phone's Bluetooth is enabled and the ring is sufficiently charged.




Image: A visual guide for Step 3 of the setup process, showing the Milavan 2.0 ring being worn and the H RING app icon for download.


## 5. OPERATING THE MILAVAN SMART HEALTH RING 2.0

### 5.1 Health Monitoring

The ring continuously monitors various health parameters using its advanced 5.0 sensor and AB2026B3 chip:

- **Heart Rate & HRV:** Automatically measures heart rate and Heart Rate Variability (HRV) every 30 minutes.
- **Blood Oxygen (SpO2):** Continuously tracks blood oxygen levels.
- **Blood Pressure:** Monitors blood pressure at regular intervals.
- **Stress Levels:** Provides real-time stress level tracking.

All data is synced to the app for comprehensive wellness insights.



## Activity Tracking with 187 Sport Modes

Effortlessly monitor your daily activities with Milavan's 187 sport modes. Simply select your preferred workout from the device, and the ring will track key data including duration, heart rate, calories burned, pace, speed (km/h), cadence, steps, and activity intensity. Gain deeper insights into your exercise performance and achieve better fitness results with ease.

Image: The Milavan 2.0 app interface showing detailed graphs for heart rate and Heart Rate Variability (HRV) monitoring.


## 5.2 Intelligent Sleep Tracking

The Milavan Smart Health Ring 2.0 offers 24-hour automatic sleep monitoring, providing detailed analysis of your sleep stages:

- Deep Sleep
- Light Sleep
- Rapid Eye Movement (REM) stages
- Wake times

This comprehensive tracking helps you understand and improve your sleep quality.

# SLEEP MONITORING



01:34, 06/24

09:25, 06/24



Midday Nap



Nighttime Sleep



Resting on the Grass



Morning Commute

**New: 24-hour** all-day sleep tracking

Old: Fixed sleep time tracking  
(8 PM – 12 PM next day)



Image: The Milavan Smart Health Ring 2.0's sleep monitoring capabilities, illustrating how it tracks sleep patterns across various times of the day.

## 5.3 Sports Modes Tracking

The ring supports 187 sports modes, including running, walking, boxing, ball games, water sports, and ice sports. It automatically records exercise data and syncs it to the app. Activate sport mode within the app for more effective workout guidance.

### Note:

- Turn off sport mode after use to prevent excessive battery drain.



## MULTIPLE SPORTS MODES



Heart Rate



Calories



Distance



Pace



Time

Image: The Milavan Smart Health Ring 2.0 demonstrating its ability to track various sports activities, including running, golf, skiing, climbing, and yoga.

### 5.4 Remote Shooting & Friends Feature

- **Remote Camera Control:** The app supports remote control of your smartphone camera.
- **Device Finding:** Use the app's device-finding function to locate your ring if misplaced.
- **Health Data Sharing:** Share your health data with friends and family through the app, promoting collective wellness.

# KEEP TRACK OF HEALTH

Link family or friends accounts to access health data anytime, anywhere



Image: The Milavan Smart Health Ring 2.0 app interface, illustrating the feature for sharing health data with family and friends.

## 5.5 Women's Health Management

When the female gender option is selected in the app, you can access an exclusive women's health management feature. This allows manual recording of menstrual cycles, pregnancy preparation plans, and pregnancy status, aiding in tracking body changes and managing daily life more effectively.



## WATERPROOF PERFORMANCE

---

With an IP68 rating and up to 50 meters of water resistance, the Milavan 2.0 Smart Ring is built to handle your daily routine effortlessly. Wear it confidently every day while maintaining its outstanding performance.

Image: The Milavan 2.0 app's women's health management feature, showing a calendar interface for tracking menstrual cycles.

## 6. MAINTENANCE

---

### 6.1 Cleaning and Care

The Milavan Smart Health Ring 2.0 features a micro-matte finish designed to resist fingerprints and scratches. To clean the ring, gently wipe it with a soft, damp cloth. Avoid using harsh chemicals or abrasive materials that could damage the finish or sensors.

# REFINED CRAFTSMANSHIP



Stylish Exterior Design



Micro-Matte Finish



Fingerprint-Resistant



Scratch-Resistant

Image: The Milavan Smart Health Ring 2.0, showcasing its micro-matte finish and resistance to scratches and fingerprints.

## 6.2 Waterproof Performance

The ring has an IP68 rating and is water-resistant up to 50 meters, making it suitable for daily routines such as washing hands, showering, and swimming. However, it is not recommended for diving or high-pressure water activities.

## Permanent Data Storage

As long as you don't clear the app's history cache, your data remains accessible after the ring disconnects and reconnects.



The app performs regular firmware and hardware updates to ensure continuous optimization and accurate data tracking.




Image: The Milavan 2.0 Smart Ring demonstrating its waterproof capabilities by being partially submerged in water.

### 6.3 Battery Life

The ring's battery lasts 3-5 days on a single charge. The upgraded charging case, with a 500mAh battery, provides additional charges, extending the total battery life to over 75 days when used together. Charge the ring regularly to ensure continuous tracking.

## 7. TROUBLESHOOTING

### 7.1 Ring Not Connecting to App

- Ensure the ring is charged.
- Verify Bluetooth is enabled on your smartphone.
- Restart the app and your smartphone.
- Try unpairing and re-pairing the device in your phone's Bluetooth settings and the app.

### 7.2 Inaccurate Readings

- Ensure the ring is worn correctly and snugly on your finger.
- Clean the sensors on the inside of the ring.
- Avoid extreme temperatures or rapid changes in environment during measurements.
- Ensure the app is updated to the latest version.

### 7.3 Short Battery Life

- Ensure all sports modes are turned off when not in use.
- Reduce the frequency of continuous monitoring if possible (check app settings).
- Ensure the ring is fully charged each time.
- The charging case should also be charged to provide extended power.

## 8. SPECIFICATIONS

Feature	Detail
Model Number	TK20

Dimensions	2.99 x 2.95 x 1.73 inches (Package)
Item Weight	3.52 ounces
Battery Type	Lithium Polymer (15 mAh)
Charging Case Battery	500 mAh
Operating System Compatibility	Android, iOS
Connectivity Technology	Bluetooth
Water Resistance	IP68 (up to 50 meters)
Special Features	Activity Tracker, Heart Rate, Blood Oxygen, Blood Pressure, Sleep, Stress, 187 Sports Modes

## 9. WARRANTY AND SUPPORT

---

For warranty information and technical support, please refer to the documentation included with your purchase or contact Milavan customer service. You may also visit the official Milavan website for frequently asked questions and support resources.

If purchased through a retailer, please contact the seller directly for specific warranty details and return policies.

