

[Manuals.plus](#) /

› [VERPAL](#) /

› VERPAL Ultra Lightweight Transport Wheelchair Instruction Manual

VERPAL B0FNCDKYX5

VERPAL Ultra Lightweight Transport Wheelchair Instruction Manual

Model: B0FNCDKYX5

1. INTRODUCTION


This manual provides essential information for the safe and effective use of your VERPAL Ultra Lightweight Transport Wheelchair. Please read it thoroughly before initial use and retain it for future reference. This wheelchair is designed to provide portable mobility assistance for individuals requiring transport.

2. IMPORTANT SAFETY PRECAUTIONS

- **Maximum Weight Capacity:** Do not exceed the maximum weight capacity of 220 lbs (100 kg).
- **Brakes:** Always engage the parking brakes when the wheelchair is stationary, especially when transferring the user in or out.
- **Seatbelt:** Always ensure the user is securely fastened with the built-in seatbelt during transport.
- **Supervision:** This is a transport wheelchair and requires an attendant to push and maneuver it. Never leave a user unattended in the wheelchair.
- **Uneven Surfaces:** Exercise extreme caution on uneven surfaces, inclines, or declines. Avoid sudden stops or turns.
- **Maintenance:** Regularly inspect all parts for wear, damage, or loose components. Do not use if any part is damaged.
- **Transferring:** Ensure the wheelchair is stable and brakes are engaged before assisting a user with transfers.



Brake: Hold Up Brake, Fix and Park Wheelchair



Unlock: Press The Trigger Release The Brake




Image: Close-up view of the brake handles and triggers. The top panel shows the 'Brake Handle' and 'Trigger' labels. The bottom left panel illustrates how to 'Hold Up Brake, Fix and Park Wheelchair' by pulling the levers upwards. The bottom right panel shows how to 'Press The Trigger Release The Brake' by squeezing the red triggers.

3. UNFOLDING AND SETUP

To prepare your VERPAL Transport Wheelchair for use, follow these steps:

1. Place the folded wheelchair on a flat, stable surface.
2. Gently pull the backrest upwards while simultaneously pushing the seat downwards until the frame locks into place.
3. Ensure the automatic locking buckle is engaged. You should hear a "click" sound.
4. Lower the footrests into position.
5. Adjust the armrests if necessary.

Automatic Locking Buckle

Never worry about forgetting to secure it. Our buckle locks automatically, preventing loosening and ensuring a safe ride, every time.



Image: A diagram illustrating the automatic locking buckle. The main image shows the wheelchair partially unfolded with arrows indicating 'Unfold' and 'Down'. An inset circle provides a magnified view of the orange locking buckle, highlighting its mechanism for securing the frame in the open position.

Video: Unfolding and Setup Demonstration

Video: This official seller video demonstrates the unfolding process of the VERPAL Ultra Lightweight Transport Wheelchair, showing how to prepare it for use.

4. OPERATING INSTRUCTIONS

Brakes

The wheelchair is equipped with reliable brakes integrated into the push handles. To engage the parking brake, press down firmly on the brake levers until you hear a click. To release the brake, squeeze the red triggers on the brake handles.

Seatbelt

For the user's safety, always fasten the built-in seatbelt. Insert the male buckle end into the female buckle receiver until it clicks securely. To release, press the button on the female buckle receiver.

Comfortable Breathable




33cm/13inch Long Armrests for Extra Comfort



Built-in Foam



Seat Belt Design



Ergonomic Design



Image: This image highlights key comfort and safety features. It shows the 'Comfortable Breathable' seat and backrest, '33cm/13inch Long Armrests for Extra Comfort', 'Built-in Foam' for cushioning, 'Seat Belt Design' with a hand fastening the buckle, and an 'Ergonomic Design' illustration.

Armrests and Footrests

The armrests are rotatable, waterproof, and non-slip, designed for user comfort and ease of transfer. The footrests are foldable, allowing for easier entry and exit from the wheelchair.

5. FOLDING AND PORTABILITY

The VERPAL Ultra Lightweight Transport Wheelchair is designed for quick and easy folding, making it highly portable and convenient for travel and storage.

1. Ensure the front wheels are aligned parallel and facing forward.

2. Step on the pedal located near the rear wheels, press the buckle, and simultaneously hold the back handle while folding the backrest down.
3. Fold up the handrails and footrests to achieve a compact folded size.
4. Utilize the retractable handle for easy pulling and transport, similar to a suitcase.

Lightweight - 3-Step Folding

1

Align front wheels parallel and forward.



2

Step on pedal, press buckle and hold back handle, folding backrest down simultaneously.



3

Fold up the handrails and footrests.



Image: A visual guide titled 'Lightweight - 3-Step Folding'. Step 1 shows aligning front wheels. Step 2 illustrates pressing a pedal and folding the backrest. Step 3 shows folding the handrails and footrests. A larger image on the right depicts a person pulling the folded wheelchair with its telescopic handle.



Image: An infographic titled 'Ultra Light 22.4 LBS'. It shows the wheelchair in its open state and then in three smaller panels: 'Foldable' showing the compact folded form, 'Quality Handle' demonstrating the retractable handle being pulled, and 'Easy to Carry' with the folded wheelchair stored in a car trunk.

Video: Folding and Portability Demonstration

Video: This official seller video demonstrates the folding process of the VERPAL Ultra Lightweight Transport Wheelchair, showcasing its compact storage and portability features.

6. MAINTENANCE

Proper maintenance ensures the longevity and safe operation of your wheelchair.

- **Cleaning:** Clean the frame with a damp cloth and mild soap. Avoid abrasive cleaners. The fabric seat and backrest can be wiped clean.
- **Wheel Inspection:** Regularly check the wheels for debris, wear, and proper rotation. Ensure all wheel fasteners are tight.
- **Brake Check:** Test the brakes regularly to ensure they engage and hold securely.
- **Frame Inspection:** Inspect the carbon steel frame for any signs of damage, cracks, or bending.

- **Fasteners:** Periodically check all screws, bolts, and other fasteners to ensure they are secure.
- **Storage:** Store the wheelchair in a dry, clean place away from extreme temperatures.

7. TROUBLESHOOTING

If you encounter any issues with your wheelchair, refer to the following common problems and solutions:

- **Wheelchair does not unfold or fold smoothly:** Ensure all locking mechanisms are disengaged before folding/unfolding. Check for any obstructions or debris in the folding joints.
- **Brakes are not holding:** Inspect the brake cables for damage or looseness. Ensure the brake levers are fully engaged. If issues persist, contact customer support.
- **Wheels are not rolling freely:** Check for hair, lint, or other debris wrapped around the wheel axles. Ensure the wheels are properly inflated (if applicable, though these are solid PVC wheels).
- **Unusual noises during use:** Inspect all moving parts for loose connections or wear. Lubricate joints if necessary with a silicone-based lubricant.

8. SPECIFICATIONS

Feature	Specification
Product Dimensions (D x W x H)	33"D x 20"W x 35"H
Seat Width	16 Inches
Seat Depth	15 inches
Wheel Size	Front: 6 Inches, Rear: 7.5 Inches
Handle Height	10.2 Inches
Maximum Weight Recommendation	220 Pounds
Frame Material	Carbon Steel
Wheel Type	Honeycomb Wheel (PVC)
Brake Style	Push-to-Lock
Color	Black

9. WARRANTY AND SUPPORT

Warranty Information

This VERPAL Ultra Lightweight Transport Wheelchair comes with a standard manufacturer's warranty covering defects in materials and workmanship under normal use. Please refer to the product packaging or contact customer support for specific warranty terms and duration.

Customer Support

For any questions, assistance, or to report issues, please contact VERPAL customer support through the retailer

where the product was purchased or visit the official VERPAL website for contact details.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0FNCDKYX5>

More documents and user manuals: <https://manuals.plus/>