

[Manuals.plus](#) /

› [Hovadova](#) /

› Hovadova HV-JLSMQDZ Pneumatic Height Adjustable Overbed Table Instruction Manual

Hovadova HV-JLSMQDZ

Hovadova HV-JLSMQDZ Pneumatic Height Adjustable Overbed Table Instruction Manual

MODEL: HV-JLSMQDZ

Introduction

This manual provides essential information for the safe and efficient use of your Hovadova HV-JLSMQDZ Pneumatic Height Adjustable Overbed Table. Please read these instructions carefully before assembly and operation to ensure proper functionality and longevity of the product.

Product Features

- **Gas Spring Pneumatic Lift System:** Allows for easy, one-handed height adjustment from 28.74 inches to 43.5 inches.
- **Tilt Panel Design:** Desktop can be adjusted from 0° to 70° to suit various activities like reading, writing, or laptop use.
- **Multi-Functional Use:** Suitable as a bedside table, sofa side table, medical tray, standing desk, or portable workstation.
- **Easy to Move and Space-Saving:** Features 4 hidden casters for mobility and an H-shaped base for stability. The desktop folds for compact storage.
- **Durable Construction:** Heavy-duty welded tubular steel structure ensures stability and long-term use.

PRODUCT DIMENSIONS

Take up a little space

Load capacity: **88 lbs**



Image: Product dimensions and load capacity of the Hovadova overbed table.

Alternate Sitting & Standing

Height ranges from 28.74 inches to 43.5 inches



Image: The overbed table adjusted for both sitting and standing use.

Ergonomic Tilting Desktop Design

Stop at any angle with in 0-70°



0 - 30°
Writing/Working



30 - 50°
Reading



50 - 70°
Watching Video

Image: Ergonomic tilting desktop design for various activities.

Saving Space

90° foldable tabletop, one-touch pull-ring storage, saves space when not in use



Image: Space-saving design with a 90° foldable tabletop.

Product Details



Desktop Latch



Pneumatic Tilting Desk



Anti-Slip Guard



Cup Holder

Image: Detailed view of the desktop latch, pneumatic tilting mechanism, anti-slip guard, and cup holder.



Easy to Move

Suitable for bedside, sofa edge, floor, and blanket use



Hidden casters



Visible Wheels

Image: Hidden casters for easy and smooth mobility.

Hospital Bed Table



Image: Overbed table used in a hospital setting.

Setup and Assembly

The Hovadova overbed table is designed for quick and easy assembly. All necessary tools are included in the package. One person can typically assemble the table within 10 minutes by following the provided instructions.

Assembly Steps:

1. Unpack all components and verify against the parts list in your manual.
2. Attach the base frame to the main vertical support column using the provided screws and tools.
3. Secure the tabletop assembly to the top of the vertical support column.
4. Connect the pneumatic lift mechanism and ensure all connections are secure.
5. Test the height adjustment and tilt functions to confirm proper assembly.

Operating Instructions

Height Adjustment:

To adjust the table height, locate the pneumatic lift lever (typically on the side of the main column). Press or pull the lever and gently guide the tabletop up or down to your desired height. Release the lever to lock the table in position. The table can be adjusted from 28.74 inches to 43.5 inches.

Desktop Tilt Adjustment:

The main desktop panel features a tilting mechanism. To adjust the tilt, release the desktop lock (refer to product details image for location). Gently tilt the desktop to your preferred angle (0° to 70°). The pneumatic system will hold it in place. Ensure the anti-slip stopper is engaged to prevent items from sliding off.

Moving the Table:

The table is equipped with hidden casters for easy mobility. Simply push or pull the table to move it to a new location. The low-profile H-shaped base allows it to slide easily under beds, sofas, or other furniture.

Folding for Storage:

To fold the tabletop for compact storage, pull the cord outward to release the tabletop from its horizontal position. The tabletop will slowly descend to a vertical position, saving space when not in use.

Video: Demonstration of the Hovadova Pneumatic Height Adjustable Overbed Table with Wheels, showcasing its features and adjustments.

Maintenance

To maintain your overbed table, regularly wipe surfaces with a soft, damp cloth. Avoid abrasive cleaners or harsh chemicals that may damage the finish. Ensure all screws and connections remain tight for stability. Periodically check the pneumatic lift system for smooth operation.

Safety Precaution:

Do not exceed the maximum weight recommendation of 88 pounds to ensure stability and prevent accidents.

Troubleshooting

- **Table not adjusting smoothly:** Check for any obstructions in the lift mechanism or ensure the lever is fully engaged/disengaged.
- **Desktop not tilting:** Verify that the desktop lock is fully released.
- **Table feels unstable:** Ensure all assembly screws are tightened. Avoid exceeding the maximum weight capacity.

If you encounter persistent issues, please refer to the warranty and support section for assistance.

Specifications

Feature	Detail
Brand	Hovadova
Model Number	HV-JLSMQDZ
Product Dimensions (D x W x H)	14.17"D x 27.48"W x 41.06"H
Item Weight	25 Pounds
Maximum Weight Recommendation	88 Pounds
Height Adjustment Range	28.74 inches to 43.5 inches
Desktop Tilt Range	0° to 70°
Frame Material	Metal
Top Material Type	Engineered Wood
Color	Brown
Included Components	Installation Tool

Warranty and Support

For warranty information or technical support, please contact Hovadova-US directly. Refer to your product packaging or the official Hovadova website for the most current contact details and warranty terms.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0FQTHVJ57>

More documents and user manuals: <https://manuals.plus/>