

VERPAL B0FQTPF1XB

VERPAL Electric Wheelchair User Manual

MODEL: B0FQTPF1XB

Brand: VERPAL

1. Introduction

This manual provides essential information for the safe and effective operation, maintenance, and troubleshooting of your VERPAL Electric Wheelchair. Please read this manual thoroughly before using the wheelchair to ensure proper function and user safety. Keep this manual for future reference.

2. Safety Information

- Always read and understand all instructions before operating the wheelchair.
- Ensure the user's weight does not exceed the maximum capacity of 265 pounds.
- Always wear the seat belt provided for enhanced safety during operation.
- Familiarize yourself with the electromagnetic braking system for immediate stopping.
- Exercise caution on inclines and uneven terrain. The wheelchair can climb inclines up to 9 degrees.
- Avoid operating the wheelchair in extreme weather conditions or on slippery surfaces.
- Regularly inspect the durable, puncture-resistant tires for wear and tear.
- Ensure all components are securely fastened before each use.

3. Product Overview

The VERPAL Electric Wheelchair is designed for adults and seniors, offering powerful dual motors and all-terrain capabilities. Key features include:

- **500W Dual Motors:** Provides powerful propulsion for various terrains.
- **16-inch All-Terrain Wheels:** Offers excellent grip and stability.
- **Electromagnetic Braking System:** Ensures immediate and secure stopping.

- **360° Joystick Control:** Allows for effortless and precise navigation.
- **Foldable Design:** Compact for easy transport and storage.
- **Long Travel Range:** Up to 15 miles on a single charge.
- **Durable Construction:** Made from aerospace-grade alloys with a robust four-layer frame.



Figure 3.1: Powerful 500W dual motors for efficient movement.

**Folds up easily
without taking up space**



Compatible with sedan; SUV; MPV two-seater

Figure 3.2: Detailed view of the seat belt, storage bag, armrest, and seat cushion.

360° Control Joystick



Intelligent
Control



Immediately
Master



Small Turning
Radius

Figure 3.3: The 360-degree control joystick for intuitive navigation.

4. Setup and Assembly

Your VERPAL Electric Wheelchair comes largely pre-assembled. The package includes the power wheelchair, user manual, tool kit, charger, and storage bags. Follow these steps for final setup:

1. Carefully remove the wheelchair and all accessories from the packaging.

2. Unfold the wheelchair by gently pulling the backrest upwards until it locks into place.
3. Attach the joystick controller to the desired armrest (left or right) and secure it.
4. Connect the battery pack securely. Ensure all electrical connections are firm.
5. Adjust the footrest to a comfortable position.

For a visual guide on assembly and setup, please refer to the following videos:

Your browser does not support the video tag.

Video 4.1: Detailed steps to install the electric wheelchair.

Your browser does not support the video tag.

Video 4.2: Visual guide on the installation process of the electric wheelchair.

5. Operating Instructions

Operating your VERPAL Electric Wheelchair is intuitive with its responsive joystick control.

5.1. Power On/Off

- Press the power button on the joystick to turn the wheelchair on or off.
- The battery indicator lights will show the current charge level.

5.2. Driving

- Use the 360° joystick to control direction and speed. Push forward to move forward, pull back to move backward, and move left or right to steer.
- Adjust speed settings using the acceleration (+) and deceleration (-) buttons on the joystick.
- The electromagnetic braking system engages automatically when the joystick is released, ensuring a smooth stop.

5.3. Folding and Transport

The wheelchair features a compact, foldable design for easy transport.

- Ensure the power is off before folding.
- Follow the folding mechanism instructions (typically involves releasing a latch and collapsing the frame).
- Once folded, the wheelchair can easily fit into a car trunk.



Figure 5.1: The electric wheelchair in its folded state, ready for transport.

Your browser does not support the video tag.

Video 5.2: Demonstration of the comfortable and easy folding mechanism.

5.4. Charging

- Connect the charger to the charging port on the battery pack.
- Plug the charger into a standard electrical outlet.
- The charging indicator on the joystick or battery will show charging status.
- A full charge provides up to 15 miles of travel range.

6. Maintenance

Regular maintenance ensures the longevity and optimal performance of your wheelchair.

- **Cleaning:** Wipe down the frame and seating surfaces with a damp cloth. Avoid harsh chemicals.
- **Tires:** Periodically check the puncture-resistant tires for debris or damage.
- **Battery:** Keep the battery charged. Avoid fully discharging the battery frequently to prolong its lifespan. Store the battery in a cool, dry place if not in use for extended periods.
- **Frame:** Inspect the aerospace-grade alloy frame for any signs of damage or loose connections.
- **Joystick:** Keep the joystick clean and free from obstructions.

7. Troubleshooting

If you encounter any issues with your wheelchair, refer to the following common problems and solutions:

- **Wheelchair does not power on:**
 - Check if the battery is properly connected and charged.
 - Ensure the power button on the joystick is pressed firmly.
- **Wheelchair does not move:**
 - Verify that the joystick is fully engaged in a direction.
 - Check for any obstructions around the wheels.
 - Ensure the manual brake is disengaged and clutches are pressed down on both sides for power operation.
- **Unusual noises during operation:**
 - Inspect wheels and motors for any foreign objects.
 - Contact customer support if the noise persists.

For more detailed troubleshooting, please consult the full user manual included with your product or contact customer support.

8. Specifications

Feature	Detail
Product Dimensions	38.2"D x 26.4"W x 37.8"H
Item Weight	66 Pounds
Weight Capacity	265 Pounds
Motor Power	500W Dual Motors
Wheel Size	16-inch All-Terrain Wheels
Max Incline	Up to 9°
Travel Range	Up to 15 miles
Folded Dimensions	29"D x 16"W x 28"H

Obstacle Clearance Height	1.6 inches
---------------------------	------------

9. Warranty and Support

Your VERPAL Electric Wheelchair comes with a 30-day return policy for refund or replacement. For any questions, concerns, or technical assistance, please contact VERPAL customer support. The package includes a tool kit and charger for your convenience.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question