

ELENKER HFK-9269KD

ELENKER HFK-9269KD 2-in-1 All Terrain Rollator Walker & Transport Chair Instruction Manual

Model: HFK-9269KD

1. INTRODUCTION

Thank you for choosing the ELENKER HFK-9269KD 2-in-1 All Terrain Rollator Walker and Transport Chair. This product is designed to provide mobility assistance, functioning both as a rollator walker for independent use and a transport chair for assisted mobility. Please read this manual thoroughly before assembly and use to ensure safe operation and proper maintenance.

2. SAFETY INFORMATION

- **Weight Capacity:** Do not exceed the maximum weight capacity of 300 pounds (136 kg).
- Ensure all parts are securely assembled and locked before each use.
- Always engage the brakes when sitting on the seat or when the unit is stationary.
- When using as a transport chair, ensure the user is secured with the safety belt and footrests are properly positioned. An attendant must always be present to push the chair.
- Avoid using on uneven or slippery surfaces without proper caution.
- Regularly inspect the brakes, wheels, and frame for any signs of wear or damage.

3. COMPONENTS LIST

The ELENKER HFK-9269KD includes the following main components:

- Main Frame with Wheels
- Padded Seat
- Adjustable Handlebars with Brakes
- Extra Wide Backrest
- Removable Footrests (2)
- Under-seat Storage Bag

- Safety Belt
- Cane/Umbrella Holder

4. SETUP AND ASSEMBLY

Follow these steps to set up your rollator walker and transport chair:

1. **Unfolding the Unit:** Place the unit on a flat surface. Gently pull the handlebars outwards until the frame locks into its open position. Ensure the seat is flat and secure.
2. **Attaching Footrests:** Insert the footrests into the designated slots on the front frame. Push down until they click securely into place. Ensure they are firmly attached before use.
3. **Adjusting Handlebar Height:** Locate the adjustment pins or levers on the handlebar posts. Press the pin or open the lever and slide the handlebars up or down to a comfortable height. Ensure the handles are at wrist level when standing upright with arms relaxed. Release the pin/close the lever to lock the handlebars in place. Verify both handlebars are at the same height and securely locked.
4. **Checking Brakes:** Squeeze the brake levers to ensure they engage the wheels. Push the brake levers down to lock the wheels for parking. Test the brakes before initial use.

The Height of Handlebars is Adjustable

For a more comfortable and natural posture while walking



Figure 4.1: Adjustable Handlebar Height and Product Dimensions. The handlebars can be adjusted from 35 to 41 inches in height for user comfort.

5. OPERATING INSTRUCTIONS

5.1. Using as a Rollator Walker

- Stand behind the rollator, holding both handlebars.
- Push the rollator forward a short distance, then step forward. Maintain a comfortable distance from the rollator.
- To rest, engage the parking brakes by pushing the brake levers down. Sit on the padded seat, ensuring the backrest is in position.



Figure 5.1: Rollator Walker Mode. The unit features a cup holder, under-seat basket, and cane/umbrella holder for convenience.

5.2. Using as a Transport Chair

- Ensure the footrests are attached and swung into position for the user's feet.
- Engage the parking brakes.
- Assist the user to sit on the padded seat. Secure the user with the safety belt.
- Disengage the parking brakes. An attendant should push the chair using the handlebars.

Multi-function Rollator Walker

Converts easily from a rollator to a transport wheelchair




Figure 5.2: Transport Chair Mode. The safety belt, armrests, and footrests provide comfort and security for the user.



Figure 5.3: Safety Belt Detail. Always ensure the safety belt is fastened when using the unit as a transport chair.

5.3. Braking System

- **To Slow Down:** Squeeze the brake levers gently towards the handlebars.
- **To Lock Wheels (Parking Brake):** Push the brake levers down until they click into the locked position. Pull them up to release.



Figure 5.4: Dual Braking System. Squeeze to slow, push down to lock.

5.4. Folding for Storage/Transport

- Remove any items from the storage bag.
- Lift the seat handle or strap to initiate the folding mechanism.
- Push the sides of the frame together until the unit is in a compact, folded position.
- Secure with any provided locking mechanisms or straps to prevent accidental unfolding during transport.



Figure 5.5: Folded Unit. The rollator folds for convenient storage and transport.

6. MAINTENANCE

- **Cleaning:** Wipe the frame with a damp cloth and mild detergent. Do not use abrasive cleaners. The padded seat and backrest can be spot cleaned.
- **Wheel Inspection:** Regularly check the 10-inch non-pneumatic wheels for debris, wear, or damage. Ensure they rotate freely.
- **Brake Adjustment:** If brakes feel loose or do not engage properly, consult a qualified technician for adjustment.
- **Storage:** Store the unit in a dry, clean place away from direct sunlight and extreme temperatures.



Figure 6.1: 10-inch Non-pneumatic Wheels. These wheels require no inflation and are designed for smooth navigation.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Brakes not engaging effectively	Brake cable too loose; worn brake pads	Consult a qualified technician for brake adjustment or replacement.
Wheels not rolling smoothly	Debris in wheels; damaged bearings	Clean wheels of any debris. If problem persists, contact support.
Unit does not fold completely	Obstruction; mechanism not fully released	Check for obstructions. Ensure seat handle is fully lifted.
Handlebars feel loose after adjustment	Adjustment pins/levers not fully engaged	Ensure adjustment pins are fully extended into the holes or levers are securely closed.

8. SPECIFICATIONS

- **Model Number:** HFK-9269KD

- **Product Dimensions (L x W x H):** 29 x 25 x 41 inches
- **Item Weight:** 33 pounds
- **Maximum Weight Capacity:** 300 pounds
- **Adjustable Handlebar Height:** 35 - 41 inches
- **Seat Height:** 22 inches
- **Wheel Size:** 10 inches (Non-pneumatic)
- **Material:** Alloy Steel, Metal
- **Color:** Black

9. WARRANTY INFORMATION

ELENKER products are manufactured to high-quality standards. This product comes with a manufacturer's warranty against defects in materials and workmanship. Please refer to the warranty card included with your purchase or contact ELENKER customer support for specific warranty terms and conditions.

10. CUSTOMER SUPPORT

For any questions, concerns, or assistance with your ELENKER HFK-9269KD rollator walker and transport chair, please contact our customer support team. Contact information can typically be found on the ELENKER official website or on the packaging materials.