

[Manuals.plus](#) /

› [COTAO](#) /

› COTAO Aurivo ReHears Hearing Red Light Therapy User Manual

COTAO SPK-700

COTAO Aurivo ReHears Hearing Red Light Therapy User Manual

Model: SPK-700

INTRODUCTION

The COTÃO Aurivo ReHears device is a red light therapy tool designed for both nasal and ear application. It utilizes a 650nm wavelength for targeted therapy. Its design aims for comfortable and unrestricted use during daily activities.

PRODUCT OVERVIEW

The Aurivo ReHears device is compact and offers dual functionality. Key features include:

- **Dual Application:** Supports both nasal and ear therapy.
- **650nm Wavelength:** Emits red light at 650nm.
- **Interchangeable Silicone Tips:** Includes soft, washable tips in three sizes (S/M/L) for ear use, and a reusable nasal cover.
- **Ergonomic Design:** Lightweight and designed for comfortable wear around the neck.
- **Adjustable Light Modes:** Provides options for slow flashing, fast flashing, or constant light.



Figure 1: The Aurivo ReHears device, showcasing its earbud-like design with red light emitting from the tips.

SETUP

1. **Unpack Contents:** Carefully remove the Aurivo ReHears device, USB cable, various silicone ear tips (S/M/L), and nasal cover from the carrying case.
2. **Select Appropriate Tips:** Choose the silicone ear tips or nasal cover that best fits your needs. The device includes options for both ear and nasal application.
3. **Attach Tips:** Gently attach the selected silicone tips or nasal cover to the red light emitters on the device.
4. **Connect Power:** The device features a USB connector with a removable cap revealing a USB-C connector. Plug the device into a compatible USB power source (e.g., phone charger, computer USB port, or power bank). The red lights will illuminate upon connection.




Figure 2: Detailed view of the USB-A and USB-C connectors for flexible power options.

OPERATING INSTRUCTIONS

The Aurivo ReHears device is designed for simple and effective use. Follow these steps for optimal results:

1. **Prepare Device:** Ensure the correct silicone tips or nasal cover are securely attached and the device is connected to a power source. The red lights should be active.
2. **Positioning for Ear Therapy:** Gently insert the ear tips into your ears. The device is designed to fit comfortably around the neck, allowing for hands-free use.
3. **Positioning for Nasal Therapy:** Insert the nasal cover into your nostrils.
4. **Select Light Mode:** The device offers three light modes:
 - **Slow Flashing:** For gentle stimulation.
 - **Fast Flashing:** For more intense stimulation.
 - **Constant Light:** For continuous therapy.

Refer to the included manual for specific instructions on how to cycle through these modes, if applicable.

5. **Session Duration:** Use the device for recommended durations as specified in the full product manual. The product description suggests daily wear for optimal benefits.

Your browser does not support the video tag.

Video 1: Demonstration of the Aurivo ReHears device, showing connection to a power source, activation of red lights, and placement in the ear for therapy.



Figure 3: A user demonstrating the ear therapy application of the Aurivo ReHears device.



Figure 4: The device in use, highlighting its comfortable and discreet design.



Figure 5: An illustration depicting the 650nm red light targeting inner ear cells.

MAINTENANCE AND CARE

- **Cleaning:** The silicone ear tips and nasal cover are washable. It is recommended to wipe them with alcohol before and after each use to ensure hygiene.
- **Storage:** Store the device and its accessories in the provided carrying case when not in use to protect it from dust and damage.
- **General Care:** Avoid exposing the device to extreme temperatures, direct sunlight, or moisture. Do not attempt to disassemble the main unit.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device not turning on / Red lights not illuminating.	Not properly connected to power source; faulty USB port/cable.	Ensure the USB cable is securely plugged into both the device and a working power source. Try a different USB port or cable.
Silicone tips do not fit comfortably.	Incorrect tip size selected.	Try a different size of the provided silicone ear tips (S/M/L) or ensure the nasal cover is correctly positioned.

SPECIFICATIONS

Feature	Detail
Brand	COTAO
Model Number	SPK-700
Wavelength	650nm Red Light
Material	Silicone (tips)
Connectivity	USB-A, USB-C (via adapter)
Dimensions	10 x 10 x 10 cm (Parcel Dimensions)
Weight	200 g (Parcel Weight)
Reusability	Reusable tips
Country of Origin	China

SAFETY INFORMATION

- Always clean silicone tips with alcohol before and after use, especially if sharing the device.
- Do not look directly into the red light emitters.
- Keep out of reach of children.
- Consult a healthcare professional before use if you have any pre-existing medical conditions or concerns.
- Discontinue use if you experience any discomfort or adverse reactions.
- This device is not intended to diagnose, treat, cure, or prevent any disease.

WARRANTY AND SUPPORT

For warranty information, product support, or to inquire about replacement parts, please contact COTÃO customer service directly or refer to the contact information provided in the product packaging. Keep your purchase receipt for warranty claims.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How long am I supposed to keep the red lights in my ears when using the COTAO Aurivo ReHears (SPK-700)?

The manual for the COTAO Aurivo ReHears (SPK-700) does not specify a strict number of minutes per session, but it provides the following guidance for usage:Frequency: The product description suggests daily wear for optimal...

[Read full answer](#)

When the COTAO Aurivo ReHears SPK-700 is charging, does the red light flicker the same way it does during ear therapy?

Based on the product manual for the COTAO Aurivo ReHears (SPK-700), the device does not have an internal battery and therefore does not "charge" in the traditional sense. Instead, it must be plugged into...

[Read full answer](#)

Can you order extra silicone covers for the ears and nasal attachments for the COTAO Aurivo ReHears (SPK-700)? If so, how do I go about ordering them?

Yes, replacement parts for the COTAO Aurivo ReHears (SPK-700) are available, though they are typically handled directly by the manufacturer rather than through general retail listings.To order extra silicone ear tips (available in S,...

[Read full answer](#)

I did not receive a manual with my COTAO Aurivo ReHears SPK-700; can you provide one?

While Manuals+ is an independent repository and cannot physically mail a paper manual to you, you can access the full digital version of the COTAO Aurivo ReHears (SPK-700) manual directly on this page.Key Operating...

[Read full answer](#)

How long does it take to charge the COTAO Aurivo ReHears (SPK-700) and how do I know when it is fully charged?

Based on the available manual for the COTAO Aurivo ReHears (SPK-700), the device does not appear to have an internal battery that requires charging. Instead, it is designed to be powered directly via a...

[Read full answer](#)