

[Manuals.plus](#) /

> [VOCIC](#) /

> VOCIC Electric Seat Lift AX16 PRO Instruction Manual

VOCIC AX16 PRO

VOCIC Electric Seat Lift AX16 PRO Instruction Manual

Model: AX16 PRO

INTRODUCTION

The VOCIC Electric Seat Lift AX16 PRO is designed to provide assistance for individuals transitioning from a sitting to a standing position. This device features a dual remote system, offering both rechargeable cordless and wired operation for convenience and flexibility. Its robust construction supports up to 330 lbs, and the ergonomic design ensures comfort during extended use. The patented 15° Golden Arc Technology aims to minimize knee strain and enhance stability during the lift.




Figure 1: Before and after using the VOCIC Electric Seat Lift, demonstrating ease of standing and fall prevention.

SAFETY INFORMATION

- Do not exceed the maximum weight recommendation of 330 pounds.
- Be aware of tripping hazards when using the cordless option; ensure cables are managed safely.
- Always ensure the seat lift is placed on a stable, flat surface or chair.
- Do not use the seat lift if any part is damaged or appears to be malfunctioning.
- Keep children and pets away from the device during operation.

PRODUCT OVERVIEW

The AX16 PRO Electric Seat Lift is designed for portability and ease of use. It includes a main cushion unit, detachable armrests, and two types of controllers: a rechargeable wireless remote and a wired remote. The unit is compact and lightweight, making it suitable for various indoor settings.

Components Included:



Figure 2: Included components of the VOCIC Electric Seat Lift AX16 PRO.

- **AX16 Main Unit:** The primary seat cushion and lifting mechanism.
- **Armrests:** Two detachable armrests for added support.
- **Rechargeable Remote:** For cordless operation.
- **Wired Remote:** For continuous power operation.
- **Power Adapter:** For charging the rechargeable remote and powering the wired remote.
- **Hooks:** For cable management or securing the remote.

SETUP

The VOCIC Electric Seat Lift AX16 PRO requires minimal setup and is ready to use out of the box. No complex assembly is needed.

1. **Unpack:** Carefully remove all components from the packaging.
2. **Placement:** Place the main AX16 unit onto a stable chair, sofa, or other seating surface. Ensure the surface is firm and level.
3. **Attach Armrests:** Insert the armrests into the designated slots on the sides of the main unit until they click securely into place.
4. **Connect Remote:**
 - For **Wired Operation:** Connect the wired remote to the main unit and plug the power adapter into a wall outlet.
 - For **Cordless Operation:** Ensure the rechargeable remote is charged. Connect the rechargeable remote to the main unit.
5. **Initial Charge (for rechargeable remote):** Before first use, fully charge the rechargeable remote using the provided power adapter.



Figure 3: The seat lift is designed for immediate use without complex installation.

OPERATING INSTRUCTIONS

The VOCIC Electric Seat Lift AX16 PRO offers simple, intuitive operation using either the wired or rechargeable remote control.

Using the Remote Control:

- **Sitting Down:** Gently sit on the cushion. The lift mechanism will remain in its lowered position.
- **Standing Up:**
 - Hold the remote control firmly.
 - Press and hold the "UP" button (usually indicated by an upward arrow).
 - The seat will slowly and smoothly begin to rise, assisting you to a standing position.
 - Release the button once you have reached a comfortable standing height or are fully standing.
- **Lowering the Seat:**
 - While standing, press and hold the "DOWN" button (usually indicated by a downward arrow).
 - The seat will slowly lower to its resting position.
 - Release the button when the seat is fully lowered or at your desired height.



Figure 4: The seat lift offers both rechargeable cordless and wired plug-in operation.




Figure 5: The 15° Golden Arc Technology provides a balanced and secure lift.

Official Product Video:

Your browser does not support the video tag.

Video 1: Official VOCIC US demonstration of the Electric Seat Lift AX16, showcasing its features and operation.

MAINTENANCE

- **Cleaning:** Wipe the cushion and frame with a damp cloth and mild soap. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** Store the seat lift in a dry place away from direct sunlight when not in use for extended periods.
- **Battery Care (Rechargeable Remote):** To prolong battery life, avoid fully discharging the remote frequently. Charge as needed.
- **Inspection:** Periodically check all connections, cables, and moving parts for any signs of wear or damage. If any issues are found, discontinue use and contact customer support.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Unit does not respond to remote.	Remote battery is low or depleted (rechargeable). Power adapter not connected or faulty (wired). Loose connection between remote and main unit.	Charge the rechargeable remote. Ensure power adapter is securely plugged into both the unit and a working wall outlet. Check all cable connections.
Lift mechanism is jerky or noisy.	Exceeding weight capacity. Unit not placed on a stable surface. Internal mechanism issue.	Ensure user weight is within the 330 lbs limit. Place the unit on a firm, level, and stable chair or surface. Discontinue use and contact customer support if the issue persists.
Armrests feel loose.	Armrests not fully engaged.	Press armrests firmly into slots until they click.

SPECIFICATIONS



Figure 6: The robust Electric Stabilized Rail-Fork System supports significant weight.

Feature	Detail
Brand	VOCIC
Model Name	AX16 PRO
Material	Carbon Steel (Frame), Foam (Seat)
Product Dimensions (D x W x H)	17.6" x 21.4" x 3.7"
Item Weight	16.5 pounds
Maximum Weight Recommendation	330 pounds
Motor Type	5000N Brushless Motor
Noise Level	Less than 50dB (Whisper-quiet operation)
Lift Angle Technology	Patented 15° Golden Arc Technology
Control Type	Rechargeable Wireless Remote, Wired Remote
Indoor/Outdoor Usage	Indoor

WARRANTY AND SUPPORT

VOCIC is committed to providing quality products and customer satisfaction. The AX16 PRO Electric Seat Lift comes with comprehensive support.

- **Lifetime Support:** VOCIC offers lifetime support for this product.
- **Professional Customer Service:** Dedicated customer service is available to assist with any inquiries or issues.
- **30-Day Free Return:** Refer to your purchase agreement for details on returns.
- **CE Certification:** This product meets European safety, health, and environmental protection requirements.

For technical assistance, troubleshooting, or warranty claims, please contact VOCIC US customer service directly. Contact information can typically be found on the product packaging or the official VOCIC website.



Figure 7: VOCIC's commitment to customer satisfaction and product quality.

© 2026 VOCIC. All rights reserved. Model AX16 PRO.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0G7YZFF1V>

More documents and user manuals: <https://manuals.plus/>