

Jinrrenyi Chair Lift

Jinrrenyi Electric Chair Lift User Manual

Model: Chair Lift

1. INTRODUCTION

This manual provides essential instructions for the safe and effective use, assembly, operation, and maintenance of your Jinrrenyi Electric Chair Lift. Designed to assist elderly individuals and those with mobility challenges, this chair lift offers support for standing, sitting, and moving between rooms. Please read this manual thoroughly before initial use and retain it for future reference.

2. SAFETY INFORMATION

Your safety is paramount. Adhere to all warnings and instructions to prevent injury or damage.

- **Maximum Weight Capacity:** Do not exceed the maximum weight capacity of 573 lbs (260 kg). Overloading can cause instability and injury.
- **Assembly:** Always ensure the chair lift is fully and correctly assembled according to the instructions before use.
- **Moving Parts:** Keep hands, fingers, clothing, and other objects clear of all moving parts during operation.
- **Safety Belts and Handrails:** Always use the integrated safety seat belts and handrails when operating the chair lift or when seated. Ensure the safety belt is securely fastened.
- **Wheel Locks:** Engage the wheel locks on the casters when the chair lift is stationary to prevent unintended movement.
- **Stable Surface:** Use the chair lift on a flat, stable, and level surface. Avoid uneven terrain or inclines.
- **Supervision:** Users with cognitive impairments or severe mobility limitations should be supervised during use.
- **Power Source:** Use only the provided power adapter and ensure it is connected to a suitable power outlet.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is a mobility aid and not a medical device for treatment.

3. PRODUCT OVERVIEW

The Jinrrenyi Electric Chair Lift is designed for comfort, stability, and ease of use, providing assistance for daily activities.

Electric Chair Lift with Wheels



Image: Multi-functional Electric Chair Lift. This image illustrates the chair lift being used in both a living room setting, with an elderly woman standing beside it, and a bedroom setting, with an elderly man using it to get up from a bed.

Key Features:

- **High Weight Capacity:** Supports up to 573 lbs, suitable for a wide range of users.
- **Padded Armrests:** Designed to assist with standing and provide support, reducing strain on joints.
- **PU Cushion:** Soft, breathable, and comfortable seating material.
- **Sturdy Frame:** Constructed from alloy steel for durability and stability.
- **Anti-slip Wheels:** Four high-performance 360-degree swivel casters for effortless movement across various floor types, with integrated brakes for stability.
- **One-Button Operation:** Simple remote control for lifting and lowering functions.
- **Safety Strap:** Built-in safety strap to secure users firmly during movement.
- **Fold-Flat Design:** Seat and armrests fold down to 10.6 inches height for compact storage.



Image: 4 Swivel Casters for Effortless Movement. An elderly woman is shown pushing the chair lift, highlighting the 360-degree swivel casters, two of which include brakes for secure positioning.



Image: More Secure Load-bearing. An elderly woman is seated in the chair lift with a blue safety belt fastened, demonstrating how the belt secures the user. Front and back views of the safety belt are also shown.



Image: Comfort & Support. This image highlights the soft, cloud-like material of the PU cushion and illustrates how the armrests fold 90 degrees when idle, offering stable support during operation.

4. SETUP INSTRUCTIONS

The Jinrenyi Electric Chair Lift is designed for tool-free, simple assembly, allowing you to set it up quickly.

Assembly Steps:

1. **Unpack Components:** Carefully remove all components from the packaging. Ensure all parts are present: main frame, seat cushion, four wheels, and remote control.
2. **Attach Wheels:** Insert the four swivel casters into the designated slots on the base of the chair lift frame. Ensure they click securely into place. Two of the casters have brakes.
3. **Install Seat Cushion:** Place the PU seat cushion onto the chair frame. Ensure it is properly aligned and secured.
4. **Connect Remote Control:** Connect the remote control cable to the designated port on the chair lift's control panel.
5. **Initial Charge:** Before first use, fully charge the chair lift's battery. Connect the power adapter to the charging port and plug it into a standard electrical outlet. An 8-hour charge provides approximately 80 lifting operations.

Electric Chair Lift with Wheels



Image: Product Dimension. This image displays the chair lift's dimensions (35.4"D x 20.5"W x 59"H) and its folded state (10.6" height), indicating its compact storage capability.

5. OPERATING INSTRUCTIONS

Operating your Electric Chair Lift is straightforward with its user-friendly remote control.

5.1. Lifting and Lowering

1. **Positioning:** Ensure the chair lift is positioned correctly and stably before use. If moving the chair, unlock the casters.
2. **Seating:** Carefully sit on the chair lift, ensuring you are centered and stable.
3. **Secure Safety Belt:** Fasten the safety belt securely around your waist or chest as appropriate.
4. **Engage Wheel Locks:** If remaining stationary, engage the wheel locks on the casters to prevent any movement.
5. **Using the Remote Control:**
 - Press the **UP arrow button** to slowly raise the seat to a standing position.
 - Press the **DOWN arrow button** to slowly lower the seat to a sitting position.
6. **Standing/Sitting:** Use the padded armrests for additional support as the chair lifts or lowers. Once at the desired height, carefully stand up or sit down.

Your Go-To for Multiple Usage Scenarios

Where to Use



Home



Medical Institutions



Elderly Care Facilities

Who Needs to Use



Elderly people



Postoperative



Caregivers

Image: Two Buttons for Lifting and Lowering Operation. An elderly man is shown using the remote control to adjust the chair's height, demonstrating both the lifting and lowering functions.

5.2. Mobility and Transfer

The chair lift's swivel casters allow for easy movement between rooms.

1. **Unlock Casters:** Ensure all wheel locks are disengaged before attempting to move the chair.
2. **Safe Movement:** Gently push or pull the chair lift to the desired location. For user safety, it is recommended that a caregiver assist with movement if the user is seated.
3. **Re-engage Locks:** Once at the new location, engage the wheel locks to secure the chair lift in place.

5.3. Charging the Battery

The chair lift is powered by a rechargeable battery.

- **Charging Port:** Locate the charging port on the control panel (refer to the image below).
- **Connect Charger:** Plug the power adapter into the charging port and then into a standard wall outlet.
- **Charging Time:** A full charge takes approximately 8 hours and provides over 80 lifting operations.
- **Indicator:** The remote control or control panel may have an indicator light to show charging status.

Your Go-To for Multiple Usage Scenarios

Where to Use



Home



Medical Institutions



Elderly Care Facilities

Who Needs to Use



Elderly people



Postoperative



Caregivers

Image: User-friendly Remote Control. An elderly man is seated in the chair lift with the remote control. A detailed inset shows the motor connection port, remote control port, and charging port on the chair's control panel.

6. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your chair lift.

- **Cleaning:** Wipe down the PU cushion and frame with a damp cloth and mild soap. Avoid abrasive cleaners or solvents. Ensure all surfaces are dry before use.
- **Caster Inspection:** Periodically check the wheels for any debris, wear, or damage. Clean any hair or lint from the caster axles to ensure smooth movement.
- **Structural Check:** Regularly inspect the frame, armrests, and seat for any signs of damage, loose screws, or structural integrity issues. Tighten any loose fasteners.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Charge the battery regularly, even if the chair lift is not in constant use.
- **Remote Control:** Keep the remote control clean and dry. Replace batteries if necessary (if applicable, for wireless remotes).

7. TROUBLESHOOTING

If you encounter issues with your chair lift, refer to the following common problems and solutions.

Problem	Possible Cause	Solution
Chair does not lift or lower.	Battery is low or depleted. Power cable is disconnected. Remote control issue.	Charge the battery for 8 hours. Ensure all cables are securely connected. Check remote control connection.

Problem	Possible Cause	Solution
Chair moves slowly or inconsistently.	Low battery. Obstruction in lifting mechanism.	Recharge the battery. Inspect for any obstructions and remove them carefully.
Wheels do not move smoothly.	Wheel locks engaged. Debris in casters.	Disengage all wheel locks. Clean any hair or debris from the caster axles.
Unusual noises during operation.	Loose components. Motor issue.	Inspect and tighten any loose screws or connections. If noise persists, discontinue use and contact support.

If the problem persists after attempting these solutions, please contact customer support.

8. SPECIFICATIONS

Detailed technical specifications for the Jinrenyi Electric Chair Lift.

Feature	Detail
Model Name	Chair Lift
Product Dimensions (Unfolded)	35.4"D x 20.5"W x 59"H (90cm D x 52cm W x 150cm H)
Folded Height	10.6 inches (27 cm)
Item Weight	72 Pounds (32.6 kg)
Maximum Weight Recommendation	573 Pounds (260 kg)
Frame Material	Alloy Steel
Seat Material	Polyurethane (PU) Cushion
Fill Material	Polyurethane Foam
Casters	4 Swivel Casters (2 with brakes)
Lifting Speed	0.15 m/s
Charging Time	Approximately 8 hours
Uses per Full Charge	Over 80 operations
Required Assembly	Yes (Tool-free)
Indoor/Outdoor Usage	Indoor

Electric Chair Lift with Wheels



Image: Dimensions. This image provides a clear visual representation of the chair lift's key measurements, including its height, width, depth, and folded height for storage.

9. WARRANTY AND SUPPORT

9.1. Warranty Information

The Jinrenyi Electric Chair Lift comes with a limited manufacturer's warranty. For specific details regarding coverage, duration, and terms, please refer to the warranty information provided with your product packaging or contact the seller directly.

9.2. Customer Support

If you have any questions, require assistance with assembly, operation, or troubleshooting, or need to claim warranty service, please contact your retailer or the manufacturer, Jinrenyi, through their official customer service channels. Please have your product model number (Chair Lift) and purchase information ready when contacting support.

10. USAGE SCENARIOS

The Jinrenyi Electric Chair Lift is versatile and suitable for various environments and user needs.



Image: Your Go-To for Multiple Usage Scenarios. This image illustrates various places where the chair lift can be used (Home, Medical Institutions, Elderly Care Facilities) and who can benefit from its use (Elderly people, Postoperative patients, Caregivers).

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question