

- Manuals.plus /
- › MEGELIN /
- › MEGELIN Multi-Function Therapy Belt User Manual

MEGELIN Multi-Function Therapy Belt

MEGELIN Multi-Function Therapy Belt User Manual

Model: Multi-Function Therapy Belt

IMPORTANT SAFETY INFORMATION

Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult with a healthcare professional before using this device, especially if you have pre-existing medical conditions, are pregnant, or have implanted medical devices. Do not use on broken or irritated skin. Discontinue use if you experience discomfort or adverse reactions.

PRODUCT OVERVIEW

The MEGELIN Multi-Function Therapy Belt is designed for at-home use, incorporating heat, Electrical Muscle Stimulation (EMS), and light therapy to support various body areas. This device is intended to assist in muscle conditioning and skin elasticity improvement.



Image: The Megelin Multi-Function Therapy Belt in grey, shown in use around the waist. The belt features a central control unit and is designed for body sculpting and recovery.

Key Features:

- **Multi-functional Therapy:** Combines heat, EMS, and light therapy.
- **Adjustable Fit:** Designed with adaptable straps for various body areas.
- **Wireless Operation:** Provides convenience and portability without the need for gels.
- **Targeted Support:** Assists in strengthening core and glute muscles and improving skin elasticity.

Sculpt & Recover Your Body At-Home

Achieve a slimmer waist with Heat + EMS + Light

Boost Fat Metabolism

Reduce Cellulite

Sculpt Waist & Core

Improve Skin Elasticity



Image: An illustration highlighting the belt's functions, including heat, EMS, and light therapy for metabolism, cellulite reduction, waist sculpting, and skin elasticity.

COMPONENTS

The Megelin Multi-Function Therapy Belt system includes the following main components:

- **Function Unit:** The central control module containing the electronics and interface.
- **Long Strap:** A longer strap for larger body areas or extended reach.
- **Medium Strap:** A shorter strap for more compact body areas.

Adjustable Fit for Every Body Area




Image: A visual breakdown of the product components, including the main function unit and two different sized straps (long and medium) for versatile application.



Image: The therapy belt with its function unit separated from the strap, illustrating how the components connect.

SETUP

1. **Charge the Function Unit:** Before first use, fully charge the main function unit using the provided charging cable. Ensure the unit is powered off during charging.
2. **Attach the Strap:** Select the appropriate strap (long or medium) for the desired body area. Connect the function unit securely to the designated area on the strap. Ensure a firm connection.
3. **Position the Belt:** Wrap the belt around the target body area (e.g., waist, abdomen, thigh, arm). Adjust the strap for a snug but comfortable fit. The belt is designed to be used without gels.

OPERATING INSTRUCTIONS

Follow these steps to operate your Megelin Multi-Function Therapy Belt:

1. **Power On:** Press and hold the power button on the function unit until the indicator lights illuminate.

2. **Select Mode:** Use the mode selection button to cycle through available therapy modes (Heat, EMS, Light, or combinations). Refer to the unit's display for current mode indication.
3. **Adjust Intensity:** Use the intensity adjustment buttons (+/-) to set the desired level for EMS and heat functions. Start with a low intensity and gradually increase to a comfortable level.
4. **Begin Session:** Once settings are configured, the therapy session will begin. The duration of each session may be pre-set or adjustable depending on the mode.
5. **Power Off:** After your session, press and hold the power button to turn off the device.



Image: A detailed view of the belt's light therapy function, showing multiple red LEDs activated on the control panel.

Usage Examples:



Image: The therapy belt positioned around the waist, demonstrating a common application area.



Image: The therapy belt applied to the thigh, showcasing its adaptability for different body parts.

MAINTENANCE

- **Cleaning:** After each use, wipe the function unit and straps with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water. Ensure the device is completely dry before storage.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and extreme temperatures. Keep it out of reach of children.
- **Battery Care:** To maintain battery life, charge the device regularly, even if not in frequent use. Avoid fully depleting the battery for extended periods.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not power on.	Low battery; improper power button press.	Charge the device fully. Press and hold the power button for a few seconds.
No sensation from EMS.	Low intensity setting; poor contact with skin.	Increase intensity gradually. Ensure the belt is snug against the skin.

Problem	Possible Cause	Solution
Heat function not working.	Heat mode not selected; low battery.	Select the heat mode. Ensure the device is sufficiently charged.

If the problem persists after attempting these solutions, please contact customer support.

SPECIFICATIONS

Model: Multi-Function Therapy Belt
ASIN: B0GJ4F7XFL
Package Dimensions: 15.67 x 6.89 x 2.17 inches
Weight: Approximately 2.09 Pounds
Manufacturer: SHENZHEN ZHENXING RUITONG TECHNOLOGY CO., LTD.
First Available: January 22, 2026

WARRANTY AND SUPPORT

This MEGELIN product comes with a **365-day warranty** from the date of purchase. This warranty covers manufacturing defects and malfunctions under normal use. It does not cover damage caused by misuse, accidents, unauthorized modifications, or improper maintenance.

For warranty claims, technical support, or any questions regarding your Multi-Function Therapy Belt, please contact MEGELIN customer service.

Contact Information: Please refer to the contact details provided with your product packaging or visit the official MEGELIN website for support.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0GJ4F7XFL>

More documents and user manuals: <https://manuals.plus/>