

dreame Chrona Photon Anti-Aging Skincare Mask

DREAME Chrona Photon Anti-Aging Skincare Mask User Manual

Model: Chrona Photon Anti-Aging Skincare Mask

- Introduction
- Safety Information
- Package Contents
- Setup
- Operation
- Light Modes
- Maintenance
- Troubleshooting
- Specifications
- Warranty & Support

1. INTRODUCTION

The DREAME Chrona Photon Anti-Aging Skincare Mask is an advanced light therapy device designed for facial skin care. It utilizes a combination of 4+1 synergistic light modes to address various skin concerns, including firming, brightening, soothing, and spot reduction. This manual provides essential information for the safe and effective use of your device.

2. SAFETY INFORMATION

Please read all safety instructions before using the DREAME Chrona Photon Anti-Aging Skincare Mask.

- Consult a physician before use if you have any medical conditions, are pregnant, or are taking photosensitizing medications (e.g., Vitamin A derivatives like retinol).
- Do not use on open wounds, broken, or severely irritated skin.
- Avoid direct eye exposure to the LED lights. The mask is designed with eye protection, but it is recommended to keep your eyes closed during treatment.
- Keep the device out of reach of children.
- Do not immerse the device in water or other liquids.
- Discontinue use if you experience any discomfort, redness, or irritation.
- Use only the provided charging cable and batteries.

3. PACKAGE CONTENTS

Verify that all items are present in your package:

- DREAME Chrona 3D LED Mask
- Remote Controller
- Spare Battery (x2)
- Quick Start Guide

- Charging Cable (USB-C)
- User Manual
- Storage Bag



Chrona 3D LED Mask×1



Controller×1



Spare Battery×1



Quick Start Guide×1



Charging Cable×1



User Manual×1



Storage Bag×1

Package List

Image: All components included in the DREAME Chrona Photon Anti-Aging Skincare Mask package, including the mask, remote, batteries, cables, and storage bag.

4. SETUP

4.1 Charging the Batteries

The DREAME Chrona mask comes with two swappable batteries. While they may arrive partially charged, it is recommended to fully charge them before first use. Each battery takes approximately 30 minutes to fully charge via the provided USB-C cable.

How to Use?



Step1 Cleanse

Remove makeup and wash your face.



Step2 Wear

Secure the mask with the adjustable strap.



Step3 Power On

Long press the power button.



Step4 Set

Use the remote to choose mode and time based on what your skin needs today.



Step5 Relax

Enjoy a hands-free treatment for 10-15 minutes.




Step6 Skincare

Remove the mask and follow your regular skincare routine.

Image: Illustration of the dual swappable batteries, highlighting their convenience for continuous use.

4.2 Inserting the Battery

Locate the battery slot on the side of the mask. Gently insert one of the charged circular batteries into the slot. The battery connects magnetically and should fit securely.



Full face coverage No Light Loss

Image: A hand demonstrating how to insert the battery into the DREAME Chrona mask.

4.3 Pairing the Remote Controller

To pair the remote controller with the mask, ensure the mask has a charged battery inserted and is powered on. Hold the power button on the remote until the lights on the mask sync, indicating successful pairing. The remote allows for easy adjustment of light modes and intensity.



Salon-Grade Glow At Home

Image: The DREAME Chrona mask shown alongside its remote controller, illustrating wireless operation.

5. OPERATING INSTRUCTIONS

The DREAME Chrona mask is designed for a 10-minute full-face light care session. Follow these steps for optimal use:

1. **Cleanse:** Thoroughly clean your face to remove makeup, dirt, and oils.
2. **Wear:** Secure the mask comfortably on your face using the adjustable strap. Ensure the mask fits snugly but not too tightly.
3. **Power On:** Press and hold the power button on the mask or remote controller to turn on the device.
4. **Set Mode & Time:** Use the remote controller to select your desired light mode and session duration (5, 10, or 15 minutes). The mask automatically recalls your last settings for convenience.
5. **Relax:** Enjoy a hands-free treatment session. You can perform other activities during this time.
6. **Skincare:** After the session, remove the mask and continue with your regular skincare routine.

Package List




Image: A visual guide demonstrating the step-by-step process of using the DREAME Chrona mask.

5.1 Video Guide: Daily Facial Care

Video: An official DREAME video demonstrating the Chrona Photon Skincare Mask for daily facial care, showing its ease of use and hands-free operation.

5.2 Video Guide: Easy Beauty Routine

Video: An official DREAME video showcasing the modern LED Skincare Mask and how it integrates into an easy beauty routine.

6. LIGHT MODES AND BENEFITS

The DREAME Chrona mask features 4+1 synergistic light modes, each targeting specific skin needs:

- **Rejuvenate Mode (Red 630nm + NIR 850nm):** Firms skin, reduces fine lines, and provides anti-aging support.
- **Calm Mode (Blue 465nm):** Soothes deep inflammation, boosts microcirculation, and helps clear blemishes.
- **Brighten Mode (Yellow 590nm):** Promotes an even skin tone and natural glow, enhancing radiance.
- **Balance Mode (Green 530nm):** Calms visible redness and balances complexion.
- **Restore Mode (Multi-Light):** Combines all four lights for comprehensive skin renewal.

Skin redness
24.73%

Skin oil
6.03%



Real participant images from a 28-day SGS study (n=32; Report No. QDCPCH25001541-01_CN). Results vary.

Image: Visual representation of the five distinct light modes (Red, Blue, Yellow, Green, and Multi-Light) available on the DREAME Chrona mask.

High-power LEDs for luminous skin.

40 mW/cm²

Optimal power. Superior results.

280 LEDs

4-core design. Greater density.

Gold-Wire Technology

Stable energy that lasts.





Image: An illustration detailing the targeted benefits of each light mode for different skin concerns.

7. MAINTENANCE

7.1 Cleaning the Mask

After each use, it is recommended to wipe down the interior and exterior of the mask with an alcohol wipe or a soft cloth dampened with a mild, non-abrasive cleaner. This helps prevent buildup and maintains hygiene.

7.2 Storage

Store the DREAME Chrona mask in the provided storage bag when not in use. This protects the device from dust and damage, ensuring its longevity.



just
338g



Wireless Freedom Ultra-Soft & Breathable

Image: The DREAME Chrona mask neatly stored in its protective bag.

8. TROUBLESHOOTING

If you encounter any issues with your DREAME Chrona mask, refer to the following common solutions:

Problem	Possible Cause	Solution
Mask has a plastic odor	New product odor from packaging materials.	Leave the mask out in a well-ventilated area for 24-48 hours to allow the odor to dissipate.
Light is visible around the eyes	Minor light spillage is normal due to ergonomic design.	The mask includes eye protection. For added caution, keep your eyes closed during treatment.
Mask feels heavy or pulls downward	Strap may be too loose or mask not positioned correctly.	Adjust the head strap to ensure a snug but comfortable fit. The mask is designed to be lightweight (approx. 338g).

Problem	Possible Cause	Solution
Mask does not power on	Battery is not charged or not inserted correctly.	Ensure the battery is fully charged and securely inserted into the mask. Try using the spare battery.
Remote not controlling mask	Remote not paired or battery low.	Re-pair the remote by holding its power button while the mask is on. Ensure the remote's battery is charged.

9. SPECIFICATIONS

- **Model:** Chrona Photon Anti-Aging Skincare Mask
- **Light Modes:** 4+1 Synergistic Lights (Red 630nm + NIR 850nm, Blue 465nm, Yellow 590nm, Green 530nm, Multi-Light)
- **LEDs:** 280 high-energy LEDs
- **Optimal Power:** 40 mW/cm²
- **Session Duration:** 5, 10, or 15 minutes (recommended 10 minutes)
- **Fit:** 3D Ergonomic Fit, soft liquid silicone
- **Weight:** Approximately 338g
- **Power Source:** Dual Swappable Nonstandard Batteries (included)
- **Control:** Wireless Remote Controller
- **Dimensions (Package):** 12.05 x 10.98 x 5.71 inches
- **UPC:** 850079600300
- **Manufacturer:** dreame

10. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the official DREAME website or contact their customer service directly. Details can typically be found in the Quick Start Guide or on the product packaging. Keep your purchase receipt for warranty claims.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0GJZTTZ3C>

More documents and user manuals: <https://manuals.plus/>