



bruin BIOMETRICS OTH-SEM250-0344 Provizio Sem Scanner S Sensor Instruction Manual

[Home](#) » [bruin BIOMETRICS](#) » bruin BIOMETRICS OTH-SEM250-0344 Provizio Sem Scanner S Sensor Instruction Manual 



bruin BIOMETRICS OTH-SEM250-0344 Provizio Sem Scanner S Sensor Instruction Manual



Contents

- [1 Features](#)
- [2 INSTRUCTIONS FOR USE](#)
- [3 Support](#)
- [4 Documents / Resources](#)
 - [4.1 References](#)
- [5 Related Posts](#)

Features



INSTRUCTIONS FOR USE

Intended Use: The single use sensor is to be used with Provizio® SEM Scanner S. The Provizio® SEM Scanner S is intended to be used by healthcare professionals as an adjunct to the standard of care when assessing the heels and sacrum of patients who are at increased risk for pressure ulcers.



Instructions: To open pouch, split apart at tabs. Remove the single use sensor and attach to the Provizio® SEM Scanner S.

Warnings and General Precautions:

This product must be used in accordance with the Provizio® SEM Scanner S instructions for use and related labeling.

References:

Provizio® SEM Scanner S Manual available at www.sem-scanner.com

Support



Bruin Biometrics, LLC

10877 Wilshire Blvd.,
Ste 1600
Los Angeles, CA 90024



Emergo Europe B. V.

Prinsessegracht 20
2514 AP, The Hague
THE NETHERLANDS
Tel: +31 70 345 8570



02637



100 Sensors

Australian Sponsor:

Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia



Documents / Resources



[bruin BIOMETRICS OTH-SEM250-0344 Provizio Sem Scanner S Sensor](#) [pdf] Instruction Manual
OTH-SEM250-0344 Provizio Sem Scanner S Sensor, OTH-SEM250-0344, Provizio Sem Scanner S Sensor, Provizio Sem Scanner, Sem Scanner, Scanner

References

- [scanner.com](#)
- [Provizio® SEM Scanner - risk assessment device pressure ulcer](#)

[Manuals+.](#)