

# **DETECTO 0065-0824-08 Portable Stretcher Scale Instruction Manual**

Home » DETECTO » DETECTO 0065-0824-08 Portable Stretcher Scale Instruction Manual



### **Contents**

- 1 DETECTO 0065-0824-08 Portable Stretcher Scale
- **2 ASSEMBLY**
- 3 Documents / Resources
- **4 Related Posts**

# Detecto

## **DETECTO 0065-0824-08 Portable Stretcher Scale**



Before starting to assemble, please unpack and remove all wrappings and fillers. Carton contains:

Item Description	Part Number	Qty.
8500/8550 Guide Rail, Wheel Side	0065-0824-08	2
8500/8550 Guide Rail, Handle Side	0065-0825-08	2
Screw Truss-Head Thread Cutting, 8-32 x .375	6021-1257	10

#### **TOOLS REQUIRED**

• No. 2 Philips screwdriver (power screwdriver recommended)

### **ASSEMBLY**

Refer to Figures No. 1, 2, and 3, then follow the instructions for attaching the Guide Rails.

- 1. Place the platform on the floor in the storage position.
- 2. Begin by holding a Wheel Side Guide Rail (with a caution label toward you) against the side of the platform (this is the longer of the two sets of guide rails). Install (2) thread-cutting screws through the holes in the guide rail, and into the holes of the platform. Securely tighten the screws. Refer to Figure No. 1.
- 3. Repeat step 2 for the other Wheel Side Guide Rail.
- 4. Next, hold a Handle Side Guide Rail (with a caution label toward you) against the side of the platform (this is the shorter of the two sets of guide rails). Install (3) thread-cutting screws through the holes in the guide rail, and into the holes of the platform. Securely tighten the screws. Refer to Figure No. 2.
- 5. Repeat step 4 for the other Handle Side Guide Rail.
- 6. The finished assembly should look like Figure No. 3







Figure No. 3

Printed in USA 0065-0832-0M Rev B 12/19

# **Documents / Resources**



**DETECTO 0065-0824-08 Portable Stretcher Scale** [pdf] Instruction Manual 0065-0824-08, Portable Stretcher Scale, 0065-0824-08 Portable Stretcher Scale, Stretcher Scale, Scale

Manuals+,